Agilent 8850

PDD detektor

LUMA detektor

8890 analyzáry

SRA Instruments

Daniel Sander

NPI

12.11.2024

Altium

Agilent 8850

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Application Note
Energy and Chemicals

The Analysis of Monocyclic Aromatic Hydrocarbons by ASTM D7504 on the Agilent 8850 GC System

Conventional and high-speed approaches using helium or hydrogen carrier gas

Author

Scott Hoy
Agilent Technologies, Inc.

Abstract


This application note highlights the speed and precision of the Agilent 8850 single-channel GC applied to both a conventional separation (39 minutes) and fast separation (6.05 minutes) of common aromatic contaminants described in ASTM D7504.¹ Also shown is the use of the Agilent method translator software to convert the conventional method into the fast method for users interested in exploring additional variations of ASTM D7504. The precision of both methods was evaluated using three different purchased standard mixtures analyzed over 20 consecutive replicate injections each. The system yielded excellent precision for both methods with calculated concentration RSDs below 1.0% for all non-trace-level compounds. The fast method yielded a benzene peak in the p-xylene standard with an average (n = 20) RMS signal-to-noise ratio (S/N) of 60 at a concentration of 14 ppmw, indicating an approximate benzene MDL of less than 1.0 ppmw.

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Agilent 8850

Method Translator
(OL CDS Workstation, MassHunter)

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Změna kolony
(20m x 0.18mm x 0.18µm místo 60m x 0.32mm x 0.5µm)

The screenshot displays the 'Method Translator' software interface. It compares 'Original Method Parameters' (Gas: He) with 'Calculated Method Parameters' (Gas: H₂). The interface includes sliders for various parameters and a table for temperature ramps.

Parameter	Original Method Parameters (He)	Calculated Method Parameters (H ₂)
Length (m)	60 m	20 m
Inner Diameter (µm)	320 µm	180 µm
Film Thickness (µm)	0.50 µm	0.18 µm
Phase Ratio	159.25	249.25
Inlet Pressure (gauge)	17.587 psi	18.311 psi
Outlet Flow (mL/min)	2.1 mL/min	1.4766 mL/min
Average Velocity (cm/s)	29.248 cm/sec	63.72 cm/sec
Outlet Pressure (abs)	34.896 psi	34.896 psi
Holdup Time	5.419 min	0.52313 min
Outlet Velocity (cm/s)	48.628 cm/sec	108.06 cm/sec

#	Ramp Rate (°C/min)	Final Temp (°C)	Final Time (min)
Init		60	10
1	5	150	0
2	50	200	10

Total Run Time: 99.00 min

#	Ramp Rate (°C/min)	Final Temp (°C)	Final Time (min)
Init		60	0.98
1	51.147	150	0
2	511.47	200	0.98

Total Run Time: 3.82 min

Pressure Units: PSI


Original Column Capacity: 7.01

Translated Column Capacity: 0.61

The column capacity of the translated method is 9% of the original column capacity. You may need to adjust your injection volume.

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