

Specifications

Item	Description
Analysis functions	Univariate analysis t-test Mann-Whitney U-test ANOVA (analysis of variance)
	Multivariate analysis PCA (Principal Component Analysis) PLS-DA
	Discriminant analysis Support Vector Machine (SVM) Random Forest
	Other Dynamic grouping
Display functions	Multivariate analysis Peak Matrix Box Plot ROC AUC Score/Loading Plot Dendrogram
	Discriminant analysis Discriminant analysis results (group and score) Unknown samples superimposed on Score Plot
Input-output data	Input Chromatogram summary file ASCII format peak lists JCAMP format peak lists mzML format peak lists (only supports Centroid mode, 1 scan/per file, data uncompressed) Note: File name must only contain half-width alphanumeric characters.
	Output Peak lists (txt format) Analysis results (xslm format) Graph screenshot
Operating system	OS Windows® 11 Professional 64bit (English/Chinese) Windows® 10 Professional 64bit (English/Chinese)



eMSTAT Solution, DPiMS, Smart Metabolites Database, Smart Aroma Database, and GCMS-TQ are trademarks of Shimadzu Corporation or its affiliated companies in Japan and/or other countries.
 Windows is either a registered trademark or a trademark of Microsoft Corporation in the United States and/or other countries.

 Statistical Analysis Software for Analytical Instruments
eMSTAT Solution


Enabling statistical analysis of chromatography and direct ionization mass spectrometry data regardless of user experience

eMSTAT Solution™



Sample group differentiation—such as when comparing against benchmark products in the product development phase or in sample pass/fail assessment during quality control—can be difficult when relying on conventional qualitative and quantitative analytical techniques. The process typically involves multivariate analysis using large amounts of chromatogram or MS spectrum data to search for markers that can help differentiate between sample groups, and to visualize sample characteristics using these markers. It also requires modeling using training data and discriminant analysis to scientifically determine the characteristics of unknown samples, including their condition (e.g., deterioration), type, and place of origin. The eMSTAT Solution statistical analysis software is equipped with a range of features that make multivariate analysis more accessible and user-friendly for everyone—even when analyzing large volumes of various analytical data (chromatograms and MS spectra).



Features



- Statistical Analysis Mode for Intuitive Data Analysis Operations**
 The software's Statistical Analysis Mode contains extensive, easy-to-use functions for univariate and multivariate analysis that enable even users without statistical analysis expertise to identify marker peaks for sample differentiation and discrimination.
- Dynamic Grouping Function for Flexible Sample Grouping**
 Leveraging the software's quality information library, the dynamic grouping function facilitates marker discovery by enabling flexible sample grouping in various combinations.
- Discriminant Analysis Mode Supports Grouping of Unknown Samples**
 Discriminant models for discrimination of unknown samples can be created using marker peaks identified in Statistical Analysis Mode. With single-step mode switching, statistical data obtained in Statistical Analysis Mode can be used in Discriminant Analysis Mode.

Workflow Using eMSTAT Solution



Statistical Analysis Mode for Intuitive Data Analysis Operations

Proper statistical analysis typically requires a certain degree of user expertise and poses a steep learning curve. The eMSTAT Solution statistical analysis software eliminates these challenges with a suite of intuitive operating features. The software enables a wide range of analyses from sample differentiation and identification of useful marker peaks in Statistical Analysis Mode to discriminant modeling of unknown samples in Discriminant Analysis Mode.




Dynamic Grouping Function for Flexible Sample Grouping

In sample group differentiation, analysis samples can be grouped according to attribute (raw materials, reaction conditions, normal/abnormal, manufacturer, type, etc.) and then searched for characteristic group markers. The eMSTAT Solution software's Dynamic Grouping Function effectively enables the creation of new sample groups while retaining existing group information.




Analytical Conditions
Instrument: LCMS-8045/8050/8060NX
Sample: Commercially available yogurt
Using the ion-pair-free LC/MS/MS method included in Shimadzu's LC/MS/MS Method Package for Primary Metabolites Ver. 3, the system can comprehensively analyze a total of 153 compounds, including sugars, in addition to hydrophilic metabolites essential for metabolomic analysis, such as amino acids, organic acids, nucleosides, and nucleotides.



Dynamic grouping


Rearranged Groups

Group Name	A+B	C+D
Company X	+	-
Company Y	-	+



Discriminant Analysis Mode Supports Grouping of Unknown Samples

Identifying Degraded Synthetic Polymers



Polystyrene (heated and unheated) polymers were separated into two groups (heated or unheated) by multivariate analysis (PLS-DA) of the MALDI mass spectrum (Score Plot). A Loading Plot can be used to confirm which peaks (marker peaks) affect the differences between the two groups.

Experimental Conditions
 • Instrument: MALDI-8020
 • Sample: Polystyrene 1,000 (heated for 2 hours at 200 °C or unheated)
 • Matrix: Dithranol
 • Cationizing Agent: AgTFA
 • Sample Preparation: The polystyrene (heated and unheated) dissolved in THF, matrix solution, and cationizing agent were spotted onto a sample plate.

Using marker peaks identified by multivariate analysis to create a discriminant model and discriminate between heated and unheated polymers in a polystyrene mass spectrum, obtained separately, resulted in correct discrimination of all polymers. By using eMSTAT Solution in combination with MALDI mass spectrometry, which can easily measure samples with large molecular weight synthetic polymers and a wide variety of other samples, such as protein, fat, or sugar samples, can be easily differentiated.




Discriminant Analysis of Beef Meats

In addition to chromatogram and MS spectrum data, the eMSTAT Solution software can read various data formats (JCAMP, ASCII, mzML, etc.), supporting the integrated analysis and multi-faceted assessment of various types of data.



Analytical Conditions
 Instrument: GCMS-TQ8040 NX
 Sample: Beef meat
 Using the Smart Metabolites Database™ Ver. 2, this system is capable of analyzing about 500 primary metabolites, including organic acids, sugars, nucleic acids, fatty acids, and amino acids. Using the Smart Aroma Database™, it can also analyze about 500 aroma compounds with SPME Arrow.

A discriminant analysis was successfully performed on an unknown meat to determine its farm of origin based on approximately 600 variables (i.e., GC/MS, UV, meat lipid measuring device, microwave-drying moisture meter, and visual assessment). In addition to determining the origin, the discriminant model revealed the characteristic variables of Farmer 2 (e.g., the roast-smelling aroma compound benzothiazole, the primary metabolite linoleic acid).



Meat lipid measuring device
 Oleic acid and monounsaturated fatty acids were measured.

Visual assessment
 Assessment of grading score, marbling, BMS, BFS, BCS, etc.


UV-Visible spectrophotometer
 Determination of glycogen concentration

Microwave-drying moisture meter
 Analysis of crude fat content, crude protein content, and water content

Approximately 600 variables

Random forest discriminant model



Easy Differentiation for Beef Classification



Extracts from commercial beef (Tasmanian and A5/A4 grade Wagyu beef) were analyzed in a DPiMS system mass spectrometer. The resulting spectra were then analyzed by PLS discriminant analysis. A Score Plot confirms grouping into three groups and a Loading Plot confirms which metabolite peaks affect grouping.

Experimental Conditions
 • Instrument: DPiMS system
 • Sample: Edible beef
 • Ionization Solution: 50% aqueous ethanol solution
 • Sample Preparation: 10 µL of the suspension from homogenizing the beef for human consumption in the 50% aqueous ethanol solution was dripped onto a sample plate.

With eMSTAT Solution, spectra obtained by convenient metabolite analysis in a DPiMS system spectrometer can be used to easily differentiate between differences in food, plant, and other samples, and screen for information about metabolites that contribute to differentiation.



The number of characteristic variables in each farm (i.e., ANOVA p-value of < 5 %)

Select	<input type="checkbox"/>
Name	Unknown (Farm 2)
Group	Farm 2
Score	70

