

AUTOMATED ONLINE SPE/TANDEM MS ANALYSIS OF TRACE ORGANIC CONTAMINANTS IN DRINKING WATER

Agilent Ultivo Triple Quadrupole LC/MS System




Figure 1A. Agilent Ultivo Triple Quadrupole LC/MS.



Figure 1B. Agilent 1290 Infinity FlexCube.

Introduction

A variety of organic molecules originating from pesticides, pharmaceuticals, personal care products, and industrial chemicals used in daily life can be introduced into the environment through domestic and industrial wastewater sources. These molecules are a subclass of compounds collectively known as Trace Organic Contaminants (TOrcs), commonly found in a variety of potable and non-potable water resources. Although not acutely dangerous as highly dilute individual entities, the effects of long-term synergistic exposure due to mixtures of these compounds are still to be determined¹. While characterization for long-term and chronic toxicity are yet to be established, it is important to monitor their presence in current water resources².

Using the online SPE capabilities of the Agilent 1290 Infinity Flexible Cube with the analytical power of the Agilent Ultivo Triple Quadrupole LC/MS, we have demonstrated a rapid and sensitive method for automated analysis of trace organic contaminants in drinking water without tedious offline enrichment.

For more information, visit:

www.agilent.com/chem/Ultivo



Agilent Technologies

Online SPE, HPLC, and MS parameters

Table 1 lists Online SPE, HPLC, and MS instrument parameters.

Table 1. Instrument conditions.

Online SPE conditions											
SPE Cartridge	Agilent ZORBAX Extend C18, 4.6 × 12.5 mm, 5 μm										
Agilent 1290 Infinity Flexible Cube solvents	A (Load): H ₂ O with 0.1 % v/v formic acid B (Elute): 1:1 ACN/IPA										
Sample volume	750 μL										
SPE Elution gradient	0 minutes: Load - Pump for 240 seconds H ₂ O with FA at 1.5 mL/min 4 minutes: Valve Switch (Elute) - Pump for 360 seconds 1:1 ACN/IPA at 1.0 mL/min 8 minutes: Equilibrate - Pump for 240 seconds of H ₂ O with FA at 1.5 mL/min										
HPLC Conditions											
Column	Agilent ZORBAX RRHD Eclipse Plus C18, 3.0 × 50 mm, 1.8 μm										
Mobile phase	A) H ₂ O with 0.2 mM NH ₄ F B) ACN with 0.2 mM NH ₄ F										
Flow rate	0.350 mL/min										
HPLC Elution gradient	<table border="1"> <thead> <tr> <th>Time (min)</th> <th>Mobile phase</th> </tr> </thead> <tbody> <tr> <td>0</td> <td>5 %B</td> </tr> <tr> <td>4</td> <td>5 %B</td> </tr> <tr> <td>11</td> <td>100 %B</td> </tr> <tr> <td>12.5</td> <td>5 %B</td> </tr> </tbody> </table>	Time (min)	Mobile phase	0	5 %B	4	5 %B	11	100 %B	12.5	5 %B
Time (min)	Mobile phase										
0	5 %B										
4	5 %B										
11	100 %B										
12.5	5 %B										
Run time	4 minutes (sample enrichment) + 8.5 minutes (HPLC gradient) + 2.5 minutes (post run) = 16 minutes total										
MS Conditions											
Gas settings											
Drying gas	11 L/min at 250 °C										
Sheath gas	12 L/min at 375 °C										
Nebulizer	45 psi										
Source voltage											
Capillary voltage	4,000 V(+), 3,500 V(-)										
Nozzle voltage	500 V(+), 500 V(-)										

Quantitative linearity over 0.5–200 ng/L concentration range

Figure 4 illustrates highly linear calibration curves ($R^2 > 0.99$) for a quantitation range of 0.5, 1, 2, 5, 10, 20, 50, 100, and 200 ng/L after Online SPE sample enrichment.




Figure 4. Selected calibration curves demonstrating a high degree of sensitivity and linearity for quantitative analysis.

Chromatographic peak reproducibility over various concentration levels

Chromatographic peak overlays in Figure 5 demonstrate the robust retention time stability and peak shape of the dual cartridge Online SPE setup over a large concentration range (0.5–200 ng/L) for a total of 54 injections.




Figure 5. Chromatographic peak overlays and retention times of selected TOxCs.

Limits of detection

Figure 6 illustrates that, of the 31 TOrCs analyzed, 29 compounds have limits of detection that easily fall into the range of 5 ppt or less with only 750 μL of sample.




Figure 6. Distribution of the limits of detection of 31 common TOrCs found in water with 750 μL injection.

Conclusions

- The Agilent 1290 Infinity Flexible Cube enables automated Online SPE sample enrichment for trace organic contaminants.
- Enhanced throughput is achieved using a dual-cartridge setup for rapid sample turnaround.
- The Agilent Ultivo Triple Quadrupole LC/MS enables enhanced trace quantitative analysis with a minimized footprint.
- Technical innovations contained in Ultivo provide optimal sensitivity and robustness, without sacrificing the performance of larger instruments.

References

1. C. G. Daughton, T. A. Terhes. Pharmaceuticals and personal care products in the environment: agents of subtle change. *Environ. Health Perspect.* **107**, 907-938 (1999).
2. T. Anumol, S. Snyder. Rapid analysis of trace organic compounds in water by automated online solid-phase extraction coupled to liquid chromatography-tandem mass spectrometry. *Talanta* **132**, 77-86 (2015).

www.agilent.com
Information, descriptions and specifications in this
publication are subject to change without notice.

© Agilent Technologies, Inc. 2017
Published in USA, June 22, 2017
5991-8155EN
Rev. 1.0



Agilent Technologies