


Introduction

What are Pt-based drugs & cisplatin?

- Used to treat various cancers (e.g., ovarian, lung, colorectal, bladder, etc)
- Prominent Pt-based drugs include:




- But, clinical usefulness decreased due to high incidence of **chemoresistance**

Why study accumulation patterns?


- Resistance mechanisms include:
 - Cellular drug accumulation*
 - Detoxification system (↑ deactivation by thiols)
 - DNA repair
 - Apoptosis regulation
 - Autophagy (self-digestion of drug)

How is ICP-MS relevant?

- Quantification of Pt + other elements
- Agilent ICP-MS modes include:
 - Conventional bulk analysis
⇒ measures Pt mass in overall cell population → limited



- Single-cell (SC) analysis
⇒ measures #cells with Pt + single-cell Pt mass
⇒ reveals heterogeneity + intercellular distributions




Project-at-a-Glance

Motivation & Significance

Aims

- Investigate accumulation patterns of 3 Pt-drugs in ovarian cancer cell models
- By applying **emerging single-cell ICP-MS technology***
- Laterally capitalize on quick bulk analysis of other elements
- To ↑ understanding of (1) Pt-drugs' behavior & (2) downstream/metabolic effects in cells with/without acquired resistance

Conclusion



Great potential for multi-purpose & logistically efficient technique for cancer drug research

Results & Discussion

(a) Validation of SC-ICP-MS analytical workflow

Table 1. Comparison of normalised Pt masses of cisplatin-treated A2780 and A2780cis cells measured by bulk vs. SC-ICP-MS.

Cell line	Cisplatin treatment	Normalised Pt mass (pg / 1000 cells)		Recovery, SC / bulk (%)
		SC-ICP-MS	Bulk ICP-MS	
A2780	10 μM	4.69 ± 1.30	7.71 ± 0.01	60.9
	20 μM	10.06 ± 1.96	16.31 ± 0.02	61.7
A2780cis	10 μM	0.43 ± 0.14	0.70 ± 0.00	62.2
	20 μM	1.60 ± 0.23	2.76 ± 0.01	58.0

- Recovery accounted for 2 analytical considerations:
 - Correction of ionic background in cell samples
⇒ To account for any leaked Pt in suspension
 - Correction of cell count based on CTE
⇒ CTE ~15% (via Gd-DTPA-treated cells)
- High SD for SC-ICP-MS is not a limitation, but indicator of tumour cell heterogeneity¹

(b) Comparison of accumulation patterns of 3 Pt- Drugs (SC)




Figure 1. Pt levels accumulated in A2780 and A2780cis cell lines at varying concentrations of Pt-drugs: (A) shows number of cells with Pt, and (B) shows mass of Pt in single cells.


- Comparing number of cells for 3 Pt-drugs in Figure (A):

⇒ A2780: cisplatin was 2x of carboplatin or oxaliplatin
⇒ More A2780 cells accumulated Pt than resistant A2780cis cells (~10-fold, except at 10μM carboplatin)

- Comparing distribution of Pt levels in Figure (B):

⇒ Pt mass was not significantly lowered in A2780cis cells, compared to sensitive cells.
⇒ Cisplatin & Oxaliplatin: ~ similar SC Pt masses
⇒ Carboplatin: lower single-cell Pt masses

Methods



Details on SC-ICP-MS

- ↓ Interferences
⇒ Use PFA syringe/tubings
- Low, stable 12 μL/min flow rate
⇒ SC nebulization requirements
- ↓ Cell loss (surface adherence)
⇒ Resuspend samples in BSA-containing solution
- ↑ Cell transport efficiency
⇒ Low Ar, high efficiency concentric nebulizer
⇒ Low vol. on-axis spray chamber
- Sample aerosol neb. efficiency
⇒ 1 μm SiO₂ reference
- Cell transport efficiency
⇒ Cell count marker (Gd-DTPA)
- * Image of spray chamber obtained from Glass Expansion (GE)'s website
- Suitable CRC for diff. elements
⇒ Pt & Gd (no gas), Si (H₂)
- Appropriate signal acquisition
⇒ 0.1 ms (for time-resolved single-cell detection)
- Software-enabled data analysis
⇒ MassHunter for 8900 ICP-QQQ (with sNP module)

- ⇒ A2780's low sensitivity to carbo/oxaliplatin; exhibited through both (a) # of cells accumulating Pt, & (b) the levels at which Pt were accumulated
- ⇒ Cross-resistance of A2780cis towards carboplatin & oxaliplatin
- ⇒ Different cellular uptake/efflux of 3 Pt-drugs especially for carboplatin^{2,3}
- ⇒ Uptake unlikely limited to passive diffusion
- ⇒ Ctr1 & OCT2 transporters with different Pt-drug affinities may be implicated^{4,5}

(c) Pt-drugs & diff. endogenous element profiles (bulk)




Figure 2. PCA score plot of the quantified levels of S, P, Mn, Mg, Fe, and S in A2780 and A2780cis cells.

- Clustering of carboplatin & cisplatin-treated samples
⇒ Attributed mainly to diff. in S & Fe levels (based on PCA loadings plot; not shown here)

- ⇒ Implicates S metabolism (possibly thiols & link with Pt detoxification)
⇒ Possibly linked to cisplatin-specific inhibition of Fe regulatory protein⁶

Further SC-ICP-MS work

- Single-nucleus analysis via nuclei extraction ⇒ may study nuclear DNA-platination
- Label metal-tagged antibody biomarkers (imitate flow/mass cytometry)
⇒ e.g., anti-cleaved poly(ADP-ribose) polymerase (PARP) to study level of apoptosis

References

- Liu, T et al. (2021). *Anal Chim Acta* 1177: 338797.
- Garmann, D et al. (2008). *J Control Release* 131(2): 100-106.
- Zisowsky, J et al. (2007). *Biochem Pharmacol* 73(2): 298-307.
- Holzer, AK et al. (2006). *Mol Pharmacol* 70(4): 1390.
- Yonezawa, A et al. (2006). *J Pharmacol Exp Ther* 319(2): 879.
- Miyazawa, M et al. (2019). *Cell Chem Biol* 26(1): 85-97.