

A Novel Strategy to Screen and Profile Steviol Glycosides of Natural Sweeteners in Food Using the ionKey/MS System's Ion Mobility Mass Spectrometry

¹ M. McCullagh, ¹ D. Douce, ¹ R. Rao, and ² S. Gosciny

¹ Waters Corporation, Wilmslow, UK

² Wetenschappelijk Instituut Volksgezondheid-Institut Scientifique de Santé Publique, Brussels, Belgium

APPLICATION BENEFITS

- Leverages the ionKey/MS™ System and the selectivity of ion mobility collision cross sections (CCS) for enhanced ionization/transmission efficiency, providing greater sensitivity for profiling stevioside isomers in complex food commodities.
- Separation of chromatographically co-eluting isomeric species can be achieved using ion mobility.
- Enables unequivocal stevioside isomer identification to authenticate route of origin in food products, i.e. natural, synthesized, semi-synthesized, and manufacturing processes.

WATERS SOLUTIONS

[ionKey/MS System](#)

[iKey™ PCA Separation Device](#)

[ACQUITY UPLC® M-Class System](#)

[Waters® Ion Mobility
Mass Spectrometry Systems](#)

[UNIFI® Scientific Information System](#)

KEY WORDS

Isomer identification, stevioside, isomers, ion mobility, CCS, sweeteners, natural food products, plant profiling, stevia, microflow

INTRODUCTION

Stevia rebaudiana Bertoni is a perennial shrub of the Asteraceae (Compositae) family native to regions of South America. Stevia is of significant economic value due to the high content of natural, dietetically valuable sweeteners in its leaves. It is referred to as “the sweet herb of Paraguay”. Currently, the stevia plant or its extracts are used as sweeteners in North America, South America, Asia, and some European countries. The Joint FAO/WHO Expert Committee on Food Additives (JECFA) has established regulations for steviol glycosides demanding a purity level of at least 95% of the seven chemically defined steviol glycosides.¹

Complete residue analysis of steviol glycosides in food, is a challenge due to nature of the task, i.e. the detection of low steviol glycoside isomer concentrations in complex food commodities, where generic extraction procedures have been used. The challenge is two-fold where low concentrations are present, because they may be used in conjunction with other sweeteners or sugars. In addition, detection at low concentrations is difficult where one or two of the glycosides are at very high concentrations (e.g. rebaudioside A). Isomers of substances may have different chemical properties – they can have different flavor, as well as possible variability in their absorption, distribution, metabolism, elimination, and toxicity. Hence it is necessary to have information on the make up of substances that can contribute to flavor. An initial highly selective, sensitive screening method could be used, where the focus is only aimed towards qualitative, but much more specific detection is required to determine the sweetener's purity, as the purity can impact taste.²⁻⁴

Full scan high resolution mass spectrometry (HRMS) offers high specificity with theoretically no limitation in the number of compounds that can be detected. The continued technology advances of time-of-flight (ToF) mass spectrometry have brought higher sensitivities, resolution, and mass accuracy (typically sub-2 ppm). ToF MS is used in combination with time tolerances, isotopic matching, fragment ions/ratios, and response thresholds to help reduce false positive and false negative detections in screening assays.

Advances in mass spectrometry have vastly improved sensitivity for full spectral analysis, but further enhancements would improve the mass spectral data quality. This is especially important to avoid compromised precursor ion or fragment ion information, and ensure high mass accuracy at low levels. Despite these MS enhancements it can still be a challenge to rapidly and efficiently identify targeted isomeric compounds present in a sample, particularly with large numbers of co-extracted matrix components. Improvements in sensitivity, and the other benefits of using the ionKey/MS System have been described previously,⁵ including enabling sample dilution to reduce matrix suppression, and subsequently increasing the overall analyte signal-to-noise values that can be achieved.

In this application note, we illustrate the selectivity of collision cross section (CCS) measurements used in combination with other recent MS technology enhancements for profiling complex mixtures. A combination of high resolution mass spectrometry and high efficiency ion mobility based measurements and separations is used. Ion mobility is a rapid orthogonal gas separation phase technique that allows another dimension of separation to be obtained within an LC timeframe. Compounds can be differentiated based on size, shape, and charge.

A CCS value is a robust and precise physicochemical property of an ion. It is an important distinguishing characteristic that is related to its chemical structure and three-dimensional conformation.⁶ Using nitrogen-based travelling wave collision cross sections (^{TW}CCSN₂) measurements can increase non targeted screening specificity. Previously generated ^{TW}CCSN₂ measurements have been entered into the UNIFI Scientific Library. This allows the expected and determined ^{TW}CCSN₂ values to be utilized for screening and confirming the presence of steviol glycosides. Here we present a unique approach to screen food products for steviol glycosides using the ionKey/MS System and ion mobility mass spectrometry (IM-MS), which provides unequivocal specificity and sensitivity. The structures of steviol and the steviol glycosides of interest are shown in Figure 1.




Figure 1. Structures of steviol and steviol glycoside used in the study.

EXPERIMENTAL

LC conditions

LC system: ACQUITY UPLC M-Class System
 Mobile phase A: 100% Water (0.1% formic acid)
 Mobile phase B: 100% Acetonitrile (0.1% formic acid)
 Gradient:

Time (min)	Flow rate	%A	%B
0.00	2	99.0	1.0
1.00	2	99.0	1.0
3.00	2	90.0	10.0
5.00	2	70.0	30.0
13.00	2	1.0	99.0
15.00	2	1.0	99.0
15.10	2	99.0	1.0
17.00	2	99.0	1.0

Flow rate: 2.0 $\mu\text{L}/\text{min}$
 Injection volume: 1 μL (full loop)
 Separation device: iKey BEH C_{18} PCA Separation Device, 130 \AA , 1.7 μm , 150 μm x 50 mm ([p/n 186007580](#))
 Separation device temp.: 40 $^{\circ}\text{C}$

MS conditions

MS system: SYNAPT G2-Si
 Ionization mode: ESI-
 Capillary voltage: 2.6 kV
 Sample cone voltage: 20 V
 Lock mass and CCS: Leucine enkephalin, $[\text{M-H}]^{-}$ =554.2620
 Acquisition range: 50 to 1,200 m/z
 Acquisition rate: 10 spectra/sec
 Collision energy ramp: 30 to 70 eV
 Resolution: 20,000 FWHM (Res mode)

Default IMS parameters

IMS T-Wave™
 velocity ramp: Start: 1,000 m/s
 End: 300 m/s
 IMS T-Wave
 pulse height: 40 V
 IMS gas flow: 90 mL
 IMS duty cycle: 10.8 ms

Extraction conditions

Blank chocolate spread sample was purified in two steps. First, fat removal was performed by liquid-liquid extraction. The fat-free extract was subjected to a C_{18} solid phase extraction to remove other matrix components. The final extract was dissolved in acetonitrile. At the end of the sample preparation step, the matrix content in the extract was 0.1 g/mL.

Steviol glycosides

10 steviol glycosides + steviol at ≤ 1 mg/mL in methanol (Table 1) were used to prepare dilutions for matrix fortification post cleanup.

Compound	Solvent standard stock solutions (H_2O)				Expected rt
	Stock conc. mg/mL*	Dilution 1 conc. ng/ μL	Dilution 2 conc. ng/ μL	Dilution 3 conc. ng/ μL	
Rebaudioside E	0.92	9.2	0.92	0.092	6.60
Rebaudioside D	0.80	8.0	0.80	0.080	6.68
Rebaudioside F	0.76	7.6	0.76	0.076	7.31
Rubusoside	0.68	6.8	0.68	0.068	7.58
Rebaudioside B	0.70	7.0	0.70	0.070	7.78
Steviolbioside	0.72	7.2	0.72	0.072	7.82
Stevioside	1.00	10.0	1.00	0.10	9.50
Rebaudioside A	1.00	10.0	1.00	0.10	9.50
Steviol	0.64	6.4	0.64	0.064	9.51
Rebaudioside C	1.00	10.0	1.00	0.10	9.69
Dulcoside A	1.00	10.0	1.00	0.10	9.71

Table 1. Concentration of solvents standards used to prepare stevioside spiked chocolate spread extract samples.

* Indicated initial concentrations are sub 1 mg/mL for 7 out of 11 analytes used.

Spiking protocol

Fortification was performed on blank processed matrix after cleanup to avoid any recovery problems (Table 2). The strategy undertaken in conjunction with CCS profiling, was to determine detection levels in a complex matrix, hence acquisitions were performed with a constant matrix concentration of 10 mg/mL (Table 3).

Starting mass of commodity sample (g/mL)	Volume of matrix	Volume of solvent standards Dilution 1	Volume of H ₂ O	Final concentration
0.1	100 µL	10 µL	790 µL	Steviosides ≤100 pg/µL Matrix=10 mg/mL
Starting mass of commodity sample (g/mL)	Volume of matrix	Volume of solvent standards dilution 2	Volume of H ₂ O	Final concentration
0.1	100 µL	10 µL	790 µL	Steviosides ≤10 pg/µL Matrix=10 mg/mL
Starting mass of commodity sample (g/mL)	Volume of matrix	Volume of solvent standards dilution 3	Volume of H ₂ O	Final concentration
0.1	100 µL	10 µL	790 µL	Steviosides ≤1 pg/µL Matrix=10 mg/mL

Table 2. Method of preparation of diluted end-spiked chocolate spread extracts at 10 mg/mL, and stevioside concentrations at ≤100 pg/µL, 10 pg/µL and 1 pg/µL.

Parameter	Sample loading details
	iKey Separation Device
Injection solvent composition	H ₂ O/Acetonitrile (90:10)
Spiking concentration of matrix (mg/kg)	≤100 µg/kg, 1 mg/kg, 10 mg/kg
Stevioside solution concentration in the final extract (ng/mL)	≤1 ng/mL, 10 ng/mL, 100 ng/mL
Matrix load (µg/mL)	10 mg/mL
Injection volume (µL)	1 µL
Loop size (µL) and injection mode	1 µL full loop
On column steviol + stevioside mass (pg)	Variable see Table 2

Table 3. Sample and matrix loadings for spiked stevioside using the iKey Separation Device.

A Post Column Addition iKey PCA Separation Device (p/n 186007580), shown in Figure 2, incorporates a 1.7 µm, ACQUITY UPLC BEH C₁₈, stationary phase in a 150 µm diameter separation channel. The iKey Separation Device temperature was set to 40 °C, and the eluent from the separation channel flows directly to an integrated ESI emitter. All microfluidic, gas, and electrical connections are automatically engaged when the iKey Separation Device is inserted into the source enclosure and locked into place. The device incorporated an additional channel, enabling post column addition of IPA solvent. The make up solvent was configured to be delivered from channel A of the MS system fluidics for this feasibility study.



Figure 2. Post Column Addition iKey PCA Separation Device.

RESULTS AND DISCUSSION

In this feasibility study, unique sensitivity and selectivity for screening steviol glycosides in complex mixtures has been achieved. Nitrogen based travelling wave collision cross sections ($^{TW}CCSN_2$), accurate mass, fragment ions, and retention time have been obtained to profile the steviol glycosides rebaudioside A to F, rubusoside, steviol, dulcoside A, steviolbioside, and stevioside. CCS measurements were readily obtained for the marker standards at ≤ 100 pg/ μ L, and this information was used to create a scientific library within UNIFI incorporating the expected steviol glycoside $^{TW}CCSN_2$ values. Previous studies have shown the benefits of $^{TW}CCSN_2$ screening, including spectral cleanup, avoidance of false positives, and discovery of pesticide protomers.⁷⁻⁹ The challenge of the assay is further complicated by the requirement to unequivocally identify three pairs of isomers. The steviol glycosides readily fragment and can be prone to insource fragmentation, resulting in isomeric fragments that can result in false positive identifications.¹⁰

Compound	Formula	[M-H] ⁻	Expected Rt	Expected $^{TW}CCSN_2$ (\AA^2)
Rebaudioside E	C44H70O23	965.4230	6.60	289.2
Rebaudioside D	C50H80O28	1127.4758	6.68	321.75
Rebaudioside F	C43H68O22	935.4124	7.32	293.18
Rubusoside	C32H50O13	641.3173	7.56	241.31
Rebaudioside B	C38H60O18	803.3701	7.77	261.19
Steviolbioside	C32H50O13	641.3173	7.81	235.78
Stevioside	C38H60O18	803.3701	7.20	269.64
Rebaudioside A	C44H70O23	965.4230	7.17	299.48
Steviol	C20H30O3	317.2117	9.48	173.38
Rebaudioside C	C44H70O22	949.4280	7.37	299.49
Dulcoside A	C38H60O17	787.3752	7.40	270.75

Table 4. Expected retention times and expected CCS values for steviol and steviosides.

A chocolate spread extract (10 mg/mL) was spiked with the analytes and analyzed using the ionKey/MS System combined with ion mobility, then screened against the stevioside $^{TW}CCSN_2$ library within UNIFI. The $^{TW}CCSN_2$ assignment for glycosides isomer pairs (rubusoside 241.31 \AA^2 /steviolbioside 235.78 \AA^2), (rebaudioside B 261.19 \AA^2 /stevioside 269.64 \AA^2), and (rebaudioside A 298.9 \AA^2 /rebaudioside E 289.2 \AA^2), shown in Table 4, have been determined using stevioside standards. For steviol and the remaining steviosides, CCS measurements were also determined: steviol (173.38 \AA^2), dulcoside A (270.75 \AA^2), rebaudioside F (293.18 \AA^2), rebaudioside C (299.49 \AA^2), and rebaudioside D (321.75 \AA^2).

The expected $^{TW}CCSN_2$ and measured $^{TW}CCSN_2$ values for steviosides spiked into chocolate spread extract are presented in Table 4. The UNIFI Component Summary results obtained (Figure 3) clearly show the benefits of using CCS measurements and the ionKey/MS System with ion mobility. When comparing the expected and measured collision cross sections, the $^{TW}CCSN_2$ measurement errors were typically <0.4%, and the mass measurement error RMS=1.85 ppm obtained for ≤ 100 pg on column loading (actual concentrations have been added (Figure 3 in red text).

Figure 3. UNIFI Component Summary for steviol and profiled steviosides ≤ 100 pg/ μ L (actual on column concentration for each analyte (pg/ μ L) shown in red text) spiked into chocolate spread extract. Three pairs of isomers are highlighted with colored stars, where each pair has been differentiated using CCS.

Component name	Observed m/z	Mass error (ppm)	Observed RT (min)	Observed CCS (\AA^2)	Expected CCS (\AA^2)	Collision cross section error (%)	Response	Adducts
1 Steviol	317.2135	3.93	9.48	172.89	173.38	64pg -0.29	38288	-H
2 Rubusoside	641.3179	0.05	7.54	241.28	241.31	68pg -0.01	8095	-H
3 Steviolbioside	641.3182	0.60	7.80	235.63	235.78	72pg -0.07	36821	-H
4 Dulcoside A	787.3778	2.52	7.39	271.06	270.75	100pg 0.11	20633	-H
5 Stevioside	803.3704	-0.38	7.19	269.48	269.64	100pg -0.06	15680	-H
6 Rebaudioside B	803.3713	0.82	7.75	261.33	261.19	70pg 0.05	24073	-H
7 Rebaudioside F	935.4154	2.65	7.30	292.71	293.18	76pg -0.16	22873	-H
8 Rebaudioside C	949.4299	1.34	7.35	299.69	299.49	100pg 0.07	29190	-H
9 Rebaudioside E	965.4242	0.68	6.60	288.12	289.20	92pg -0.37	47499	-H
10 Rebaudioside A	965.4253	1.87	7.17	298.83	298.90	100pg -0.03	25921	-H
11 Rebaudioside D	1127.4780	1.43	6.67	321.91	321.75	80pg 0.05	21232	-H

In Figure 4, the results obtained at ≤ 1 pg/ μ L are presented in the UNIFI Component Summary, where an overall RMS error of 2.72 ppm was obtained. However the mass error for rubusoside, was -7.46 ppm, where the peak area response was just 63. At this low level the mass measurement error has been affected by other more abundant ions from the matrix background. Screening with a mass measurement tolerance of 10 ppm and CCS% tolerance of 2% has enabled rubusoside to be identified with confidence at trace levels (680 fg) and false negative detections have been avoided.

Component name	Observed m/z	Mass error (ppm)	Expected RT (min)	Observed RT (min)	Observed CCS (\AA^2)	Expected CCS (\AA^2)	Collision cross section error (%)	Response	Adducts
1 Steviol	317.2120	-0.74	9.48	9.47	174.23	173.38	640fg	0.49	679 -H
2 Rubusoside	641.3131	-7.46	7.56	7.55	242.34	241.31	680fg	0.43	63 -H
3 Steviolbioside	641.3173	-0.91	7.81	7.81	235.28	235.78	720fg	-0.21	372 -H
4 Dulcoside A	787.3739	-2.44	7.40	7.39	271.89	270.75	1pg	0.42	146 -H
5 Stevioside	803.3688	-2.32	7.19	7.19	269.47	269.64	1pg	-0.06	163 -H
6 Rebaudioside B	803.3688	-2.31	7.77	7.76	261.84	261.19	700fg	0.25	347 -H
7 Rebaudioside F	935.4106	-2.47	7.32	7.30	292.37	293.18	760fg	-0.28	300 -H
8 Rebaudioside C	949.4277	-0.92	7.37	7.36	299.86	299.49	1pg	0.12	334 -H
9 Rebaudioside A	965.4230	-0.53	7.17	7.17	299.47	298.90	1pg	0.19	370 -H
10 Rebaudioside E	965.4240	0.46	6.61	6.61	288.71	289.20	920fg	-0.17	1099 -H
11 Rebaudioside D	1127.4755	-0.72	6.68	6.67	322.51	321.75	800fg	0.24	229 -H

Figure 4. Component Summary for steviol and profiled steviosides ≤ 1 pg/ μ L (actual on column concentration for each analyte (pg/ μ L shown in red text) spiked into chocolate spread extract. Three pairs of isomers are highlighted with colored stars, where each pair has been differentiated using CCS.

The combined extracted mass chromatogram for steviol and profiled steviosides ≤ 1 pg/ μ L spiked into chocolate spread extract are shown in Figure 5. It can be seen that stevioside and rebaudioside A coelute at 7.18 minutes.




Figure 5. Extracted mass chromatogram for steviol and profiled steviosides ≤ 1 pg/ μ L spiked into chocolate spread extract.

However Figure 6 reveals using ion mobility, that there are two isobaric species at m/z 803.3707 present at retention time 7.18 mins. The retention time aligned multicomponent spectrum at 7.19 mins for rebaudioside A (brown m/z 965.4), and coeluting stevioside (green m/z 803.37) is shown in Figure 7.




Figure 6. UNIFI ion mobility 3D Data Viewer showing two mobility separated compounds at m/z 803.3707, calculated $^{TW}CCSN_2$ values 261.15 Å² and 269.48 Å² at retention time 7.18 mins.




Figure 7. Retention time aligned multicomponent spectrum at 7.19 mins for rebaudioside A (brown m/z 965.4) and coeluting stevioside (green m/z 803.37). Also shown (lower) is the retention time aligned fragmentation spectrum.

Using the “cleaned up” ion mobility product ion spectra, that are retention and drift time aligned, it is possible to obtain the single component precursor and ion mobility product ion spectra for stevioside (Figure 8) and rebaudioside A (Figure 9) resolved from co-eluting components. Use of ion mobility reveals coelution of isobaric species, which would not have been observed without ion mobility separation. As Figure 6 shows, rebaudioside A has produced an insource fragment ion with a CCS of 261.15 \AA^2 , compared to 269.48 \AA^2 for stevioside.




Figure 8. Retention time (7.19 mins) and drift (7.24 ms) aligned spectrum for stevioside (green m/z 803.37). Also retention/drift time aligned ion mobility product ion spectrum.




Figure 9. Single component spectrum at retention time (7.17 mins) and drift (7.83 ms) aligned spectrum for rebaudioside A (green m/z 965.4). Also retention/drift time aligned ion mobility product ion spectrum.

For the first time unique CCS measurements, precursor ion, and corresponding isomer fragmentation spectra for steviol glycosides have been obtained using microfluidic chromatography $^{TW}CCS_N_2$ ion mobility screening. This approach reduces the quantity of two expensive commodities, i.e. high purity standards and solvent, and it has the potential to negate the requirement to repeatedly purchase expensive high purity standards, (€2500 for the standards used in this study) for future screening assays.

CONCLUSIONS

- The ionKey/MS System with ion mobility offers some unique advantages for profiling complex matrices:
 - Sensitivity in combination with high resolution full spectra acquisition.
 - Spectral clean up.
 - Collision cross section to provide unique selectivity and added confidence in identification.
- Ion mobility selectivity has been illustrated, where accurate mass measurement and ^{13}C CCSN₂ measurement have been used to successfully detect and differentiate stevioside isomer residues at trace levels.
- Chromatographically coeluting isomeric species (stevioside and rebaudioside fragment) have been separated.
- The iKey/MS System provides increased sensitivity in order to profile steviosides at low concentrations.
- For the analyst, the ionKey/MS System, brings the benefits of microfluidic chromatography to the required “routine use” platform, in combination with routine ion mobility screening.
- The application of the ionKey/MS System ion mobility offers the potential to reduce analysis costs:
 - Reduced solvent consumption and waste disposal.
 - Reduced requirement to purchase expensive analytical standards.

References

1. J Pól, B Hohnová, T Hyötyläinen. Joint FAO/WHO Expert Committee on Food Additives (JECFA) (2010) *Compendium of Food Additive Specifications, Monograph 5, Steviol glycosides*. <http://www.fao.org/docrep/011/i0345e/i0345e00.htm> (accessed November 2010) 73rd JECFA (2010) published in FAO JECFA Monographs 10 (3).
2. *EFSA Journal* 2011; 9(5):2181.
3. C Well, O Frank, T Hofmann.* Quantitation of Sweet Steviol Glycosides by Means of a HILIC-MS/MS-SIDA Approach. *J Agric Food Chem*. 61: 11312–11320, 2013.
4. R Shah, L S De Jager, T H Begley. Simultaneous determination of steviol and steviol glycosides by liquid chromatography-mass spectrometry. *Food Additives & Contaminants: Part A*. 29,12: 1861–1871, No. 12, December 2012.
5. M McCullagh, S Gosciny, V Hanot, R Tyldesley Worster, and D Douce. Exploring the Benefits and Potential of iKey Microfluidic Chromatography and Time-of-Flight Mass Spectrometry for Pesticide Residue Analysis. Waters Application Note [720005195en](#), 2014.
6. The use of collision cross section measurements (CCS) in food and environmental analysis. Waters Technical Note [720005374en](#), 2015.
7. M McCullagh, G Cleland, V Hanot, J Williams, S Gosciny. Collision Cross Section a New Identification Point for a “Catch All” Non Targeted Screening Approach. Waters Application Note [720005055en](#), 2014.
8. M McCullagh, V Hanot, S Gosciny. The Use of Ion Mobility Spectral Cleanup and Collision Cross Section Values to Increase Confidence and Efficiency in Pesticide Residues Screening Strategies. [Waters Application Note 720005080en](#), 2014.
9. M McCullagh, D Eatough, V Hanot, S Gosciny. Discovery of Pesticide Protomers Using Routine Ion Mobility Screening. Waters Application Note [720005028en](#), 2014.
10. B F Zimmermann. Tandem mass spectrometric fragmentation patterns of known and new steviol glycosides with structure proposals. *Rapid Commun. Mass Spectrom.* 25: 1575–1582, 2011.

Waters

THE SCIENCE OF WHAT'S POSSIBLE.®

Waters, ACQUITY UPLC, SYNAPT, High Definition Mass Spectrometry, HDMS, UNIFI, UPLC, and The Science of What's Possible are registered trademarks of Waters Corporation. T-Wave, ionKey/MS, iKey, and ionKey are trademarks of Waters Corporation. All other trademarks are the property of their respective owners.

©2015–2016 Waters Corporation. Produced in the U.S.A. February 2016 720005421EN AG-PDF

Waters Corporation
 34 Maple Street
 Milford, MA 01757 U.S.A.
 T: 1 508 478 2000
 F: 1 508 872 1990
www.waters.com