6.1050.4x0 Metrosep C 4 - xxx/4.0

6.1050.410	Metrosep C 4 - 100/4.0
6.1050.420	Metrosep C 4 - 150/4.0
6.1050.430	Metrosep C 4 - 250/4.0
6.1050.450	Metrosep C 4 - 50/4.0

DE

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EN

Column material

Silica gel with carboxyl groups, particle diameter 5 µm

Dimensions

6.1050.410	100 x 4.0 mm
------------	--------------

6.1050.420	150 x 4.0 mm
6.1050.430	250 x 4.0 mm
6.1050.450	50 x 4.0 mm

pH range

2 - 7

Temperature range

recommended standard temperature: room temperature

Maximum pressure

25 MPa (250 bar)

Maximum flow rate

2.0 mL/min

Standard flow rate: 0.9 mL/min

Application

Determination of monovalent and divalent cations.

Standard eluent

1.70 mmol/L nitric acid

0.70 mmol/L dipicolinic acid

Preparation

- Rinse the column with eluent for 1 to 2 h.
- To avoid high backpressure we recommend rinsing the column after installation at a low flow rate (0.4 mL/min) for approx. 20 min until the working temperature is reached.

Guard column

Metrosep C 4 Guard/4.0 (6.1050.500)

Metrosep C 4 S-Guard/4.0 (6.1050.510)

Metrosep C 4 S-Guard - 50/4.0 (6.1050.530)

Storage

Store the column in eluent.

Regeneration

NOTE

Ensure that the maximum pressure is never exceeded during regeneration.

If the pressure becomes too high, reduce the flow rate.

Clean the column as follows:

1. Disconnect the column outlet from the detector inlet.
2. Rinse the column with ultrapure water before and after regeneration.

If the column is contaminated with organic components, clean the column as follows:

1. Regenerate the column against flow direction with the following solutions in succession (flow rate 0.90 mL/min, room temperature):
 - 1 h ultrapure water
 - 1 h acetonitrile/water (40/60)
 - 1 h ultrapure water

If the column is contaminated with inorganic components, clean the column as follows:

1. Regenerate the column against flow direction with the following solution (flow rate 0.90 mL/min, room temperature):
 - 1 h 10 mmol/L nitric acid and 4 mmol/L dipicolinic acid

Organic modifiers

- **In the eluent**
0 - 100% acetone and acetonitrile, no alcohol.
- **In the sample matrix**
0 - 100% acetone, acetonitrile and alcohols.

General notes

- Sample solutions must be microfiltered (0.45 µm).
- We recommend producing fresh eluent regularly.

- To avoid high backpressure when changing from/to organic modifiers, adjust the flow rate in small increments from 0.1 mL/min to the measuring conditions within 30 minutes. The maximum pressure of the column must not be exceeded.
- The column must not dry out.
- Regularly replace the eluent aspiration filter and the inline filter because oils, fats, bacteria and other residues may deposit on them.
- The column is not resistant to alcohols.
- To protect the separation column, use the pulsation absorber (6.2620.150) to reduce the injector pressure surges.

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