

Scale-Up of Fast Analytical to Fast Preparative Separations

This article describes a simple, reproducible, and predictable process for scaling up from fast analytical to fast preparative HPLC separations. Chromatographers performing lab scale preparative or larger process scale purifications may use the recommendations discussed. Major considerations for preparative work include ease of scale-up, cost, time, and throughput. The former affects the latter points and is the major focus of this report. Scale-up is simplified by using preparative columns that contain the same smaller diameter packing used in the analytical scout columns. In turn, the time and cost required for method development is reduced. Throughput is increased since loading capacity is higher on the more efficient, smaller diameter packing. Backpressure from smaller diameter packings is not excessive for the shorter, wider columns used in fast preparative separations.

Simple Scale-Up

Smaller, 5µm packings offer distinct advantages over 10µm or larger packings for preparative applications. Successful scale-up from analytical to preparative work is quickly and reliably obtained by using two simple formulas. The resulting chromatography for the preparative column is essentially the same as the analytical column (Figure A). These formulas do not necessarily hold true if the preparative column contains larger diameter packing than the analytical packing:

Loading Capacity: $I_p = I_a \times (D_p/D_a)^2 \times L_p/L_a$

Where, I_p = Injection load of preparative column
 I_a = Injection load of analytical column
 L_p = Length of preparative column
 L_a = Length of analytical column
 D_p = Internal diameter of preparative column
 D_a = Internal diameter of analytical column

Flow Rate: $F_p = F_a \times (D_p/D_a)^2$

Where, F_p = Flow rate of preparative column
 F_a = Flow rate of analytical column
 D_p = Internal diameter of preparative column
 D_a = Internal diameter of analytical column

Cost and Time

Methods developed on 5µm analytical columns can only be easily transferred to preparative scale columns, using the formulas above, when the same 5µm packing is used in the preparative column. Column costs for 10µm preparative columns can be about 40% less expensive than the corresponding 5µm columns. However, scale-up development time and, therefore, expense is often significantly more for 10µm columns. This is due to variances in

(continued on page 4)

Figure A. Same Chromatography Obtained in Fast Analytical to Fast Preparative Column Scale-Up.





Figure B. Efficiency and Backpressure on 5µm vs. 10µm Particles in Preparative Columns.



CONTENTS

Page 1:

- Scale-Up of Fast Analytical to Fast Preparative Separations

Page 2:

- New Products
Preparative HPLC Columns
Injectors

Page 3:

- Literature
Custom Resin & Media
Processing Services

2001 Supelco Chromatography Catalog

- LC Performance Tip
Turn to the Supelco Website First for All Your Application Needs
- NEW SERVICE!
LC Answers by email

Page 4:

- Scale-Up of Fast Analytical to Fast Preparative Separations: (contd.)
- Case Study
Detection Wavelength Adjustment for Preparative HPLC

LITERATURE

Good News From Supelco! Custom Resin and Media Processing Services

Supelco supplies most types of chromatographic resin from the major manufacturers of adsorbents and ion exchangers. We carry an extensive inventory of these media for your convenience. Now you can purchase small quantities of the leading resins without having to stock large amounts of resin due to minimum order levels required by resin producers.

But there's more good news...Supelco will also process your resin into the exact form for your application. You define the exact process to follow in producing your low pressure chromatography media and Supelco takes it from there. Our capabilities include: GMP manufacture, sieving and grinding, blending of different resins, cleaning (aqueous and organic), acid/base wash, ionic conversion, drying, sterilization, cartridge packaging (from 1mL to 4L range) and unit-use packaging.


And there's more good news...Supelco will test your processed resin to assure that it meets your requirements and provide a Certificate of Analysis in most cases. We can perform the following tests on the final resin: USP or EPA extractables, sterility/endotoxigenicity, moisture content, exchange capacity, conductivity and pH.

When you factor in the costs for raw resin and expense of handling and disposing of large volumes of solvents, the value of operator time for processing the media and the convenience of purchasing your resin in the exact form you need it makes sense to consider Supelco's Custom Resin and Media Processing Service.

Here's how it works:

1. We work with you to define the product specification, including testing required, the expected demand, and other relevant information.
2. The Supelco Resin Team reviews the project, conducts a feasibility study and determines if we can produce the resin to the specification.
3. If feasible, we give you an estimate of price and delivery.
4. We refine the process until the product is exactly what you want.
5. We assign a catalog number (at your request) to make it easier to purchase the resin in the future.


S To see how this service can work for you please contact Supelco at 800-652-2572, ext. 9.



T400049

The 2001 SUPELCO Chromatography Catalog

Be sure to check the box on the literature request form to reserve your copy of this valuable resource. Catalogs are scheduled to be mailed in January, 2001.




LC PERFORMANCE TIP

Turn to the Supelco Web Site First for All Your Application Needs

If you are trying to develop a method, need to know how best to analyze for a particular compound or have applications questions, visit the Sigma-Aldrich web site at www.sigma-aldrich.com. At the Supelco section of the site we have hundreds of chromatograms on the many applications we have done, including applications on our Discovery family of columns and SPE materials. You can search the applications section by alphabetical name, technology, compound class or industry type. You will also have access to all of Supelco's other technical information such as Bulletins, Trouble Shooting Guides and Reporter articles. Having difficulty finding the information you need? You can contact our Technical Service Department right from the web site through email. Finally, you can sign up for free newsletters, view the other Sigma-Aldrich brands' web sites and even order directly on line.

S Visit www.sigma-aldrich.com and meet your needs!




NEW SERVICE!

LC Answers by email

If you are trying to develop a method, looking for state-of-the-art column recommendations and separation strategies, or seeking trouble-shooting advice, contact Keith Duff at kduff@sial.com for a fast, accurate response.

Win a free Barlow Multitool, or Executive Discovery Pen if your question is chosen to appear in a future "LC Performance Tip."



Scale-Up of Fast Analytical...

(continued from page 1)

Trademarks and Registered Trademarks:

Barlow - Barwin Limited
Tefzel, Vespel - E.I. DuPont de Nemours and Company
LabPro, Rheodyne - Rheodyne, L.P.
Discovery, SUPELCOSIL, Suplex - Sigma-Aldrich
PEEK - Victrex plc

chromatographic performance between 5 μ m and 10 μ m columns.

Throughput

Peak efficiency is superior on 5 μ m columns compared to 10 μ m columns (Figure B). Preparative columns packed with 5 μ m packings can withstand higher sample loads than 10 μ m columns since the peaks are innately narrower on the 5 μ m columns. This allows more crude sample to be purified in fewer injections. Also, faster flow rates do not broaden peak widths as much on 5 μ m columns as they do on 10 μ m columns. Each individual preparative run can be run faster on 5 μ m than on 10 μ m preparative columns.

Normal Operating Backpressure

Short, wide columns are commonly used in fast preparative separations. A typical fast preparative column size is 5cm x 21.2mm. Backpressure is directly proportional to the length of the column and to the square of the particle radius. The short length of fast preparative columns helps to compensate for the higher backpressure associated with the smaller 5 μ m particle size. This allows fast preparative separations to be run under normal operating backpressures.

For more information, request T495004, and see Keith J. Duff and Richard Ludwig, "Packing and Evaluation of Small Particle Preparative Columns", *American Laboratory*, 32KK-33MM, March (1994).

CASE STUDY 5

Detection Wavelength Adjustment for Preparative HPLC

Detector response is much higher in preparative HPLC than standard analytical scale runs due to the large amount of sample that is injected. To compensate, small volume preparative HPLC detector cells are available from instrument manufacturers, which reduces the response. Still, peaks can (and often do) go off scale, even with a setting at the highest attenuation. The peaks appear as shown in Figure C, where they are "gated" or run off of the top of the chromatogram. Adjustments should be made to keep the response within the detector linear range of about 2AU or less. The best solution is to use a wavelength setting that is somewhat away from the lambda maxima of the peaks of interest, particularly of the largest peak. Most modern multi-wavelength and diode array detectors are capable of scanning the UV absorbance of analytes. Alternatively, a scanning UV instrument is helpful in choosing a good preparative HPLC detection wavelength. Keep in mind that peak heights/areas may change disproportionately to each other when switching the detection wavelength. Some peaks of interest might disappear entirely if a proper wavelength is not chosen. In the example shown, the detector wavelength was switched from 205nm to 230nm (Figure C and Figure D, respectively) to keep the largest peak within a reasonable detection limit.

For more information, request T400183.

Figure C. "Gated" Peak due to Too Strong of a Detector Response @ 205nm






Figure D. Good Detector Response Obtained @ 230nm



SUPELCO • Supelco Park, Bellefonte, PA 16823-0048 • 800-247-6628 • www.sigma-aldrich.com

The
SIGMA-ALDRICH
Family

 **SIGMA**
Biochemicals and
Reagents for Life
Science Research

 **ALDRICH**
Organics and
Inorganics for
Chemical Synthesis

 **Fluka**
Specialty Chemicals
and Analytical
Reagents for Research

 **SUPELCO**
Chromatography
Products for Analysis
and Purification

 **Riedel-deHaen**
Laboratory Chemicals
and Reagents for
Research and Analysis

© 2000 Sigma-Aldrich Co. Printed in USA. Supelco brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.