



Polycyclic aromatic hydrocarbons

Analysis of polycyclic aromatic hydrocarbons to EPA 8100

Application Note

Environmental

Authors

Agilent Technologies, Inc.

Introduction

Gas chromatography with an Agilent CP-Sil 8 CB column separates 20 polycyclic aromatic hydrocarbons according to EPA 8100 in 45 minutes.




Agilent Technologies

Conditions

Technique : GC-capillary
Column : Agilent CP-Sil 8 CB, 0.32 mm x 25 m fused silica
WCOT CP-Sil 8 CB (df = 0.12 µm) (Part no. CP7741)
Temperature : 100 °C (2 min) → 300 °C, 5 °C/min
Carrier Gas : H₂, 45 kPa (0.45 bar, 6.5 psi)
Injector : Split, 100 mL/min
T = 265 °C
Detector : FID
T = 325 °C

Peak identification

1. naphthalene
2. acenaphthylene
3. acenaphthene
4. fluorene
5. phenanthrene
6. anthracene
7. fluoranthene
8. pyrene
9. benz[a]anthracene
10. chrysene
11. benzo[b]fluoranthene
12. benzo[k]fluoranthene
13. benzo[a]pyrene
14. 3-methylcholanthrene
15. indeno[1,2,3-cd]pyrene
16. dibenz[a,h]anthracene
17. benzo[g,h,i]perylene
18. dibenzo[a,e]pyrene
19. dibenzo[a,i]pyrene
20. dibenzo[a,h]pyrene



www.agilent.com/chem

This information is subject to change without notice.

© Agilent Technologies, Inc. 2011

Printed in the USA

31 October, 2011

First published prior to 11 May, 2010

A00801



Agilent Technologies