



# Fatty acid ethylene glycol esters, $C_{16} - C_{36}$

## Analysis of montan wax

### Application Note

Energy & Fuels

#### Authors

Agilent Technologies, Inc.

#### Introduction

The Agilent CP-SimDist UltiMetal column allows high temperature analysis (HT-GC) of the boiling components of refined/modified and raw montan wax. Montan wax is a plant wax derived from lignite by solvent extraction. It is similar to carnauba wax.



**Agilent Technologies**

## Conditions

Technique : GC-wide-bore  
Column : Agilent CP-SimDist UlitiMetal, 0.63 mm x 10 m,  
WCOT CP-SimDist UlitiMetal (df = 0.17  $\mu$ m)  
(Part no. CP7542)  
Temperature : 90 °C (2 min)  $\rightarrow$  440 °C, 10 °C/min; 440 °C (60 min)  
Carrier Gas : N<sub>2</sub>, 10kPa (0.1 bar 1.4 psi)  
Injector : on-column  
Detector : FID, T = 440 °C  
Sample Size : 1.0  $\mu$ L  
Concentration Range : 0.3 g/100 mL  
Solvent Sample : toluene

Courtesy : Dr. L. Matthies and F. Preusser, Völpker  
Montanwachs GmbH, Völpke, Germany


## Peak identification chromatogram A

1. ethylene glycol mono-esters of C<sub>16</sub> - C<sub>36</sub> saturated fatty acids ("montanic acids") hydrocarbons, traces of C<sub>16</sub> - C<sub>36</sub> saturated fatty alcohols ("montanic alcohols")
2. ethylene glycol di-esters of montanic acids.


## Peak identification chromatogram B

1. montanic alcohols, hydrocarbons
2. montanic alcohol esters of montanic acids

**Chromatogram A: Refined and modified montan wax**



**Chromatogram B: Raw montan wax**



[www.agilent.com/chem](http://www.agilent.com/chem)

This information is subject to change without notice.

© Agilent Technologies, Inc. 2011

Printed in the USA

31 October, 2011

First published prior to 11 May, 2010

A01386



**Agilent Technologies**