



# Amino alcohols

## Analysis of ethanol amines

### Application Note

Environmental

#### Authors

Agilent Technologies, Inc.

#### Introduction

Ethanol amines are highly polar compounds and are difficult to elute from a standard capillary column. Only with very thick film is an acceptable peak obtained. However, elution temperatures are high and bleed will cause detection problems at low level.

The Agilent CP-Sil 8 CB for Amines phase elutes the ethanol amines at nanogram level with good peak shape at lower temperatures. The temperature stability of 325 °C allows a quick bake out and a wide range of components to be analyzed.




**Agilent Technologies**

## Conditions

Technique : GC-capillary  
Column : Agilent CP-Sil 8 CB for Amines, 0.32 mm x 30 m  
fused silica WCOT (df = 1.0  $\mu$ m) (Part no. CP7596)  
Temperature : 60 °C (5 min)  $\rightarrow$  220 °C, 6 °C/min  
Carrier Gas : H<sub>2</sub>, 50 kPa (0.5 bar, 7 psi)  
Injector : Split,  
T = 270 °C  
Detector : FID  
T = 300 °C  
Concentration Range : 5-10 ng per component on the column  
Solvent Sample : methanol

## Peak identification

1. MEA (mono-ethanol amine)
2. DEA (di-ethanolamine)
3. TEA (tri-ethanolamine)



[www.agilent.com/chem](http://www.agilent.com/chem)

This information is subject to change without notice.

© Agilent Technologies, Inc. 2011

Printed in the USA

31 October, 2011

First published prior to 11 May, 2010

A01538



**Agilent Technologies**