

Application Bulletin 181/5

Automated potentiometric titration of aluminum and magnesium in the same solution

Branch

General analytical laboratories; pharmaceutical industry; metals, electroplating

Keywords

Titration; automation; aluminum; magnesium; Cu ISE; complexometric titration; back-titration; DCTA; branch 1; branch 4; branch 10; 6.0502.140

Summary

Mixtures of aluminum and magnesium ions can be analyzed automatically by potentiometric titration. After addition of 1,2-diaminocyclohexanetetraacetic acid (DCTA) and complex formation the DCTA excess is back-titrated with copper(II) sulfate. The ion-selective copper electrode is used as indicator electrode. First the aluminum is determined in acidic solution, then the magnesium in alkaline solution.

Instruments

- Titrator with MET mode
- 10 mL buret, 2x
- Sample changer

Electrodes

Cu ISE	6.0502.140
LL ISE Reference	6.0750.100

Reagents

- Copper sulfate pentahydrate, $\text{CuSO}_4 \cdot 5 \text{H}_2\text{O}$
- 1,2-Diaminocyclohexanetetraacetic acid monohydrate, DCTA
- Sodium acetate
- Glacial acetic acid
- Sulfuric acid, concentrated, H_2SO_4
- Ammonia, $w(\text{NH}_3) = 10\%$
- Calcium standard, 6. 2303.070

- $\text{Cu}(\text{NH}_4)_2\text{EDTA}$, Cu-EDTA
- NaOH

Solutions

Titrant	$c(\text{CuSO}_4) = 0.1 \text{ mol/L}$ If possible this solution should be bought from a supplier. Dissolve 24.968 g $\text{CuSO}_4 \cdot 5 \text{H}_2\text{O}$ in approx. 500 mL deion. H_2O , add 0.5 mL conc. H_2SO_4 and make up to 1 L with deion. H_2O .
DCTA solution	$c(\text{DCTA}) = 0.1 \text{ mol/L}$ Dissolve 36.463 g DCTA in 400 mL $c(\text{NaOH}) = 0.5 \text{ mol/L}$ and make up to 1 L with deion. H_2O .
Buffer solution	Dissolve 123 g sodium acetate and 86 mL conc. acetic acid in deion. H_2O and make up to 1 L.
Cu-EDTA	$c(\text{Cu}(\text{NH}_4)_2\text{EDTA}) = 0.1 \text{ mol/L}$ If possible this solution should be bought from a supplier.
NaOH solution	$c(\text{NaOH}) = 1 \text{ mol/L}$

Standard solution

Calcium standard	$c(\text{Ca}^{2+}) = 0.1 \text{ mol/L}$ 6.2303.070
------------------	---

Sample preparation

No sample preparation is required.

Analysis

Titer DCTA

3 to 5 mL calcium standard is dosed into a titration vessel. 50 mL deion. H_2O and 10 mL $w(\text{NH}_3) = 10\%$ and 1 mL Cu-EDTA are added. The solution is then titrated with $c(\text{DCTA}) = 0.1 \text{ mol/L}$ until after the first equivalence point.

Titer CuSO_4

3 to 5 mL $c(\text{DCTA}) = 0.1 \text{ mol/L}$ is dosed into a titration vessel. 50 mL deion. H_2O and 10 mL $w(\text{NH}_3) = 10\%$ are added. The solution is then titrated with $c(\text{CuSO}_4) = 0.1 \text{ mol/L}$ until after the first equivalence point.

Sample

Strongly acidic sample solutions (e.g. from acid digestions) are pre-neutralized to $\text{pH} = 2 \dots 3$ with $c(\text{NaOH}) = 1 \text{ mol/L}$.

Pipet a sample volume containing not more than 12 mg Al^{3+} and 20 mg Mg^{2+} into a titration vessel and dilute with 50 mL deion. H_2O . First the aluminum is determined. Add 6.00 mL $c(\text{DCTA}) = 0.1 \text{ mol/L}$ and 5 mL buffer solution and allow to react for 1 min under stirring. Afterwards the DCTA excess is back-titrated with $c(\text{CuSO}_4) = 0.1 \text{ mol/L}$ until after the first equivalence point.

Now the magnesium can be determined. Add another 6.00 mL $c(\text{DCTA}) = 0.1 \text{ mol/L}$ and 20 mL $w(\text{NH}_3) = 10\%$ to the titrated sample solution and back-titrate the DCTA excess with $c(\text{CuSO}_4) = 0.1 \text{ mol/L}$ until after the first equivalence point.

Parameters

Titer DCTA and CuSO_4

Mode	MET U
Stirring rate	8
Signal drift	50 mV/min
Min. waiting time	5 s
Max. waiting time	26 s
Volume increment	0.1 mL
EP criterion	30 mV
EP recognition	greatest

Al determination

Mode	MET U
Stirring rate	8
Signal drift	50 mV/min
Min. waiting time	5 s
Max. waiting time	26 s
Volume increment	0.1 mL
Stop EP	1
Volume after EP	1 mL
EP criterion	30 mV
EP recognition	greatest

Mg determination

Mode	MET U
Stirring rate	8
Signal drift	50 mV/min
Min. waiting time	5 s
Max. waiting time	26 s
Volume increment	0.1 mL
EP criterion	15 mV
EP recognition	greatest

Calculation

Titer DCTA

$$f_1 = \frac{V_{\text{Std}} \times c_{\text{Std}}}{V_{\text{EP1}} \times c_{\text{DCTA}}}$$

f_1 : Titer of DCTA

V_{Std} : Added volume of standard solution in mL

c_{Std} : Concentration of the standard solution in mol/L

V_{EP1} : Titrant consumption until the first equivalence point in mL

c_{DCTA} : Concentration of DCTA in mol/L

Titer CuSO_4

$$f_2 = \frac{V_{\text{DCTA}} \times f_1 \times c_{\text{DCTA}}}{V_{\text{EP1}} \times c_{\text{CuSO}_4}}$$

f_2 : Titer of CuSO_4

V_{DCTA} : Added volume of DCTA in mL

f_1 : Titer of DCTA

c_{DCTA} : Concentration of DCTA in mol/L

V_{EP1} : Titrant consumption until the first equivalence point in mL

c_{CuSO_4} : Concentration of CuSO_4 in mol/L

Sample

$$\beta_{\text{Al}} = \frac{(V_{\text{DCTA}} \times f_1 \times c_{\text{DCTA}} - V_{\text{EP1.1}} \times f_2 \times c_{\text{CuSO}_4}) \times M_{\text{Al}}}{V_{\text{S}}}$$

$$V_{\text{Ex}} = V_{\text{End}} - V_{\text{EP1.1}}$$

$$\beta_{\text{Mg}} = \frac{(V_{\text{DCTA}} \times f_1 \times c_{\text{DCTA}} - (V_{\text{EP1.2}} + V_{\text{Ex}}) \times f_2 \times c_{\text{CuSO}_4}) \times M_{\text{Mg}}}{V_{\text{S}}}$$

β_{Al} : Aluminum content in g/L

β_{Mg} : Magnesium content in g/L

V_{DCTA} :	Added volume of DCTA in mL
$V_{\text{EP1.1}}$:	Titrant consumption until the first equivalence point of the aluminum titration in mL
$V_{\text{EP1.2}}$:	Titrant consumption until the first equivalence point of the magnesium titration in mL
V_{Ex} :	Excess of titrant added in the aluminum titration in mL
V_{End} :	End volume of the aluminum titration in mL
f_1 :	Titer of DCTA
f_2 :	Titer of CuSO_4
c_{DCTA} :	Concentration of DCTA in mol/L
c_{CuSO_4} :	Concentration of CuSO_4 in mol/L
M_{Al} :	Molecular mass of aluminum; 26.982 g/mol
M_{Mg} :	Molecular mass of magnesium; 24.305 g/mol
V_{S} :	Sample size in mL

Comments

- The surface of the Cu ISE has to be polished from time to time with aluminum oxide powder (6.2802.000 polishing set).
- A waiting time of at least 1 min before the titration of the aluminum is recommended, due to the slow reaction of aluminum with DCTA.

Author

Competence Center Titration
Metrohm International Headquarters

Example determination




Fig. 1: Titration curve of the aluminum determination




Fig. 2: Titration curve of the magnesium determination