

IC Application Note No. C-123

Title: Cations including strontium in brine

Summary: Determination of lithium, sodium, ammonium, potassium, calcium, magnesium and strontium in brine using cation chromatography with direct conductivity detection.

Sample: Brine


Sample Preparation: Dilution 1:500 with 2 mmol/L nitric acid, direct injection after filtration (0.45 µm)

Column: 6.1050.420 Metrosep C 4 – 150

Eluent: 1.7 mmol/L nitric acid
0.7 mmol/L dipicolinic acid

Flow: 0.9 mL/min

Injection Volume: 20 µL



| Results: | Li ⁺ mg/L | Na ⁺ mg/L | NH ₄ ⁺ mg/L | K ⁺ mg/L | Ca ²⁺ mg/L | Mg ²⁺ mg/L | Sr ²⁺ mg/L |
|-------------------|-------------------------|-------------------------|--------------------------------------|------------------------|--------------------------|--------------------------|--------------------------|
| Standard solution | 0.02 | 100 | 0.2* | 20 | 1.0 | 100 | 0.1 |
| Brine (diluted) | 0.01 | 84.1 | 0.1* | 15.7 | 0.55 | 93.7 | n.d. |
| Brine | 5 | 43300 | 54* | 8090 | 287 | 48200 | n.d. |

* ammonium concentration is only an approximation as the separation is not sufficient