# **BHIMADZU**

## Automated PY-GCMS Workflow for the Qualitative and Quantitative Analysis of Microplastics in Environmental Samples <u>Alan Owens, Andy Sandy, Yoshiyuki Okamura; Ruth Marfil-Vega</u>

Shimadzu Scientific Instruments, Columbia, MD

### 1. Introduction

Millions of tons of plastic from food packaging or other household and commercial products are dumped into landfills, or simply thrown out to find their way into various water bodies. The implications of MP pollution in our environment are currently not well understood. In this study, we demonstrated the suitability of a Frontier Lab multi-shot pyrolyzer with a Shimadzu gas chromatography/mass spectrometer (Py-GC/MS) for identification and quantitative measurement of major polymeric components of microplastics. Evolved gas analysis (EGA) analysis followed by flash pyrolysis is used to qualitatively identify these polymers. A quantitative method is also used to determine the number of analytes in the sample.

### 2. Experimental

The configuration consisted of a Shimadzu GCMS, model QP2020 NX, a Frontier Lab multi-shot pyrolyzer, model EGA/PY-3030D, an auto-shot sampler (Shown in figure 1). The PY-GCMS was configured both in EGA and single shot modes. In EGA mode, evolved gases from the polymer appear from heating the sample and were injected into a short inert tube, in the absence of a traditional GCMS column. The temperature program range was from 100 to 700 °C. Leading to a thermogram, a plot of detector response of analytical signal versus furnace temperature, being generated.



Figure 1. Shimadzu GCMS-QP2020 NX and Frontier Multi-Shot Pyrolyzer EGA/PY-3030D.

Gas Chromatography	Nexis GC-2030
Injection port mode	Split mode, 100:1 split ratio
Carrier gas	Helium
Injection port temperature (°C)	300
Column	Single Shot: SH-Rxi-5 MS, 30 m x 0.25 mmlD x 0.25 µm; EGA : EGA
	deactivated tube 2.5 m x 0.15mm
Flow control mode	Linear velocity, 36.1 cm/sec
Oven Temperature	Single Shot :40 °C (4.0 mins.), 20 °C/mins. to 280 °C (7mins); EGA :
	isothermal 300 °C
Mass Spectrometer	QP2020 NX
Interface Temperature (°C)	280
Ion Source Temperature (°C)	230
Detector Voltage (kV)	Relative to Tune -0.03
Threshold	100
Scan Range	m/z 29 to 400
	Scan Speed 1666
Pyrolyzer	EGA/PY-3030D
EGA Furnace Temperature	100 °C, 20 °C/mins. to 700 °C (Total mins 30)
Single Shot Furnace Temp (°C)	600
Interface Temp (°C)	300
Selective Sampler gas	Helium
Selective Sampler Pressure Stability	20 sec
Time	
Auto sampler Purge Time	10 sec

### **3a. Analytical Method**





#### **3b. Sample Prep**

- 1. Polymer Preparation
- Solid polymers are sliced or ground
- Placed into an eco cup with quartz wool
  - 2. EGA-MS Analysis
  - Thermal zone established
  - Optimum PY furnace temp determined for Single Shot
    - 3. Single Shot GCMS
    - Pyrolyzates formed by flash pyrolysis
    - Pyrogram obtained
      - 4. Data Analysis and Comparison
      - All pyrolyzates identified via F-Search library
      - Characteristic pyrolyzate determined via data comparison

Figure 2. Analytical flow of sample preparation and analysis

For qualitative analysis, seven polymers including: acrylonitrile-butadiene-styrene copolymer (ABS); Polycaproamide (Nylon-6); Polyhexamethylene adipamide (Nylon-6,6); Polyethylene terephthalate (PET); Polypropylene (PP); Polystyrene (PS) and Polyvinyl chloride (PVC) were selected. The solid samples were sliced into fine particles (< 200 ug); samples were then placed in a PY eco sample cup for analysis.

Regarding quantitative analysis, twelve solid polymers in CaCO3 were obtained using Frontier MP kit. The homogeneous mixture of the polymers with the CaCO3 diluent made weighing easier as well as prevent reaction within the PY micro furnace. Standards were weighed in the range 0.20 – 4.0 mg. The mass of each polymer in the 4mg of the standard range from  $2.10 - 145.8 \ \mu g$ .

#### 4. Results

#### **Qualitative Analysis**

Figure 2. shows EGA results of the two selected polymers analyzed in this study. EGA was used for identifying a thermal zones



**Figure 3**. Average Mass spectrum and EGA thermogram for ABS and Nylon 6,6.

From the thermal zone, a single shot temperature of 600 °C was determined to use for all analytes. Figure 3 shows pyrogram of 2 selected polymers analyzed by single shot analysis. Table 2 lists the pyrolyzates found in the two selected standards shown in figure 3.



#### **5. Results**

From the thermal zone, a single shot temperature of 600 °C was determined to use for all analytes. Figure 3 shows pyrogram of 2 selected polymers analyzed by single shot analysis. Table 2 lists the pyrolyzates found in ABS shown in figure 4.



For a precision test, eight cups were analyzed in a straight order within the sequence, i.e. injection 1 through 8. ranged from 2.634-13.672 (Table 4).

additional chromatograms, please refer to Shimadzu published application notes.

#### 6. Conclusion

The study demonstrated the satisfactory performance of the Shimadzu GCMS-QP2020NX coupled to a Frontier PY-3030D pyrolizer in quantitation of selected plastics. In this application a fast, robust, and precise workflow was developed for quantitation of twelve polymers. Calibration curve results showed linearity for all compounds, coefficient of determination (r2) ranging from 0.9947-0.9999. Using 8 replicates of standards at 3 mg sample weight, a precision experiment was conducted. Percent RSD results for the targeted compound ranged from 2.634-13.6718.



 Table 3. Summary of method performance