



# Analysis of Pinot Noir Wines by HS-SPME GC/Q-TOF: Correlating Geographical Origin with Volatile Aroma Profiles

---

Philip L. Wylie<sup>1</sup>, Anna K. Hjelmeland<sup>2</sup>, Ron  
Runnebaum<sup>3</sup> & Susan E. Ebeler<sup>3</sup>

- 1) Agilent Technologies, Wilmington, DE 19808  
phil\_wylie@agilent.com
- 2) Agilent Technologies, Santa Clara, CA
- 3) Department of Viticulture and Enology  
University of California, Davis, CA

# Outline

Larger Pinot Noir study objectives

HS-SPME GC/Q-TOF sub-study

HS-SPME method

GC/Q-TOF method

HS-SPME GC/Q-TOF Results

Statistical results




# Larger Study Supported by Jackson Family Wines



- Obtain *Pinot noir* grapes from 15 different vineyards (2015)
  - Same grape clones
  - Same root stock (10 vineyards)
  - Different soils and microclimates
- Deliver grapes to UC Davis winery
- Make four replicate wines from each vineyard
  - Same enological practices used for all wines
- Analyze all wines to see how soil & microclimates affect wine
  - **Volatile Analysis- HS-SPME GC-qTOF (comparing to GC-MS acquired data)**
  - Elemental Analysis- ICP-MS
  - Sensory Analysis- Descriptive Analysis
  - Polyphenolic Analysis- LC-DAD



# American Viticultural Areas in CA and OR Providing Grapes






# Macro weather data available. Need climate in the vineyard where the grapes are grown

One degree day per degree [Fahrenheit](#) over 50 °F.  
Summed from April 1 – Oct. 31

Other considerations:  
*Pinot noir* grown in cooler regions, some with coastal/marine influence

Altitude from near sea level to 2000 ft.




Winkler Heat Index regions in California

Recording weather stations will be installed to get microclimate data




# Solid Phase Microextraction (SPME) Steps



# We have used HS-SPME GC/MS for ultra-trace analysis of haloanisoles in wine


- GC conditions: initial 40C, ramp @ 30C/min to 280, hold for 3 Min, flow rate 1.2 mL/min
- Extraction conditions: SPME headspace, 100 µm PDMS, pre-extraction agitation @ 500 rpm & 40 C for five minutes, extract 10 minutes at 500 rpm & 40 C
- Injection: Splitless, desorb at 280 for 11 min
- Internal standards: d5-TCA, d5-TBA & C13-6 PCA; – for TeCA, C13-6 PCA was used as the internal standard



Solid Phase Micro Extraction

# TCA at 0.1 ng/L 210→195 m/z S/N 13.1


100 ppq is equivalent to 1 second in 320,000 years





# TCA in Customer Complaint Wine – Measured TCA = 2.3 ng/L (ppt)

Figure 1: Chromatogram of customer complaint wine A. Calculated TCA level 2.3 ng/L.



# HS-SPME GC/Q-TOF Method



7200 Accurate Mass High Res. GC/Q-TOF

TOF mode @ 5Hz  
30 m X 0.25 mm X 0.25  $\mu$ m DB-WAXETR  
40°C (5 min); 3°C/min  $\rightarrow$  180°C (min);  
30°C/min  $\rightarrow$  240°C (10 min)





PAL 3 Autosampler for SPME,  
Liquid or HS Injections

100  $\mu$ m, 1 cm Fiber  
Pre-extraction sample equilibration = 5 min @ 30°C  
Headspace extraction = 45 min @ 30°C  
Fiber desorption in MMI inlet = 2 min @ 240°C  
Fiber conditioning = 10 min @ 250°C

# Typical HS-SPME GC/Q-TOF chromatogram of *Pinot noir* wines in the study


Zoom Abundance



# Identification of *Pinot noir* wine volatiles

- Many of the compound were identified by running authentic standards (using RT Locked method on a different GC/MS).
- **For other compounds, we:**
  - Used spectral matching of high resolution accurate mass spectra to NIST 14 unit mass library
  - Calculated Retention Index (RI) values
  - Compared observed RI value to other published values (polar column)
  - Used knowledge of characteristic red wine volatiles with aroma impact
- **Compounds found by spectral searching and RI comparison are *tentatively identified*.**

# Most significant compounds labeled (some identities are tentative)





# Statistical Analysis

Use MassHunter Quant to produce table of Analyte response/ISTD response for 65 identified and tentatively identified compounds


ANOVA used to determine which compounds differed significantly by vineyard

PCA Scores and Loadings plots done using Mass Profiler Professional




# PCA Scores Plot Averaging 3 GC/MS replicates for three wine replicates

Each data point represents 9 measurements





# PCA Scores Plot for 10 vineyards with same combination of grape clone and root stock

First 3 components account for 69 % of the variance



# North CA Coast wines cluster away from Central CA Coast wines






X-Axis Component 1 (29.53%)

Y-Axis Component 2 (23.1%)

Z-Axis Component 3 (16.11%)



# PCA Loadings Plot of significant compounds



# Conclusions

- ❑ Four batches of wine were made from Pinot Noir grapes harvested from 15 different vineyards in California and Oregon (3 batches used for GC/Q-TOF analysis)
- ❑ Three replicate HS-SPME injections made for each wine batch  
15 vineyards X 3 wine batches X 3 replicates = 135 analyses  
9 analyses for each vineyard
- ❑ Volatile profile differs between vineyards in Oregon, CA north coast and CA central coast
- ❑ All vineyards could be separated in PCA
- ❑ This approach could be useful in evaluating regional differences in botanicals

# What still needs to be done

- Obtain climate information for each AVA and each vineyard
  - Placing recording weather station at each vineyard location
- Correlate GC/Q-TOF results with
  - Vineyard microclimate
  - Low resolution GC/MS results
  - Metals analysis
  - Polyphenolic analysis by LC/DAD
  - Sensory Analysis
- Continue investigation over multiple years
- Add more vineyards with same grape clone and same root stock



Thanks to my coauthors: Anna Hjelmeland, Ron Runnebaum & Susan Ebeler

Thanks to Jackson Family Wines for Support and for contributing the grapes

# Thank You!





# UC Davis Department of Viticulture and Enology Teaching & Research Winery

152 research  
Fermenters

Highly  
automated  
e.g. temperature  
control with  
minimal gradients




14 500-gallon  
fermenters




# UC Davis Departments of Viticulture and Enology and Food Science Share the Robert Mondavi Institute – LEED Platinum building Complex



# Solid Phase Microextraction (SPME) Steps




# PCA Loadings Plot of the significant compounds



Not very discriminating peak in the TIC (left). EIC of  $\beta$ -Damascenone (m/z = 121) shows more discrimination

TIC



EIC m/z =121

