

Introduction

Helium (He) Kinetic Energy Discrimination (KED) is a widely used technique to resolve spectra interference of ICP-MS. He KED is, however ineffective against doubly charge ion (M^{++}) overlap. In general, interference by M^{++} is not a big issue though, there are cases where doubly charged ion of rare earth elements (REEs) interfere with As and Se; $^{75}As^+$ suffers from overlap by $^{150}Nd^{++}$ and $^{150}Sm^{++}$, and $^{78}Se^+$ from $^{156}Gd^{++}$ and $^{156}Dy^{++}$. This reports a method to measure trace As and Se in presence of REEs by single quadrupole ICP-MS (ICP-SQ).

Figure 1-1 is a stopping curve of singly charged ion (M^+) and doubly charged ion (M^{++}) of same isotope, neodymium (Nd) in He cell. It indicates M^{++} and M^+ has almost same collision cross section area to He molecule, suggesting that He KED cannot attenuate M^{++} . As shown in Figure 4, BEC of As and Se in presence of REE is worse in He KED than in Nogas, because lighter As^+ and Se^+ lose more Kinetic Energy (KE) than heavier REE^{++} with every collision.

Figure 1-2 is a stopping curve of Nd^+ and Nd^{++} in Hydrogen (H_2) cell, showing that M^{++} drops more rapidly than M^+ . It indicates M^{++} has larger cross section area than M^+ to H_2 molecule. It is a typical ion-molecule interaction via induced dipole due to the large polarizability of H_2 molecule. Based on the finding, H_2 was studied to attenuate REE^{++} interference on As^+ and Se^+ .




Figure 1 Stopping curve of Nd^{++} and Nd^+ at octpoleBias = -18V. Signal of the ions, given in normalized signal was measured as a function of KED. KED is a difference between quadrupole bias voltage and octpole ion-guide bias voltage.

Figure 1-1 (left) stopping curve in He cell and Figure 1-2 (right) stopping curve in H_2 cell.

Experimental 1 : H_2 KED mode 1) typical tune

Figure 4 gives BEC of As and Se in presence of REE in H_2 KED using typical tuning conditions summarized in table1. While BEC of Se was dramatically improved as expected, BEC of As was Not.

This is because As^+ is a unique singly charged ion having large collision cross section area to H_2 . Figure 2 shows signal suppression of elemental ions in H_2 KED mode. Some element such as Ni, Cu and As shows significant suppression, which suggests they have large cross section area to H_2 .




Figure 2 Signal suppression in H_2 KED. Signal ratio (with KED /without KED) is given to each element. Small ratio indicates that the elemental ion has large collision cross section area toward H_2 molecule.

Experimental 2 : H_2 KED mode 2) OctpBias = -6V

There are three tuning parameters which affect cell performance ; Cell gas flow rate, OctpBias and KED. I found that high OctpBias (-6V) attenuated REE^{++} interference on As^+ . As shown in Figure 4, BEC < 0.5ppb was achieved for As and Se in presence of 1ppm REE with the condition. Applying common interference mathematical correction to data further reduced the BEC, below 50ppt.

Measured stopping curve of As^+ and Nd^{++} was given in Figure 3 for two H_2 KED conditions. With OctpBias = -6V, the cross section area of Nd^{++} gets larger (signal drops faster with KED) to H_2 . It improves the BEC of As to satisfactory level at sacrifice of sensitivity. However, with integration time of 3s, MDL < 30ppt (0.1ppb x n = 9) was achieve for both As and Se.




Figure 3 Stopping curve of As^+ and Nd^{++} in two H_2 conditions. Figure 3-1 (left) H_2 1) typical tuning and Figure 3-2 (right) H_2 2) OctpBias = -6V. KE distribution of ions, calculated from the stopping curve are also given.

Result and Summary

BEC of As and Se in presence of 1ppm of REE mix




Figure 4 Summary of the study; BEC of As and Se in presence of 1ppm REE mixture (1ppm of each of Nd, Sm, Gd and Dy).

Table 1 : Tuning conditions
Instrument ; Agilent 7900 x-lens
Plasma condition ; General purpose Plasma

	Nogas	He KED	H_2 KED 1)	H_2 KED 2)
Cell gas	NA	He	H_2	H_2
Cell gas flow rate	ml/min	NA	4.5	3.0
octpBias	V	-5	-18	-6
KED	V	3	3	4
Cell Exit	V	-50	-50	-50
Deflect	V	15	-5	7
Plate Bias	V	-40	-50	-20

Summary and Conclusions

- He KED cannot attenuate the REE^{++} interference on As and Se, since M^{++} has almost same collision cross section area as M^+ to He.
- H_2 KED has more chance to resolve the interference, since M^{++} has larger cross section area than M^+ to H_2 molecule.
- H_2 KED attenuated REE^{++} interference on Se^+ , but didn't REE^{++} interference on As^+ , since As^+ uniquely has large cross section area to H_2 .
- H_2 KED with OctpBias = -6V, attenuated the REE^{++} interference on As^+ , allowing **BEC < 0.5 ppb for As and Se in presence of 1ppm REEs.**
- Use of common mathematical correction in the H_2 KED mode further lowered the BEC to 0.03 ppb for As and Se, allowing accurate determination of 1ppb As and Se in presence of 1ppm REEs.