

RapidFire 360 System – Site Preparation Checklist

Thank you for purchasing an Agilent **solution**. To get you started and to assure a successful and timely installation, please refer to this specification or set of requirements.

Correct site preparation is the key first step in ensuring that your instruments and software systems operate reliably over an extended lifetime. This document is an **information guide AND checklist** prepared for you that outlines the supplies, consumables, space and utility requirements for your equipment for your site.

Customer Responsibilities

Make sure your site meets the following prior specifications before the installation date. For details, see specific sections within this checklist, including:

- The necessary laboratory or bench space available
- The environmental conditions for the lab as well as laboratory gases and plumbing
- The power requirements related to the product (e.g. number & location of electrical outlets)
- The required operating supplies necessary for the product and its installation
- Please consult Other Requirements section below for other product-specific information.
- For more details, please consult the product-specific Site Preparation or Pre-Installation manual.

If Agilent is delivering installation and familiarization services, users of the instrument should be present throughout these services; otherwise, they will miss important operational, maintenance and safety information.

Important Customer Information

1. If you have questions or problems in providing anything described as a Customer Responsibilities above, please contact your local Agilent or partner support/service organization for assistance prior to delivery. In addition, Agilent and/or it's partners reserve the right to reschedule the installation dependent upon the readiness of your laboratory.
2. Should your site not be ready for whatever reasons, please contact Agilent as soon as possible to re-arrange any services that have been purchased.
3. Other optional services such as additional training, operational qualification (OQ) and consultation for user-specific applications may also be provided at the time of installation when ordered with the system, but should be contracted separately.

**RapidFire 360 System –
Site Preparation Checklist**


Dimensions and Weight

Identify the laboratory bench space before your system arrives based on the table below.

Pay special attention to the **total height and total weight requirements for all system components you have ordered and avoid bench space with overhanging shelves**. Also pay special attention to the total weight of the modules you have ordered to ensure your laboratory bench can support this weight.

Special Notes


1. 8 x 8 feet (or 2.5 x 2.5 meters) of space is sufficient for the RapidFire system, mass spectrometer and computer controller.
2. Agilent recommends the laboratory physical arrangement to be such as to minimize the distance between the left wall of the RapidFire platform and the source inlet of the mass spectrometer. This precaution indeed minimizes the length of the fluidic line connecting the two systems.
3. Access to the back of the RapidFire instrument (two feet) will be needed on a regular basis to empty waste containers.

Instrument Description	Weight		Height		Depth		Width	
	Kg	lbs	cm	in	cm	in	cm	in
RapidFire 360 System	320	700	168	66	81	32	155	61




**RapidFire 360 System –
Site Preparation Checklist**

Typical RapidFire set-up – View from above




front




front

RapidFire 360 System – Site Preparation Checklist

Scheme 1



Scheme 2



**RapidFire 360 System –
Site Preparation Checklist**


Environmental Conditions

Operating your instrument within the recommended temperature ranges insures optimum instrument performance and lifetime.

Special Notes

1. Performance can be affected by sources of heat & cold e.g. direct sunlight, heating/cooling from air conditioning outlets, drafts and/or vibrations.
2. The site's ambient temperature conditions must be stable for optimum performance.
3. This system requires a source of vacuum with these specifications: 34 L/min of flow and -70 kPa pressure. The barbed-end fitting to the vacuum source will accommodate the ID 1/4" RapidFire vacuum-rated tubing.

Instrument Description	Operating temp range °C	Operating humidity range (%)
RapidFire system	4 °C to 35 °C	0 % to 95 %, non-condensing




Power Consumption

Special Notes

1. If a mass spectrometer is supplied in conjunction with your RapidFire instrument, be sure to account for those electrical outlets.

Instrument Description	Dedicated Outlet Required (A)	Line Voltage & Frequency (V, Hz)	Maximum Power Consumption (W)
RapidFire system US	12 A	100 / 120 V, 50 - 60 Hz	1440 W
RapidFire system EU	10 A	220 - 240 V, 50 Hz	2400 W
RapidFire system JP	15 A	100 V, 50 / 60 Hz	1500 W

**RapidFire 360 System –
Site Preparation Checklist**

**Software Specifications for
Workstations, Clients, and Servers**
Special Notes

1. Full integration between the RF360 front end and the MS (in order to enable sequence-based batch processing capability in concert with an Agilent or AB Sciex detector) requires one of the following software versions on the MS computer:

Specification Description	Minimum	Comments (if applicable)
MassHunter QQQ	4.1.4114	B08.01
MassHunter QTOF	4.0.4033	6.01 SP1


**Computer Hardware Specifications for
Workstations, Clients, and Servers**
Special Notes

1. The RapidFire computer (supplied by Agilent) will comply with the following requirements:

Specification Description	Minimum	Recommended (if applicable)
Operating system type(s), versions	Windows XP Pro x32 SP3	Windows 7 Pro is fully supported
Processor type and recommended speed	Intel® Xeon® 2.53GHz	
.NET	.NET 3.5 SP1 .NET 4.5	.NET 3.5 SP1 .NET 4.5
Language settings/compatibility		US English
Memory	3Gb	
Mass storage	RAID 1, at least 250Gb	
Specific drivers	1 TCP/IP, 4 COM, 3 USB	

RapidFire 360 System – Site Preparation Checklist



Networking Specifications

Special Notes

1. The RapidFire computer and the mass spectrometer (MS) computer will communicate via a TCP/IP based dedicated network; that is isolated from the MS corporate network.



Required Operating Supplies by Customer

Special Notes

1. For information on Agilent consumables, accessories and laboratory operating supplies, please visit <http://www.chem.agilent.com/en-US/Products/consumables/Pages/default.aspx>
2. Please prepare the set of solvents (HPLC grade) listed in the table below to ensure successful installation.
3. For your application, you may also need these HPLC-grade, filtered solvents: methanol, acetone, triethylamine (TEA), ammonium acetate, and isopropyl alcohol.
4. You may also wish to use 96- and 384-well plates with barcodes on their south side (standard Code 128 ASCII, < 18 characters).

Item Description (including dimensions etc)	Recommended Quantity
Water	2 L
Acetonitrile (MeCN)	2 L
Water + 0.09% formic acid (FA) + 0.01% trifluoroacetic acid (TFA)	2 L
MeCN + 0.09% FA + 0.01% TFA	2 L

Important Customer Web Links

- For additional information about our solutions, please visit our web site at <http://www.chem.agilent.com/en-US/Pages/HomePage.aspx>
- Need to get information on your product?
Literature Library - <http://www.agilent.com/chem/library>
- Need to know more?
Customer Education – <http://www.agilent.com/chem/education>
- Need technical support, FAQs? – <http://www.agilent.com/chem/techsupp>
- Need supplies? – <http://www.agilent.com/chem/supplies>