



Mestrelab Research

Gears Mpublish 1.0

STARTING GUIDE



Document Number

P/N 227 R1



COPYRIGHT

©2021 MESTRELAB RESEARCH S.L.

All rights reserved. No parts of this work may be reproduced in any form or by any means - graphic, electronic, or mechanical, including photocopying, recording, taping, or information storage and retrieval systems - without the written permission of the publisher.

Products that are referred to in this document may be either trademarks and/or registered trademarks of the respective owners. The publisher and the author make no claim to these trademarks.

While every precaution has been taken in the preparation of this document, the publisher and the author assume no responsibility for errors or omissions, or for damages resulting from the use of information contained in this document or from the use of programs and source code that may accompany it. In no event shall the publisher and the author be liable for any loss of profit or any other commercial damage caused or alleged to have been caused directly or indirectly by this document.


The sharing of primary research data is an extremely powerful way to preserve research integrity and replicability and to allow for its reuse by peers. However, preparing research data for publication requires considerable time and effort on the part of the researcher. It can be tedious and error prone.

Mpublish is a plugin for Mnova Gears that solves this problem. It allows easy and automatic preparation of supporting information for publications. Mpublish inputs the primary instrumental data files and generates both human- and machine-readable output files of the research data and metadata.

The current version of Mpublish works with NMR primary datasets (raw data files) and with processed, analyzed, and annotated spectra within the Mnova platform (.mnova files); however, the concept is extensible to many other types of spectroscopic data.

The workflow


Mpublish workflow is both simple and straightforward. So, let's launch Mgears and get started.




Input

In the **Input** section, you should select the analytical data you want to publish. Data files can be saved on a **Disk** directory or can be retrieved from a connected **database (DB)**. In this guide we will be using input files from the **Disk**. *(Please refer to the Mgears manual for further details about other input modes).*


Click on and select your data folder. Both 1D and 2D spectra files can be processed by Mpublish.



You can then choose and configure your experiment detection mode: you can either use **Manual Detection** and therefore enter **Path Masks** for each type of experiment manually, or you can use the **Automatic Detection** option and select the experiments of interests, as per the image below.




If needed, configure the **Advanced Options** or use the **Filtering Options** to refine the detection of your data files of interest.


Top Tip! Before moving to the next step, use the **Automatic inspection** button  to check if the data detected corresponds to your configuration.

Plugins

In the **Plugins** section, select and add the Mpublish plugin.




The plugin will appear on the right-hand side of the window. You should now click on **Mpublish Plugin Settings** to configure the format and information to include in the final report.





In the **Mnova Publish** Settings dialog, first choose the journal template required for your output document. Two formats are available by default with Mpublish 1.0:



- Organic Letters for Supporting Information
- ACS for Supporting Information




If you need additional customizations, please contact us to discuss your requirements.


Type the **Title** of your publication, then click on  to add the **Authors** and corresponding **Organization**. You can indicate the publication's **Main Author** simply by ticking the corresponding checkbox.



The authors added will be listed in the **Authors** table. You can select any such entry and click on  to delete it. You can also click on  to clear the table content.



Finally, you can add any relevant information, abstracts, etc., in the **Information** section.





Name	Organization	Main Author
Robert Boyle	USC	<input checked="" type="checkbox"/>
Mario Novelli	Bruker	<input type="checkbox"/>
Cristina Meyer	Mestrelab	<input type="checkbox"/>

Once all the above has been completed to your satisfaction, click on **OK**.


Output

The last step would be to configure the **Output** section. As with any Mgears plugin, you must choose the location where you would like to save your analysis results.


To save your results on the **Disk**, click on  and select the **Directory** folder.



Here you can also choose to create an Mnova document or a PDF by ticking the corresponding checkbox and indicating the desired output location.




Finally, it is possible to include data in a DB. To do so, you must enter the DB **Server**, **Port**, **User**, and **Password** then click on **Connect**. You will then be able to select the Database.

Now that everything is set, click on  **Run** to start the analysis. Mgears will run the analysis and save the results in the previously defined location.

The output folder

The output folder is named according to the date and time at which your analysis was conducted and contains all the output generated in the current evaluation.



In the Mpublish subfolder you will find:

- your final report saved as an ODT (open document text) file including the pre-defined title and author information, and where all the processed data has been templated according to the journal format selected.

SUPPORTING INFORMATION

Publication title example

Robert Boyle¹, Mario Novelli², Cristina Meyer³

1. USC
2. Bruker
3. [Mestrelab](#)

Additional information

Supporting Information

Figure S1. ¹H NMR spectrum of 3H (400 MHz, DMSO, 29 °C)
 Figure S2. ¹H NMR spectrum of Felodipine.632.fid (600 MHz, DMSO, 298 K)
 Figure S3. ¹H NMR spectrum of m-drive-ds.81.fid (400 MHz, CDCl₃, 297 K)
 Figure S4. ¹H NMR spectrum of Quinine 3H (809 MHz, CDCl₃, 25 °C)
 Figure S5. ¹³C NMR spectrum of Quinine 13C (100 MHz, CDCl₃, 25 °C)
 Figure S6. ²H₂C HSQC-EDITED NMR spectrum of Quinine HSQC (100, 399) MHz, CDCl₃, 25 °C)

Figure S1. ¹H NMR spectrum of 3H (400 MHz, DMSO, 29 °C)

Figure S2. ¹H NMR spectrum of Felodipine.632.fid (600 MHz, DMSO, 298 K)

- a folder containing the primary raw data used to generate these processed and analyzed data. The subfolders and files are renamed and reorganized so they match the titles included in the automatically generated report.

↑ documentation_2021-03-11T15.29.54.zip - archivo ZIP, tamaño descomprimido 4,031,497 bytes

Nombre	Tamaño	Comprimido	Tipo	Modificado	CRC32
..			File folder		
Figure S1	115,977	57,145	File folder		
Figure S2	885,989	550,056	File folder		
Figure S3	1,360,158	933,067	File folder		
Figure S4	96,696	58,710	File folder		
Figure S5	356,302	212,045	File folder		
Figure S6	1,209,371	548,595	File folder		
Structures	7,004	1,038	File folder	3/2/2021 8:39 ...	

The other output subfolders include the Mnova and PDF documents, and Mgears execution files, including the log file and a copy of the mgears settings.

Pro-Tip! The log file is a great way to trouble shoot any issues with the workflow.

For more detailed information about Mgears, please refer to the Mgears Manual or to our Resources section on the web.