

P004b-80A

CLARITY CHROMATOGRAPHY STATION SDK SOLUTIONS

Control Module Development

- Development by DataApex
 - Communication protocols necessary
 - Unit for testing needed
 - Price & time depend on complexity and annual volume
- Development using <u>Clarity</u> <u>Software Development Kit</u> (SDK)

Clarity Software Development Kit

Set of tools for <u>rapid development</u> of control modules <u>by third parties</u>

System Configuration		×			
Setup Control Modules	Number of Instrumen	ts: 2			
Used S/N AS IC C C C C C C C C C C C C	A Instrument 1 AA Instrument 2 Instrument Type GPC Image for Closed Instrument	X Instrument 3 X Instrument 4 Name Instrument 1 Image for Opened Instrument			
		Available Control Modules			
A Colbrick - 4 Instrument 2	> AS	± =	Filter: All	•	×
Add Remove About Setup	C	Name Status S	Vendor Aglent Aglent Aglent Angle Technology Center Dani PG Instruments Shimadzu Labio Dani Netel Bikuia YL Instrument YL Instrument	Comment 4890D, 5890 II, 5890A 6820, 6850, 6850 II, 6890, 68 GC 2010, GC 2014, GC 2014 Gas Chromatograph with UV Analyte 2900A, Chromite 3000	Developed by YL Instrument Under development. Develop
			III		
		Add Cancel			Help

Control Module / Extension Module

- <u>Control Module</u>: is a bridge between Clarity CDS and specific hardware
 - Control modules may be developed by External developers
- <u>Extension Module</u>: enhances Clarity CDS by new data, computations, outputs, etc.
 - Extensions are only developed by DataApex
- Control module may contain additional calculations based on data from Clarity


Clarity SDK

- Provides stable well-documented interface between Clarity software and the control/extension module
- SDK is implemented in C++ and uses Microsoft MFC and ATL libraries (COM)
- Designed to be used with: C++ in Microsoft Visual Studio 2010 (MFC 10) or C++ in Microsoft Visual Studio 2015 (MFC 14)

Clarity SDK Content

- Dynamically linked library (CswSDK10Mfc10.DLL or CswSDK10Mfc14.DLL)
- Set of tools
 (SDK Wizard, DependencyWalker etc.)
- Documentation
 (*.CHM help, examples, SDK Forum)

Clarity SDK Overview



Scenarios of Clarity SDK Usage

- Control module communicates directly with hardware over hardware lines such as COM, USB, TCP/IP
- Control module communicates with hardware through producer's own SDK (COM objects, dlls, etc.)

Clarity SDK

- Control modules developed by SDK will be added to Clarity's portfolio of instrument drivers (Details in DataApex SDK Policy D040)
- Exclusive drivers- a part of OEM version



Clarity SDK

- Control/Extension Modules developed by Clarity SDK will be compatible with future versions of Clarity
- Development Cooperation Agreement required
- SDK is free of charge, all services related to a SDK are charged (price as per valid pricelist)
- SDK support includes an introductory webinar, online email/phone support, Clarity Software loaner and testing of a developed driver(s)
- We offer an optional SDK training in Prague

Localization of control modules

- Clarity is localized to several languages, control module can be too
- Default language for a new control module should be English
- To allow localization, all strings (including states, error messages) should be saved in *.rc files

Released SDK Projects

Clarity already includes dozens of control modules for different instrumentation, both developed by DataApex and by external developers

- GC (Agilent 7890, Shimadzu GC-2010, ...)
- HPLC Systems (Agilent 1200, Shimadzu Prominence, ...)
- HPLC Pumps (Spark SPH1240, Sykam S2100, ...)
- Detectors (Sedere ELSD, ESA Coulochem, ...)
- Autosamplers (Spark Alias, CTC PAL, ...)
- Valves (Upchurch, Rheodyne, Valco, ...)
- PDA Detectors (Knauer S2800, Duratec DDT-3200, ...)
- MS Detectors (Advion Expression CMS, ...)
- Column Ovens (Spark Mistral, ...)
- Fraction Collectors (Buechi 660, ...)
- Special Applications (Sykam Amino Acid Analyzer, ...)

SDK Technical questions:

sdk@dataapex.com

DataApex