



Analysis of Environmental Samples Following US EPA Guidelines Utilizing a New Simultaneous CCD Detector ICP-OES System

Doug Shrader¹, Vincent Calderon², Andrew Ryan³ and Henk Inthorn³
¹ Varian, Inc., 201 Hansen Court, Suite 108, Wood Dale, IL 60191, USA, ² Varian, Inc., 679 Springvale Road, Mulgrave, VIC 3170, Australia
³ Varian B.V., Herculesweg 8, 4338 PL Middleburg, The Netherlands

Introduction

The United States Environmental Protection Agency (US EPA) Contract Laboratory Program (CLP) defines the analytical methods accepted for the isolation, detection and quantitative measurement of target analytes in both water and soil / sediment environmental samples. Data from the Statement of Work (SOW) for Multi Media, Multi Concentration Inorganic Analysis (ILM05.3) is used to define the nature and extent of contamination, and determine appropriate cleanup actions, emergency response actions and enforcement / litigation activities.¹

This poster describes the use of a new simultaneous CCD detector ICP-OES system to carry out the US EPA / CLP compliant analysis of water samples.

¹ ILM05.3, EPA Publication 540-F-04-001, 2004.

US EPA ILM05.3

Determine

- 22 target analytes in water
- NIST Certified Standard Reference Material 1643e Trace Elements in Water
- Melbourne drinking (tap) water

Purpose

1. Provide analytical data of known and documented quality
2. Environmental Regulations

- Extent of contamination
- Determine cleanup actions
- Emergency response and remedial actions
- Enforcement/litigation activities
- Hazardous waste site investigations

Three New ICP-OES Instruments

Varian 710/715-ES

- Entry level system with CCD detector
- A fully featured ICP-OES for laboratories with moderate or small sample loads

Varian 720/725-ES

- High sample throughput
- Ultimate in performance
- Flexibility to upgrade in the future for higher productivity and performance




Varian 730/735-ES

- Highest sample throughput
- Ultimate in performance

The Varian 730-ES simultaneous CCD ICP-OES was utilized for this study. The 730-ES has an axial view plasma. I-MAP technology exactly matches the detector to the echelle optical image and provides continuous wavelength coverage from 167 to 785 nm. The 730-ES features a high efficiency 40 MHz RF generator and low gas consumption. Its Cooled-Cone Interface (CCI) displaces the cooler tail of plasma leading to reduced interferences and increased linear dynamic range. The system includes MFC nebulizer gas, 4 channel pump and the SVS1 Switching Valve System as standard for maximum productivity.


700-ES Series Optical Design

Computer-optimized echelle optical design provides complete wavelength coverage utilizing a single detector and single entrance slit. All wavelengths are captured in one reading. Minimal number of optical components provide high light throughput and excellent signal-to-noise levels. The polychromator is thermostated to 35°C for stability and fast start-up.




The 730-ES features a custom-designed and patented CCD detector with pixels arranged to exactly match the 2-dimensional echelle image (I-MAP technology). Duplex readout circuitry halves readout time. Processing electronics are 40 times faster with a 1 Megahertz pixel processing speed. The detector is Peltier cooled to -35°C for low noise leading to improved detection limits. Processing electronics are off-chip for high QE.

Adaptive Integration Technology



True Simultaneous Measurement

- Simultaneous measurement of all wavelengths, background signal and internal standard irrespective of signal intensity
- AIT automatically adjusts measurement time in real time
- Achieves optimum signal-to-noise ratio
- Prevents over range signals
- 1 ms to 100 s read times provide wide dynamic range
- Ideal for time resolved / specification type applications




Less than 10 minutes of warm-up time is required

MultiCal

Manages the data collected when using alternate wavelengths:

- Combine element data of strong and weak lines for extended working range
- Measure multiple wavelengths of similar sensitivity for confirmation of results
- Take the mean, weighted mean or median result
- Take the minimum result – most likely to be least interfered result
- Automatically performs calculations on data, e.g. ratios, water hardness

Sample Labels	Nu Mean mg/L (mg/L)	Nu 588.811 mg/L (mg/L)	Nu 589.302 mg/L (mg/L)
Blank	0.00	0.00	0.00
Zppsh CS4+6	0.02	0.02	0.02
100pph CS4+6	0.10	0.10	0.10
1ppm CS4+6	1.00	1.00	1.00
10ppm CS4+6	10.00	10.00	10.00
Stipgen Ver Major	20.00	20.00	20.00
200ppm Ver Major	200.0	200.0	200.0
500ppm Ver Major	500.0	500.0	500.0
Blank Mill Powder	395.0	395.0	395.0
Blank Flour 10g	17.20	17.20	17.20
Blank Flour 10g (dup)	17.08	17.08	17.08




4 Port Switching Valve System (SVS1)

Provides shorter sample uptake and more efficient washout. SVS1 is compatible with SP3-3 and Cetac ASX-520 autosamplers and is controlled by the ICP Expert II software. Sample introduction is rinsed while next sample is introduced, leading to reduced carry over, increased sample throughput and lower cost per analysis.





Four Software Triggers

1. Start Run: **Position A (instant rise)**
2. At Start Uptake Delay + user specified time delay: **Position B (sample)**
3. End Solution Measurement: **Position A (instant rise)**
4. End Run: **Position B (autosampler rinse)**



Features and Benefits of SVS1

- 33% improvement in productivity
- 25% reduction in argon use
- 25% reduction in torch use
- Extra 460 samples/day
- Faster analysis times – but longer rinse times!
- Reduces contamination
- Improves accuracy
- Reduces sample loading
- Less torch wear and salt build-up
- Improves long term stability for high TDS sample



Sample Type and Preparation

- Certified Reference Material
- NIST SRM 1643e Trace Elements in Water
- Melbourne drinking (tap) water
- Sample matrix was 1% v/v HNO₃ and 5% v/v HCl

Calibration, QC and Buffer Solutions

Prepared from Inorganic Ventures Inc, custom-grade multi-element solutions

Superfund CLP ICP Kit for ILM05.2

- CLPP-CAL-1, CLPP-CAL-2, CLPP-CAL-3
- CLP-AES-CROL, CLPP-ICS-A, CLPP-CAL-ICS-B4
- CLPP-SPK1, CLPP-SPK-5
- QCP-CICV-1, QCP-CICV-2, QCP-CICV-3

Ionization buffer

- Merck Tracepur™ CsNO₃, 1% w/v final solution

Calibration standards, QC solutions and ionization buffer

- Diluted with >18MW/cm² deionized water (Millipore systems)
- 1% v/v HNO₃ (Merck Ultrapur™)
- 5% v/v HCl (Merck Ultrapur™)

Read Back of Blank + 5 mg/L Analyte

Run	As	Se	Sb	Tl	Pb
Run 1					
Avg. Conc.	5.25	4.81	5.46	5.02	5.39
Run 2					
Avg. Conc.	5.12	4.82	5.01	5.12	5.04
Run 3					
Avg. Conc.	5.13	5.16	4.98	4.80	5.34

Instrument Parameters

Power	1.4 kW	Plasma gas flow	15 L/min
Auxiliary gas flow	1.5 L/min	Spray chamber type	Glass cyclonic
Torch	Standard axial torch	Nebulizer type	SeaSpray
Nebulizer gas flow	0.75 L/min	Pump speed	15 rpm
Sample uptake	2.5 mL	Replicate read time(s)	30
No. of replicates	2	Sample delay time(s)	25
Switching valve delay(s)	22	Stabilization time(s)	10
Rinse time(s)	30	Peristaltic pump type	4 Channel
Background Correction	1 or 2 point off-peak		

Method Detection Limits

Element Wavelength	CRQL (µg/L)	CRQL ILM05.3 (µg/L)	MDL Required (µg/L)	MDL Obtained (µg/L)	Result
Ag 328.068	5	10	5	0.5	Pass
As 227.312	200	200	100	5	Pass
As 188.980	5	10	5	1	Pass
Ba 585.367	20	200	100	0.6	Pass
Be 313.042	1	5	2.5	0.009	Pass
Ca 315.887	5000	5000	2500	1	Pass
Ca 214.439	2	5	2.5	0.09	Pass
Co 228.615	5	50	25	0.4	Pass
Co 267.716	5	10	5	0.2	Pass
Co 324.754	5	25	12.5	0.7	Pass
Cr 199.897	100	100	50	0.3	Pass
Cr 199.897	5000	5000	2500	2	Pass
Mn 285.213	5000	5000	2500	0.4	Pass
Mn 257.610	10	15	7.5	0.06	Pass
Nb 589.592	5000	5000	2500	0.6	Pass
Ni 231.604	20	40	20	0.7	Pass
Pb 220.353	3	10	5	0.8	Pass
Sr 206.834	5	60	30	1	Pass
Se 196.206	5	25	12.5	1	Pass
Tl 190.794	5	25	12.5	1	Pass
V 292.401	10	50	25	0.3	Pass
Zn 206.200	10	60	30	0.5	Pass

LCS Analysis

Element Wavelength	NIST 1643e Measured (mg/L)	NIST 1643e Certified (mg/L)	LCS % Recovery	NIST 1643e Duplicate LCS Measured (mg/L)	Control Limit	% RPD or Difference (mg/L)	Sample Spike Measured (mg/L)	Added QC Spike Conc. (mg/L)	% Spike Recovery
Ag 328.068	0.001963	<CRQL	-	<CRQL	-	-	0.0493	0.0493	100.2%
As 227.312	0.1418	0.151	106.6	0.160	CRQL	0.009	1.03	1.72	104.3%
As 188.980	0.06045	0.0590	97.5	0.0575	20%RPD	2.42%	0.0923	0.0379	97.8%
Ba 585.367	0.0542	0.554	101.9	0.561	CRQL	0.002	2.35	1.72	106.8%
Be 313.042	0.01398	0.0160	100.0	0.0142	CRQL	0.0002	0.0585	0.0430	105.6%
Ca 315.887	32.3	32.0	99.0	32.1	20%RPD	0.50%	-	-	-
Ca 214.439	0.00568	0.00642	97.8	0.00645	CRQL	0.0003	0.0558	0.0473	103.2%
Co 228.615	0.02706	0.0280	103.5	0.0283	CRQL	0.0003	0.484	0.450	106.6%
Co 267.716	0.0204	0.0209	102.4	0.0211	CRQL	0.0003	0.203	0.172	106.8%
Co 324.754	0.02276	0.0229	100.7	0.0242	CRQL	0.0013	0.247	0.215	104.9%
Cr 199.897	0.0981	0.105	106.8	0.104	CRQL	0.001	1.03	0.859	108.5%
Cr 199.897	2.034	2.11	103.7	2.13	CRQL	0.02	-	-	-
Mn 285.213	0.027	0.35	106.6	0.65	CRQL	0.10	-	-	-
Mn 257.610	0.03897	0.0410	105.1	0.0411	CRQL	0.0001	0.503	0.430	108.2%
Nb 589.592	20.74	21.6	104.1	20.9	CRQL	0.7	-	-	-
Ni 231.604	0.06241	0.0629	100.9	0.0639	CRQL	0.001	0.516	0.430	106.1%
Pb 220.353	0.01963	0.0207	105.7	0.0202	CRQL	0.0005	0.0385	0.0189	100.7%
Sr 206.834	0.0583	0.0596	102.2	0.0608	CRQL	0.0012	0.444	0.096	93.8%
Se 196.206	0.01197	<CRQL	-	<CRQL	-	-	0.0561	0.0473	118.5%
Tl 190.794	0.00745	<CRQL	-	<CRQL	-	-	0.0542	0.0473	114.6%
V 292.401	0.02786	0.0389	102.7	0.0388	CRQL	0.0001	0.486	0.430	104.7%
Zn 206.200	0.0785	0.0803	102.3	0.0820	CRQL	0.0017	0.538	0.430	107.7%

Local Water Analysis

Element Wavelength	Melbourne Tap Water Measured (mg/L)	Duplicate Measured (mg/L)	Control Limit	% RPD or Difference (mg/L)	Sample Spike Measured (mg/L)	Added QC Spike Conc. (mg/L)	% Spike Recovery
Ag 328.068	<CRQL	<CRQL	CRQL	-	0.0484	0.0491	98.6
As 227.312	0.0914	0.0948	CRQL	0.001	2.1	1.96	103
As 188.980	<CRQL	<CRQL	CRQL	-	0.0395	0.0361	109
Ba 585.367	0.0180	0.0172	CRQL	0.00078	2.05	1.96	104
Be 313.042	<CRQL	<CRQL	CRQL	-	0.0513	0.0491	104
Ca 315.887	3.64	3.63	CRQL	0.01	-	-	-
Ca 214.439	<CRQL	<CRQL	CRQL	-	0.0488	0.0451	108
Co 228.615	<CRQL	<CRQL	CRQL	-	0.50	0.491	104
Co 267.716	<CRQL	<CRQL	CRQL	-	0.206	0.196	105
Co 324.754	0.162	0.161	20% RPD	0.40%	0.412	0.246	102
Cr 199.897	0.0935	0.0912	CRQL	0.0023	1.10	0.982	102
Cr 199.897	0.598	0.596	CRQL	0.002	-	-	-
Mn 285.213	1.118	1.112	CRQL	0.003	-	-	-
Mn 257.610	0.00617	0.00611	CRQL	0.00006	0.524	0.491	105
Nb 589.592	4.075	4.073	CRQL	0.002	-	-	-
Ni 231.604	<CRQL	<CRQL	CRQL	-	0.516	0.491	105
Pb 220.353	<CRQL	<CRQL	CRQL	-	0.0201	0.0180	112
Sr 206.834	<CRQL	<CRQL	CRQL	-	0.091	0.0901	112
Se 196.206	<CRQL	<CRQL	CRQL	-	0.093	0.091	109
Tl 190.794	<CRQL	<CRQL	CRQL	-	0.0474	0.0451	105
V 292.401	<CRQL	<CRQL	CRQL	-	0.503	0.491	102
Zn 206.200	0.05089	0.0685	CRQL	0.00096	0.530	0.491	107

Interference Check Sample A

Element	CRQL ILM05.3 (µg/L)	ILM05.3 Limit (µg/L)	CSA (µg/L)	Result
Ag 328.068	10	20	-10	Pass
As 188.980	10	20	-2	Pass
Ba 585.367	200	200	0.4	Pass
Be 313.042	5	10	0.1	Pass
Ca 214.439	5	10	0.4	Pass
Co 228.615	50	100	1	Pass
Co 267.716	10	20	0.2	Pass
Co 324.754	25	50	2	Pass
Mn 257.610	15			