

Introduction

The analysis of metals during the beer brewing process is important as certain elements can affect the quality and taste of final beer products. Metals can originate from a range of sources including the brewing water, malt grains, hops, adjuncts, fruits, and spices. They can also be introduced through the brewing and packaging process. Brewers deliberately introduce metals in the form of salts (CaSO₄, MgSO₄, ZnSO₄, CaCl₂) to control pH, adjust taste, improve efficiency, and enhance fermentation performance. Metals that can be detrimental to the overall taste of beer include iron, which can impart a metallic taste. As a result, the concentration of Fe must be kept as low as possible in the finished product.

This study focused on the determination of Ca, K, Mg, Na, Cu, Fe, Mn, and Zn in wort and finished beer samples using the Agilent 5110 VDV ICP-OES.

Experimental

Instrumentation All measurements were performed using an Agilent 5110 VDV ICP-OES. The sample introduction system consisted of a SeaSpray nebulizer, double-pass cyclonic spray chamber, and a 1.8 mm i.d injector torch. Instrument operating parameters are shown in Table 1.

Parameter	Setting
Using a combined run, some elements were measured in radial mode and others in axial.	
Read time (s)	10
Replicates	3
Sample uptake delay (s)	35
Stabilization time (s)	25
Rinse time (s)	35
Pump speed (rpm)	13
Fast pump during uptake and rinse (rpm)	On
RF power (kW)	1.20
Auxiliary flow (L/min)	1.00
Plasma flow (L/min)	12.0
Nebulizer flow (L/min)	0.70
Viewing mode*	Axial
Viewing height (mm)	8
Sample pump tubing	Black/black
Internal standard pump tubing	Orange/green
Waste pump tubing	Blue/blue
Background correction	Fitted

Table 1. Agilent 5110 ICP-OES instrument and method parameters.

Experimental

Standards and sample preparation Wort and beer samples were obtained from the New Belgium Brewing Company (Fort Collins, CO, USA). Sample details are listed in Table 2. All in-process and finished beer samples were sonicated for 20 mins to remove CO₂. Following sonication, 10 mL of beer was diluted in 10 mL of 14% HNO₃. All wort samples were filtered through Whatman paper before preparation. 8 mL of wort was diluted in 32 mL of 8.8% HNO₃ and 3.8% ethanol. All beer and wort samples were stored at 4 °C before analysis.

Sample code	Sample description
W1	Voodoo Ranger IPA wort
W2	Fat Tire wort
W3	Voodoo Ranger 8 Hop Pale Ale wort
W4	Voodoo Ranger Imperial wort
B1	Voodoo Ranger IPA finished beer
B2	Fat Tire finished beer
B3	Voodoo Ranger 8 Hop Pale Ale finished beer
B4	Voodoo Ranger Imperial finished beer

Table 2. New Belgium Brewing Company wort and finished beer samples

Multi-elemental calibration standards were prepared at 0.1 and 0.5 ppm for Cu, Fe, Mn, and Zn; and at 15.0 and 50.0 ppm for Ca, K, Mg, and Na. All standards were prepared in 7% HNO₃ and 3% ethanol. Continuing Calibration Verification (CCV) standards were prepared at 0.02 ppm for minor element checks, 15 ppm for majors, and 500 ppm for K check. A control beer sample spiked with 0.1 ppm Zn was also prepared. Internal standards (ISTD), gallium and yttrium, were prepared at 200 ppm (Ga) and 40 ppm (Y). All calibration standards, Quality Control (QC) checks, and internal standards were matrix matched with 7% HNO₃ + 3% ethanol.

Results and Discussion

Detection limits The Detection Limit (DL) for each element was based on three sigma of seven replicate measurements of the blank solution during an analytical run (Table 3). Sub-ppm (mg/L) level DLs were achieved for all elements. Method DLs (MDLs) were determined by calculating 3 x the SD of 10 replicate beer and wort samples, respectively (Table 3).

Results and Discussion




Figure 1. Calibration curves for Fe 238.204 nm, Cu 327.395 nm, Mg 279.553 nm, and Ca 422.673 nm.

Element & wavelength (nm)	DL (mg/L)	Wort MDL (mg/L)	Beer MDL (mg/L)
Major elements			
Ca 422.673	0.0126	13.35	10.41
K 769.897	0.6539	57.35	66.93
Mg 279.553	0.0091	8.12	15.28
Na 589.592	0.0301	2.07	2.16
Minor elements			
Cu 327.395	0.0005	0.027	0.014
Fe 238.204	0.0008	0.061	0.018
Mn 259.372	0.0001	0.015	0.020
Zn 213.857	0.0009	0.034	0.016

Table 3. Element wavelengths, DLs, and MDLs.

Element and wavelength (nm)	Background correction/interference correction	Calibration range (mg/L)	Correlation coefficient
Ca 422.673	Fitted + ISTD (Y 371.029)	0-50	0.9999
K 769.897	Fitted + ISTD (Ga)	0-50	0.9995
Mg 279.553	Fitted + ISTD (Ga)	0-50	0.9996
Na 589.592	Fitted + ISTD (Ga)	0-50	0.9999
Cu 327.395	Fitted + ISTD (Y 360.074)	0-0.5	0.9999
Fe 238.204	Fitted + ISTD (Y 371.029)	0-0.5	1.0000
Mn 259.372	Fitted + ISTD (Ga)	0-0.5	0.9999
Zn 213.857	Fitted + ISTD (Y 360.074)	0-0.5	1.0000
Ga 417.204	Used as ISTD		
Y 371.029, 360.074	Used as ISTD		

Table 4. Wavelength and working calibration range.

Detectability test To validate the method, a series of QC spike recovery tests were carried out during each analytical run. The tests consisted of a Continuing Calibration Blank (CCB), followed by two CCV solutions (low concentration for Cu, Fe, Mn and Zn, high concentration for Ca, Mg and Na), a K check, and control beer sample spiked with 0.1 ppm Zn. Each QC solution was analyzed six times and averaged, with all recoveries within ±10% of the target values. The QC results are given in Table 5. Analysis of wort and finished beer samples. The four wort and four finished beer samples were analyzed using the developed method. The quantified concentration results are displayed in Table 6. All the results for all elements were within the specification limits set by the manufacturer. Comparing the concentration levels of the elements in wort and beer enables the analyst to monitor the beer brewing process. The data is also useful to characterize the product.

Results and Discussion

	Element and wavelength (nm)							
	Ca 422.673	Cu 327.395	Fe 238.204	K 769.897	Mg 279.553	Mn 259.372	Na 589.592	Zn 213.857
Low conc check 0.02 ppm (mean, n=6), CCV	-	0.0215	0.0210	-	-	0.0208	-	0.0205
% Recovery	-	107	105	-	-	104	-	103
Zinc check 0.1 ppm, control sample (mean, n=6)	-	-	-	-	-	-	-	0.098
% Recovery	-	-	-	-	-	-	-	98
High conc check 15 ppm (mean, n=6), CCV	14.758	-	-	15.761	15.494	-	14.610	-
% Recovery	98	-	-	105	103	-	97	-
Potassium check 500 ppm (mean, n=6)	-	-	-	464.014	-	-	-	-
% Recovery	-	-	-	92	-	-	-	-

Table 5. QC spike recovery results of the low and high CCVs, and control samples.

Sample	Ca 422.673	Cu 327.395	Fe 238.204	K 769.897	Mg 279.553	Mn 259.372	Na 589.592	Zn 213.857
	3	5	7	7	3	2	2	7
	Concentration (mg/L)							
W1	96	0.242	0.045	857	129	0.183	30	0.163
W2	80	0.195	0.149	675	111	0.150	26	0.167
W3	85	0.204	0.054	735	124	0.241	26	0.183
W4	74	0.357	0.055	1125	178	0.263	36	0.262
B1	90	0.099	0.029	719	134	0.287	25	0.016
B2	61	0.047	0.038	464	90	0.132	17	0
B3	101	0.076	0.023	692	106	0.308	12	0
B4	62	0.148	0.032	850	145	0.295	26	0.028

Table 6. Concentration of major and minor elements in wort and finished beer samples.

IntelliQuant semiquantitative results A wort and finished beer sample were analyzed using IntelliQuant during the analytical run. A custom IntelliQuant calibration (1.0, 5.0 and 10.0 mg/L) was created using an Agilent QC standard solution. The semiquantitative results for all elements in the wort and finished beer samples are shown in Table 7. The results reveal the presence of high levels of silicon in the samples. Silicon is present in barley and is introduced at high levels during the beer brewing process. Levels are typically not monitored during the process but can be of interest in terms of dietary intake.

Element and wavelength (nm)	Wort	Finished beer
	IntelliQuant value (mg/L)	IntelliQuant value (mg/L)
Na	9.3	17.5
Mg	24.9	48.3
Si	3.35	17.6
K	136.1	309.5
Ca	20.4	30.6
Mn	0.1	0.1
Sr	-	0.1
B	0.1	-

Table 7. IntelliQuant semiquantitative results of elements in wort and finished beer samples.



The lab at New Belgium. Agilent 5110 VDV (Vertical Dual View) ICP-OES

Conclusions

- Major and minor metals in wort and finished beer samples were measured with good speed and accuracy using the Agilent 5110 VDV ICP-OES.
- The method can be used to provide valuable information to brewers at different stages of production allowing them to optimize the quality of the final product.
- In addition to the quantified method results, high levels of Si were identified in the wort and beer samples using the IntelliQuant (semiquantitative) feature of the ICP Expert software.

Acknowledgements

Agilent Technologies would like to thank the New Belgium Brewing Company for their co-operation and collaboration.

For the full Application note of this poster please visit www.agilent.com/chem

Published September 28, 2017
 Publication number: 5991-8394EN