


# Sample Preparation Method for Analysis of Environmental Pesticides (Simazine, Thiobencarb and Thiuram (Thiram)) in Water Using ISOLUTE® 101 Solid Phase Extraction Columns



**Figure 1.** Analyte Structures.

## Introduction

This application note describes a method for the extraction of pesticides from environmental water samples using ISOLUTE® 101 solid phase extraction columns.

Simazine, a triazine herbicide, Thiobencarb, a thiocarbamate herbicide, and Thiuram (alternatively named Thiram), a dithiocarbamate bactericide are widely used pesticides (Figure 1). As typical pesticides detected in the environment, they are designated as “Target Setting Items for Water Quality Management” in the Japanese Water Supply Law, “Environmental Standards for the Protection of Human Health” in the Environmental Standards for Water Pollution, and “Class III Specified Hazardous Substances” in the Soil Pollution Control Law, and compliance with the standard values set by each of them is required.

ISOLUTE® 101 is a solid phase extraction media consisting of an unmodified form of polystyrene divinylbenzene copolymer. Since its hydrophobic interactions are stronger than with ODS (C18) type media, it is possible to efficiently collect organic compounds with a wide range of hydrophobicities contained in water samples. In addition, as the column can easily be dried with nitrogen gas purge, the overall analysis time can be shortened which increases efficiency.

This method has been prepared in accordance with Japanese official law of the Ministry of the Environment, Notice No. 59, Annex 5 and Annex 6, 1971.

## Analytes

Simazine, CAS: 122-34-9

Thiobencarb, CAS: 28249-77-6

Thiuram (Thiram), CAS: 137-26-8

## Sample Preparation Procedure

### Format

ISOLUTE® 101 200 mg/6 mL, part number 101-0020-C

### Sample Pre-Treatment

To 200 mL of the sample water, add 4 mL of a solution of 10% EDTA·2Na. If thiuram is to be measured, adjust to pH 3.5 with hydrochloric acid.

### Condition

Condition the column with Methanol (5 mL).

### Equilibration

Equilibrate the column with 10% EDTA·2Na solution (4 mL) followed by purified water (15 mL).

### Sample Loading

Load the sample at a flow rate of 100 mL/min using a vacuum pump.

### Interference Elution

Wash the column with purified water (10 mL) at a flow rate of 100 mL/min using a vacuum pump, and dry the column using a flow of nitrogen for 30 minutes.

### Elution

Elute simazine and thiobencarb with acetone (5 mL)

Elute thiuram with acetonitrile (5 mL)

### Evaporation and Reconstitution

Simazine, Thiobencarb: Concentrate the resulting eluate to 1 mL under nitrogen gas in a 35 °C water bath. After evaporation, redissolve the residue in 5 mL of hexane, concentrate again to 1 mL under nitrogen gas in a 35 °C water bath. Add 20 µL of 10 mg/L phenanthrene-D<sub>10</sub> and 20 µL of 5,000 mg/LPEG (1: 1 mixture of PEG 200 and PEG 300) as the internal standard, and perform the GC/MS analysis.

Thiuram: Concentrate the resulting eluate to 1 mL under nitrogen gas in a 35 °C water bath. Transfer the concentrated eluate to vials as it is, and perform HPLC/UV measurements.

## GC/MS Conditions (Simazine, Thiobencarb)

### Instrument

Agilent 7890A/5975C GC-MSD system

### Column

HP-5 MS (30 m x 0.25 mm, membrane pressure 0.25  $\mu$ m, Agilent)

### Carrier

Helium

### Column Flow Rate

1.1 mL/min

### Vaporizing Chamber Temperature

150 °C

### Interface Temperature

280 °C

### Ion Source Temperature

230 °C

### Ionization Method

EI

### Injection Volume

1  $\mu$ L (Splitless)

**Table 1.** GC Column Heating Conditions.

Time (min)	Temperature
0 → 1	50 °C
0 → 11.5	(20 °C/min)
0 → 12.5	280 °C

**Table 2.** MS Conditions.

	For Assay (m/z)	For Detectability (m/z)
Simazine	201	173
Thiobencarb	100	72
Phenanthrene-D <sub>10</sub>	188	160

## HPLC/UV Conditions (Thiuram)

### Equipment

Agilent 1100 HPLC system

### Column

XDB-C18 (75 mm x 3.0 mm, 3.5  $\mu$ m, Agilent)

### Column Temperature

40 °C

### Mobile Phase

A: Water

B: Acetonitrile

### Detection Wavelength

272 nm

### Injection Volume

10  $\mu$ L


**Table 3.** HPLC Mobile Phase Conditions.

Time (min)	%A	%B
0 → 2	80	20
0 → 8	80 → 10	20 → 90
0 → 12	10	90


## Results

The mass chromatogram of Simazine and Thiobencarb mixed standard and the UV chromatogram of Thiuram standard are shown in Figures 3 and 4. Good peaks are obtained for each. In addition, Figure 5 shows calibration curves using the standard solutions. Good linearity is obtained in the corresponding concentration range when

the sample is concentrated 200 times. Reproducibility and recovery at  $n = 5$  were also obtained at the lower limit of the generated calibration curve concentrations (Table 4). Good recoveries and reproducibility have been obtained for all samples. We have also confirmed the validity of this method through spiking and recovery tests on real samples.




**Figure 3.** Mass chromatogram of Simazine and Thiobencarb mixed standard solution (0.2 mg/L).



**Figure 4.** Chromatogram ( $\lambda = 272$  nm) of Thiuram standard solution (0.25 mg/L).

**Table 4.** Reproducibility and recovery rate at the lower measurement limit of Simazine, Thiobencarb and Thiuram ( $n = 5$ ).

	Mean (ug/L)	RSD (%)	Recovery Rate (%)
Simazine	0.3076	2.18	102.5
Thiobencarb	0.2954	2.64	98.5
Thiuram	0.5318	2.68	106.4



**Figure 5.** Calibration curves of Simazine, Thiobencarb (0.05-0.5 mg/L) and Thiuram (0.1-1.0 mg/L).

## Summary

The use of ISOLUTE® 101 allowed the sample to pass through at a superior flow rate (100 mL/min) than with the equivalent competitor product (13 mL/min). In addition, since column drying by centrifugation after column washing, which had been essential until now, has become unnecessary, it was possible to reduce the working time by about 50 minutes per batch while maintaining a recovery rate of approximately 100%.

Further efficiency can be achieved by using automated devices such as TurboVap® LV (Automated Evaporation System) and Biotage® Horizon SmartPrep (Automated Sample Preparation System).

## Ordering Information

Part Number	Description	Quantity
<b>101-0020-C</b>	ISOLUTE® 101 200 mg/6 mL	30
<b>415000</b>	TurboVap® LV	1
<b>41964</b>	TurboVap® LV Multi Rack (48 Positions, 10–20 mm Tubes)	1
<b>SPE2-2616</b>	Biotage® Horizon SmartPrep II Extractor Base Module Configured with 6 mL Plunger Assembly	1

This application note was prepared in cooperation with Sougoukankyoubunseki Co., Ltd.

### EUROPE

Main Office: +46 18 565900  
Toll Free: +800 18 565710  
Fax: +46 18 591922  
Order Tel: +46 18 565710  
Order Fax: +46 18 565705  
order@biotage.com  
Support Tel: +46 18 56 59 11  
Support Fax: + 46 18 56 57 11  
eu-1-pointsupport@biotage.com

### NORTH & LATIN AMERICA

Main Office: +1 704 654 4900  
Toll Free: +1 800 446 4752  
Fax: +1 704 654 4917  
Order Tel: +1 704 654 4900  
Order Fax: +1 434 296 8217  
ordermailbox@biotage.com  
Support Tel: +1 800 446 4752  
Outside US: +1 704 654 4900  
us-1-pointsupport@biotage.com

### JAPAN

Tel: +81 3 5627 3123  
Fax: +81 3 5627 3121  
jp\_order@biotage.com  
jp-1-pointsupport@biotage.com

### CHINA

Tel: +86 21 68162810  
Fax: +86 21 68162829  
cn\_order@biotage.com  
cn-1-pointsupport@biotage.com

### KOREA

Tel: + 82 31 706 8500  
Fax: +82 31 706 8510  
korea\_info@biotage.com  
kr-1-pointsupport@biotage.com

### INDIA

Tel: +91 22 4005 3712  
india@biotage.com

Distributors in other regions are listed on [www.biotage.com](http://www.biotage.com)

#### Part Number: AN911

© 2019 Biotage. All rights reserved. No material may be reproduced or published without the written permission of Biotage. Information in this document is subject to change without notice and does not represent any commitment from Biotage. E&OE. A list of all trademarks owned by Biotage AB is available at [www.biotage.com/legal](http://www.biotage.com/legal). Other product and company names mentioned herein may be trademarks or registered trademarks and/or service marks of their respective owners, and are used only for explanation and to the owners' benefit, without intent to infringe.