


# Extraction of a Drugs of Abuse Panel from Human Urine Using Biotage® Mikro ABN SPE Microelution Plates Prior to UPLC-MS/MS Analysis



**Figure 1.** Example structures by class.

## Introduction

This application note describes the extraction of 49 drugs of abuse from human urine using Biotage® Mikro ABN microelution plates, prior to LC/MS-MS analysis.

The simple sample preparation procedure delivers clean extracts and analyte recoveries mostly greater than 80% with RSDs lower than 10% for most analytes. Linearity of greater than 0.99 is achieved for all analytes from 1–1000 pg/mL.

Mikro plate extraction allows for very low elution volumes and enhanced workflow efficiency.

## Analytes

Amphetamine, Methamphetamine, 3,4-Methylenedioxyamphetamine (MDA), 3,4-Methylenedioxyamphetamine (MDMA), 3,4-Methylenedioxy-N-ethylamphetamine (MDEA), Hydromorphone, Morphine, Benzoyllecgonine (BZE), Oxycodone, Dihydrocodiene, Oxycodone, Mephedrone, Norfentanyl, 7-amino-flunitrazepam, 7-amino-clonazepam, Hydrocodone, Codeine, 6-Monoacetylmorphine (6-MAM), Cocaine, Norketamine, 2-Ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine (EDDP), Zaleplon, Zopiclone, Norbuprenorphine, Ketamine, Nitrazepam, Flunitrazepam, Clonazepam,  $\alpha$ -OH-triazolam, Oxazepam, Estazolam, Temazepam, Zolpidem, Alprazolam, Methadone, Lorazepam, Bromazepam,  $\alpha$ -OH-alprazolam, 2-OH-ethyl-flurazepam, Triazolam, Nordiazepam, Diazepam, Midazolam, Fentanyl, Flurazepam, Buprenorphine, Phencyclidine (PCP), Lysergic acid diethylamide (LSD).

## Internal Standards

Amphetamine-D<sub>5</sub>, Morphine-D<sub>3</sub>, BZE-D<sub>3</sub>, 6-MAM-D<sub>3</sub>, Diazepam-D<sub>5</sub>,

## Sample Preparation Procedure

### Format

Biotage® Mikro ABN Plate, 2 mg, part number 600-0002-LVP

### Sample Pre-Treatment

Spike urine (1 mL) with internal standard solution and allow to equilibrate for 1 hour. Dilute sample with 100 mM NH<sub>4</sub>OAC pH 5 (950  $\mu$ L) and add  $\beta$ -glucuronidase (50  $\mu$ L). Incubate at 60 °C for 2 hours.

Internal standard solution consisted of a 10 pg/  $\mu$ L methanolic solution. 100  $\mu$ L of this was added to 1 mL of urine to give a 1 ng/mL spike concentration.

### Condition (Optional)

Condition wells with methanol (100  $\mu$ L)

### Equilibration (Optional)

Equilibrate wells with 0.1% ammonium hydroxide (aq) (100  $\mu$ L)

### Sample Loading

Load 400  $\mu$ L of the pre-treated urine sample

### Wash 1

Elute interferences with 0.1% ammonium hydroxide (100  $\mu$ L)

### Wash 2

Elute interferences with H<sub>2</sub>O:MeOH (90:10, v/v, 100  $\mu$ L)

### Elution

Elute analytes with DCM:MeOH (90:10, v/v, 30  $\mu$ L)

### Post Elution and Reconstitution

Dry the extract in a stream of air or nitrogen using a Biotage® SPE Dry at 40 °C, 20 to 40 L/min, or a TurboVap® at 40 °C, 1.5 L/min, for approximately 5 minutes.

Reconstitute evaporated samples with H<sub>2</sub>O:MeOH (90/10, v/v) containing 0.1% formic acid (30  $\mu$ L).

Vortex mix and cover plate with a sealing mat prior to injection.

### Processing Conditions

Biotage® Mikro plates were processed using a Biotage® Pressure+ Positive Pressure Manifold.

Settings: Condition, equilibrate, load, wash and elute steps: 7–9 psi (fine control setting).

## UHPLC Conditions

### Instrument

Shimadzu Nexera UHPLC

### Column

Restek Raptor™ Biphenyl 2.7 µm (100 x 2.1 mm)

### Mobile Phase

A: 2 mM ammonium formate (aq) containing 0.1% formic acid

B: 2 mM ammonium formate (MeOH) containing 0.1% formic acid

### Flow Rate

0.4 mL/min

### Injection Volume

5 µL

### Column Temperature

30 °C

**Table 1.** HPLC Gradient.

Time (min.)	%A	%B
0	80	20
2.00	80	20
7.50	40	60
11.25	40	60
12.75	0	100
13.50	0	100
13.51	80	20
15.00	80	20

## MS Conditions

### Instrument

Shimadzu 8060 Triple Quadrupole MS using ES interface

### Nebulizing Gas Flow

3 L/min

### Drying Gas Flow

3 L/min

### Heating Gas Flow

17 L/min

### Interface Temperature

400 °C

### DL Temperature

250 °C

### Heat Block Temperature

300 °C

### CID Gas Flow

270 kPa

**Table 2.** MS conditions for target analytes in positive mode.

Analytes	MRM Transition	Collision Energy
Morphine-D <sub>3</sub>	289.0 > 201.1 (289.0 > 152.1)	-26.0 -50.0
Morphine	286.0 > 152.1 (286.0 > 201.1)	-50.0 -25.0
Oxymorphone	302.00 > 227.1 (302.00 > 198.1)	-30.0 -45.0
Hydromorphone	286.0 > 185.0 (286.0 > 157.0)	-30.0 -40.0
Amphetamine-D <sub>5</sub>	141.0 > 93.0 (141.0 > 124.15)	-15.0 -20.0
Amphetamine	136 > 91.05 (136 > 119.1)	-15.0 -14.0
Methamphetamine	150.0>90.95 (150>119.1)	-20.0 -14.0
MDA	180 > 105 (180 > 77)	-20.0 -40.0
Dihydrocodiene	302 > 119.05 (302 >171)	-35.0 -45.0
Codeine	300.0>215.1 (300.0>165)	-25.0 -40.0
6-MAM-D <sub>3</sub>	331.0 > 165.1 (331.0 > 211.1)	-40.0 -25.0
6-MAM	328.0 > 165.1 (328.0 > 211.1)	-40.0 -25.0
MDMA	194.0 > 163.1 (194.0 > 105.0)	-15.0 -25.0
Oxycodone	316.2 > 241.2	-20.0
Mephedrone	178.00 > 145.05 (178.00 > 144.00)	-20.0 -30.0
Hydrocodone	300.0 > 199.05 (300.0 > 171.1)	-30.0 -40.0
MDEA	208 > 163.05 (208 > 105.05)	-15.0 -25.0
Nor-Ketamine	223.9 > 125 (223.9 > 179.05)	-20.0 -15.0
Nor-Fentanyl	233.0 > 84.05 (233.0 > 56.05)	-20.0 -26.0
BZE-D <sub>3</sub>	293.00 > 171.05 (293.00 > 77.00)	-20.0 -50.0
BZE	289.90 > 168.05 (289.90 > 105.00)	-20.0 -30.0
Ketamine	237.90 > 125.00 (237.90 > 207.05)	-30.0 -14.0
7-Aminoclonazepam	285.90 > 222.10 (285.90 > 121.10)	-25.0 -29.0
Cocaine	304.00 > 182.05 (304.00 > 82.05)	-20.0 -30.0
Zopiclone	388.90 > 245.05 (388.90 > 217.00)	-15.0 -35.0

Table 2. Continued.

Analytes	MRM Transition	Collision Energy
Norbuprenorphine	414.00 > 101.25 (414.00 > 187.20)	-39.0 -38.0
LSD	323.50 > 208.10 (323.50 > 223.25)	-29.0 -23.0
7-Aminoflunitrazepam	283.90 > 135.05 (283.90 > 227.05)	-30.0 -26.0
Zolpidem	308.00 > 235.10 (308.00 > 263.10)	-35.0 -25.0
Buprenorphine	468.10 > 396.25 (468.10 > 414.30)	-40.0 -35.0
Fentanyl	337.00 > 188.10 (337.00 > 105.00)	-20.0 -40.0
Flurazepam	388.00 > 315.00 (388.00 > 288.00)	-20.0 -26.0
PCP	244.00 > 91.05 (244.00 > 159.15)	-35.0 -14.0
Midazolam	325.90 > 249.10 (325.90 > 223.00)	-35.0 -40.0
Bromazepam	315.80 > 182.10 (315.80 > 209.10)	-31.0 -27.0
EDDP	278.00 > 234.00 (278.00 > 234.00)	-30.0 -45.0
Lorazepam	320.80 > 275.00 (320.80 > 229.05)	-22.0 -30.0
Oxazepam	320.80 > 229.05 (286.90 > 104.20)	-23.0 -35.0
Nitrazepam	286.90 > 104.20 (281.90 > 180.10)	-25.0 -35.0
Clonazepam	315.90 > 270.05 (315.90 > 214.05)	-25.0 -38.0
a-OH-Triazolam	358.90 > 331.10 (358.90 > 239.05)	-28.0 -44.0
2-OH-et-flurazepam	332.90 > 211.10 (332.90 > 109.00)	-37.0 -27.0
Methadone	310.50 > 265.10	-16.0
a-OH-Alprazolam	324.90 > 216.10 (324.90 > 205.10)	-39.0 -46.0
Nordiazepam	270.90 > 140.05 (270.90 > 208.10)	-26.0 -28.0


Analytes	MRM Transition	Collision Energy
Zaleplon	305.90 > 236.15 (305.90 > 264.20)	-28.0 -22.0
Flunitrazepam	313.90 > 268.10 (313.90 > 239.10)	-25.0 -35.0
Estazolam	294.90 > 267.05 (294.90 > 205.05)	-20.0 -40.0
Temazepam	300.90 > 255.05 (300.90 > 177.05)	-20.0 -39.0
Triazolam	342.90 > 308.10 (342.90 > 239.05)	-27.0 -41.0
Alprazolam	308.90 > 281.00 (308.90 > 205.05)	-25.0 -40.0
Diazepam-D <sub>5</sub>	289.90 > 193.05 (289.90 > 154.00)	-32.0 -27.0
Diazepam	285.10 > 193.05 (285.10 > 154.00)	-32.0 -27.0




Biotage® PRESSURE+ 96 Positive Pressure Manifold.

## Results

High (mostly > 80%) reproducible (RSD < 10%) recoveries were achieved using the method described in this application note using the Biotage® Mikro plate format.



**Figure 2.** Shows analyte recoveries (1 ng/mL) using the optimized Mikro ABN protocol described in this application note. The time saving Load-Wash-Elute procedure gave similar analyte recoveries. See Discussion section for recovery data generated using the Load-Wash-Elute method.




**Figure 3.** Representative chromatography for application analytes in urine spiked at 1 ng/mL.

Calibration curve performance was investigated from plasma spiked between 1–1000 pg/mL. Good linearity was observed for all analytes typically delivering  $r^2$  values greater than 0.99. Table 3. details linearity performance and associated LOQ for each analyte.

**Table 3.** Analyte calibration curve  $r^2$  and LOQ performance.

Analyte	$r^2$ 200 $\mu$ L Load	LLOQ (pg/mL) 200 $\mu$ L Load
Morphine	0.9990	100
Oxymorphone	0.9991	25
Hydromorphone	0.9996	50
Amphetamine	0.9990	50
Methamphetamine	0.9995	50
Dihydrocodiene	0.9993	1
Codeine	0.9990	25
6-MAM	0.9993	50
MDMA	0.9993	5
Oxycodone	0.9996	25
Mephedrone	0.9993	50
Hydrocodone	0.9990	50
MDEA	0.9993	10
Nor-Ketamine	0.9993	10
Nor-Fentanyl	0.9992	5
BZE	0.9995	50
Ketamine	0.9995	5
7-Aminoclonazepam	0.9998	100
Cocaine	0.9996	25
Norbuprenorphine	0.9995	100
LSD	0.9995	25
7-Aminoflunitrazepam	0.9993	100
Zolpidem	0.9992	5
Buprenorphine	0.9991	10
Fentanyl	0.9990	50
Flurazepam	0.9990	1
PCP	0.9994	10
Midazolam	0.9997	50
Bromazepam	0.9995	100
EDDP	0.9990	1
Lorazepam	0.9996	< 250
Oxazepam	0.9990	250
Nitrazepam	0.9990	500
Clonazepam	0.9994	100
a-OH-Triazolam	0.9992	50
2-OH-et-flurazepam	0.9995	< 100
Methadone	0.9990	250
a-OH-Alprazolam	0.9997	250
Nordiazepam	0.9990	50
Zaleplon	0.9991	5
Flunitrazepam	0.9991	50
Estazolam	0.9993	5
Temazepam	0.9995	< 250
Triazolam	0.9993	1
Alprazolam	0.9994	25
Diazepam	0.9996	5



**Figure 4.** Calibration curves for Fentanyl (a), 6-MAM (b), Cocaine (c) and Methamphetamine (d) using the Biotage® Mikro plate to extract human urine.


## Discussion and Conclusion

The Biotage® Mikro solid phase extraction microelution plate provided robust extraction of a large drugs of abuse panel from hydrolyzed urine samples. High, reproducible recoveries were achieved, with an overall processing time of approximately 22 minutes (including the evaporation step). An evaporation step was required in this application, as the elution solvent (DCM/MeOH) which gave the highest analyte recoveries was not compatible with the reversed phase analytical UPLC system. However, due to the low elution volume, evaporation time for 96 samples was ~ 5 minutes.

Compared to the equivalent procedure using a 10 mg EVOLUTE® EXPRESS ABN plate, a reduction in total processing time (~22 mins vs ~33 mins) was possible. In addition, higher, more reproducible recoveries for some analyte types (such as benzodiazepines) were achieved.

### Comparison with EVOLUTE® Load-Wash-Elute Methodology

Elimination of the optional conditioning and equilibration steps in this method was investigated. A reduced overall extraction time (and solvent usage) was possible, and analyte recoveries were not adversely impacted (see figure 5).



**Figure 5.** Analyte recoveries achieved using the Load-Wash-Elute methodology (where conditioning and equilibration steps are eliminated). This approach reduced the overall processing time by approximately 5 minutes (~17 mins for L-W-E method vs ~22 mins for full method).



## Chemicals and Reagents

- » Methanol (LC-MS grade), Ultra-Pure Methanol (Gradient MS), and dichloromethane (99.8%) were purchased from Honeywell Research Chemicals (Bucharest, Romania).
- » All analyte standards, deuterated internal standards, ammonium acetate, ammonium formate, formic acid and ammonium hydroxide (27–30%) were purchased from Sigma- Aldrich Company Ltd. (Gillingham, UK).
- » Water used was 18.2 MOhm-cm, drawn daily from a Direct-Q5 water purifier.
- » Mobile phase A (2 mM ammonium formate (aq), 0.1% formic acid) was prepared by adding 0.126 mg of ammonium formate to 1 L purified water with 1 mL formic acid.
- » Mobile phase B (2 mM ammonium formate (methanolic), 0.1% formic acid) was prepared by adding 0.126 mg of ammonium formate to 1 L ultra-pure MeOH with 1 mL formic acid.
- » Internal standards (100 pg/μL) were prepared from a 10 ng/μL stock solution by adding 10 μL of each to 950 μL of MeOH. 10 μL of this solution was then added to each calibration solution.
- » Hydrolysis buffer solution (100 mM ammonium acetate) was made by adding 0.3854 mg of ammonium acetate to 50 mL of water (18.2 MOhm-cm).
- » Hydrolysis enzyme β-Glucuronidase from Helix pomatia, Type HP-2; ≥100,000 units/mL purchased from Sigma- Aldrich Company Ltd. (Gillingham, UK).
- » Equilibration and wash 1 solvent (0.1% ammonium hydroxide) was made by adding 100 μL of ammonium hydroxide to 99.9 mL of water (18.2 MOhm-cm).
- » Wash 2 solvent (H<sub>2</sub>O:MeOH (90:10, v/v)) was made up by measuring out 90 mL of water (18.2 MOhm-cm) and 10 mL of methanol and adding both to a bottle.
- » Elution solvent ((DCM:MeOH (90:10, v/v)) was made up by measuring out 90 mL of DCM (18.2 MOhm-cm) and 10 mL of methanol and adding both to a bottle.
- » Reconstitution solvent was made by measuring out 90 mL of purified water (18.2 MOhm-cm) and 10 mL of MeOH and adding them to the same bottle with 100 μL formic acid.

## Additional Information

- » All data shown in this application note was generated using human urine donated by healthy human volunteers.

## Ordering Information

Part Number	Description	Quantity
<b>600-0002-LVP</b>	Biotage® Mikro ABN Plate, 2 mg	1
<b>PPM-96</b>	Biotage® PRESSURE+ 96 Positive Pressure Manifold	1
<b>SD-9600-DHS</b>	Biotage® SPE Dry Sample Concentrator System	1
<b>121-5202</b>	Collection Plate, 1 mL Square	50
<b>121-5204</b>	Pierceable Sealing Mat	50

### EUROPE

Main Office: +46 18 565900  
 Toll Free: +800 18 565710  
 Fax: +46 18 591922  
 Order Tel: +46 18 565710  
 Order Fax: +46 18 565705  
 order@biotage.com  
 Support Tel: +46 18 56 59 11  
 Support Fax: +46 18 56 57 11  
 eu-1-pointsupport@biotage.com

### NORTH & LATIN AMERICA

Main Office: +1 704 654 4900  
 Toll Free: +1 800 446 4752  
 Fax: +1 704 654 4917  
 Order Tel: +1 704 654 4900  
 Order Fax: +1 434 296 8217  
 ordermailbox@biotage.com  
 Support Tel: +1 800 446 4752  
 Outside US: +1 704 654 4900  
 us-1-pointsupport@biotage.com

### JAPAN

Tel: +81 3 5627 3123  
 Fax: +81 3 5627 3121  
 jp\_order@biotage.com  
 jp-1-pointsupport@biotage.com

### CHINA

Tel: +86 21 68162810  
 Fax: +86 21 68162829  
 cn\_order@biotage.com  
 cn-1-pointsupport@biotage.com

### KOREA

Tel: +82 31 706 8500  
 Fax: +82 31 706 8510  
 korea\_info@biotage.com  
 kr-1-pointsupport@biotage.com

### INDIA

Tel: +91 22 4005 3712  
 india@biotage.com

Distributors in other regions are listed on [www.biotage.com](http://www.biotage.com)

### Literature Number: AN941

© 2020 Biotage. All rights reserved. No material may be reproduced or published without the written permission of Biotage. Information in this document is subject to change without notice and does not represent any commitment from Biotage. E&OE. A list of all trademarks owned by Biotage AB is available at [www.biotage.com/legal](http://www.biotage.com/legal). Other product and company names mentioned herein may be trademarks or registered trademarks and/or service marks of their respective owners, and are used only for explanation and to the owners' benefit, without intent to infringe.