

Application News

Analysis of SAF Using Comprehensive Two-Dimensional Gas Chromatography (GC×GC)

Kei Usui, Yu Nagao

User Benefits

- ◆ By comparing chromatograms obtained as two-dimensional images, it can be utilized for quality control inspection of SAF.
- ◆ Changing the column combination customizes separation characteristics to the sample composition.
- ◆ Utilizing the FASST (Fast Automated Scan/SIM Type) allows for the simultaneous analysis of the composition of SAF and the quantitative analysis of antioxidants.

Introduction

SAF (Sustainable Aviation Fuel) produced from waste cooking oil, used paper, and sugarcane has gained attention in recent years from the perspective of achieving a carbon neutral society. The physical properties of SAF, such as combustibility, vary depending on its components. Therefore, component analysis of SAF is crucial; however, since SAF is a mixture of various hydrocarbons, complete separation of each component is extremely challenging. Comprehensive two-dimensional gas chromatography (GC×GC), which utilizes the separation characteristics of two different types of columns, is an effective method for analyzing SAF with such a complex composition. Additionally, the GCMS-QP2050 enables ultra-fast scanning (30,000 u/sec), which improves the separation of adjacent peaks and allows for high-precision qualitative analysis([Application News No.01-00788-EN](#)).

In this application, two types of SAF (A and B) and one intermediate product from the SAF production process (C) were measured using a combination of GC×GC with a non-polar column as the first column and a medium-polarity column as the second column. The two-dimensional images allow for a visual representation of the composition distribution, making quality control inspection. Moreover, high separation of paraffins and naphthenes using a combination of a medium-polar column as the first column and a non-polar column as the second column, as well as the quantification of antioxidants in SAF, which are regulated by concentration in ASTM D2425 will be introduced.

System Configuration

The system configuration is shown in Tables 1. In this analysis, the ZX1 modulator from Zoex Corp. (Houston, TX, USA) was used. This modulator features thermal modulation using liquid nitrogen, allowing for the attainment of sharp peaks, including those of low-boiling compounds.

Table 1 System Configuration

GC Model	: Nexis™ GC-2030 / AOC-30i
MS Model	: GCMS-QP2050 (TMP exhaust: 255 L/sec)
Injection Port	: SPL
Modulator	: ZX1 thermal modulator (ZOEX corp.)
Nitrogen generator	: MT-24F (System Instruments Co., Ltd.)
Software	: Labsolutions™ GCMS GC Image (GC Image, LLC)
Library	: NIST Mass Spectral Library

Comparison of SAF 2D images (Low-polarity → Medium-polarity)

An analysis of three different samples (A, B, and C) were conducted. The analysis conditions are shown in Table 2, and the obtained two-dimensional images are presented in Figure 1. A low-polar column, DB-5MS, was used as the first column, and a medium-polar column, BPX-50, was used as the second column. In this combination of columns, the horizontal axis represents boiling point, while the vertical axis corresponds to polarity. Sample A contains fewer aromatic compounds and naphthenes compared to sample B. Sample C is composed solely of low-boiling paraffins and naphthenes. Thus, by comparing two-dimensional images, intuitive comprehension of compounds distribution can be achieved, making it applicable for quality control of SAF.

Table 2 Analysis Conditions of SAF (Low polarity → Medium polarity)

GC	
1 st Column	: DB-5MS (30 m × 0.25 mm I.D., df= 0.25 μm)
2 nd Column	: BPX-50 (P/N : 054740) (2.5 m × 0.10 mm I.D., df= 0.10 μm) *cut the column into 2.5 m
Injection Temperature	: 280 °C
Flow Control Mode	: Pressure (He)
Inlet pressure Program	: 150 kPa → 0.6 kPa/min → 200 kPa → 1.2 kPa/min → 250 kPa
Purge Flow	: 3 mL/min
Injection Volume	: 0.2 μL
Split ratio	: 100
Column Oven Temp. Program	: 40 °C → 1.0 °C/min → 120 °C → 2.0 °C/min → 220 °C → 4.0 °C/min → 280 °C (15 min)
MS	
Ion Source Temperature	: 230 °C
Interface Temperature	: 280 °C
Measurement Mode	: Scan
Scan Range (m/z)	: 35-500
Event Time	: 0.02 s (= 50 Hz, 30,000 u/sec)
Modulator	
Modulation period	: 6 s
Hot pulse duration	: 350 ms (320 °C)
Cold Gas Flow	: 15 L/min (N ₂)




Fig. 1 TIC two-dimensional plot of the analysis of SAF samples (Low-polarity → Medium-polarity)

■ Medium-polarity → Non-polarity

In SAF that contains a high amount of paraffins and naphthenes, using a medium-polar column (SH-17) as the first column and a non-polar column (SH-1) as the second column is effective for improving their separation. In this case, the distribution differs from that of the previously mentioned combination of low-polar and medium-polar columns, with lower polar compounds being distributed towards the upper position. An example of a two-dimensional image is shown in Figure 2, where it can be seen that paraffins and naphthenes are separated over a wide range of retention times.

Additionally, clear separation was observed in the region of naphthenes and aromatic compounds. For paraffins, both straight-chain and branched types were separated vertically by the second column, making them easy to distinguish. Furthermore, the distribution of monoaromatic and diaromatic compounds was also distinguishable. In the analysis of diesel fuel, monocyclic, bicyclic, and tricyclic aromatic compounds can be distinguished from each other (Figure 3). The analysis conditions for SAF and diesel fuel are shown in Tables 3 and 4, respectively.




Fig. 2 TIC two-dimensional plot of the analysis of SAF sample (Medium-polarity → Non-polarity)




Fig. 3 TIC two-dimensional plot of the analysis of diesel sample (Medium-polarity → Non-polarity)

Table 3 Analysis Conditions of SAF (Low polarity→Medium polarity)

GC	
1 st Column	: SH-17 (P/N : 221-75907-30) (30 m × 0.25 mm I.D., df= 0.25 μm)
2 nd Column	: SH-1 (P/N : 221-75723-30) (2.5 m × 0.32 mm I.D., df= 0.25 μm) *cut the column into 2.5 m
Injection Temperature	: 250 °C
Flow Control Mode	: Pressure (He)
Inlet pressure Program	: 20 kPa → 1.0 kPa/min → 30 kPa → 2.0 kPa/min → 180 kPa
Purge Flow	: 3 mL/min
Injection Volume	: 0.1 μL
Split ratio	: 200
Column Oven Temp. Program	: 40 °C → 0.5 °C/min → 45 °C → 2.0 °C/min → 195 °C
MS	
Ion Source Temperature	: 230 °C
Interface Temperature	: 280 °C
Measurement Mode	: Scan
Scan Range (<i>m/z</i>)	: 35-500
Event Time	: 0.02 s (= 50 Hz, 30,000 u/sec)
Modulator	
Modulation period	: 11 s
Hot pulse duration	: 350 ms (350 °C)
Cold Gas Flow	: 15 L/min (N ₂)

Table 4 Analysis Conditions of diesel (Low polarity→Medium polarity)

GC	
1 st Column	: SH-17 (P/N : 221-75907-30) (30 m × 0.25 mm I.D., df= 0.25 μm)
2 nd Column	: SH-1 (P/N : 221-75723-30) (2.5 m × 0.32 mm I.D., df= 0.25 μm) *cut the column into 2.5 m
Injection Temperature	: 250 °C
Flow Control Mode	: Pressure (He)
Inlet pressure Program	: 20 kPa → 2.0 kPa/min → 220 kPa
Purge Flow	: 3 mL/min
Injection Volume	: 0.2 μL
Split ratio	: 100
Column Oven Temp. Program	: 40 °C → 2.0 °C/min → 240 °C
MS	
Ion Source Temperature	: 230 °C
Interface Temperature	: 280 °C
Measurement Mode	: Scan
Scan Range (<i>m/z</i>)	: 35-500
Event Time	: 0.02 s (= 50 Hz, 30,000 u/sec)
Modulator	
Modulation period	: 10 s
Hot pulse duration	: 350 ms (350 °C)
Cold Gas Flow	: 15 L/min (N ₂)

Quantitative Analysis of Antioxidant

Additives such as antioxidants and static dissipators are contained in SAF. While the concentrations of various compounds in typical SAF are often in the percentage range, antioxidants like BHT (Butylhydroxytoluene) are regulated to around several tens of ppm in standards such as ASTM D2425. Therefore, the concentration of BHT is significantly different from that of SAF components, making it undetectable in Scan analysis of SAF components. However, by conducting measurements using FASST (Fast Automated Scan/SIM Type), it is possible to perform SIM analysis simultaneously during the retention time of BHT in addition to Scan analysis, allowing for the detection of BHT (Figure 4). When BHT was added to sample A at a concentration of 20 mg/L and a calibration curve was created using the external standard method, it was quantified at 19 mg/L, achieving a recovery rate of 95%.




Fig. 4 FASST analysis of sample A
(a)Scan, (b)SIM, (c)Calibration curve

Analysis Condition: See Table2,
Event time : Scan (0.02s), SIM (*m/z*:205, 220, 0.01s)

Conclusion

Three SAF samples were measured using a combination of low-polar to medium polar columns and compared the two-dimensional images. The two-dimensional images allow for intuitive comprehension of the distribution of each component, making it easy to identify differences in composition and conducting quality control inspection. Additionally, in the measurement of SAF using a combination of medium-polar to non-polar columns, the separation of paraffins and naphthenes was improved. Furthermore, by utilizing the FASST, composition analysis of SAF and quantitative analysis of antioxidants were simultaneously performed, which are regulated by concentration in ASTM D2425.

<Acknowledgement>

The authors wish to thank Dr. Naoto Fujino of Biomaterial in Tokyo Co.,Ltd., Dr. Haruki Tani and Mr. Shugo Onitsuka of Environment Energy Co., Ltd. for their generous cooperation, including providing the samples used in these measurements.

<Related Applications>

1. Comprehensive Characterization of Diesel Fuel on GC×GC Utilizing Impressive High-Speed Scan Technology of GCMS-QP2050, [Application News No.01-00788-EN](#)

Nexis and LabSolutions are trademarks of Shimadzu Corporation or its affiliated companies in Japan and/or other countries.



SHIMADZU

Shimadzu Corporation

www.shimadzu.com/an/

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. See <https://www.shimadzu.com/about/trademarks/index.html> for details.

Third party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

01-00862-EN

First Edition: Aug. 2025

➤ Please fill out the survey

Related Products

Some products may be updated to newer models.



➤ GCMS-QP2050

Gas Chromatograph Mass
Spectrometer

Related Solutions

Hydrocarbon
➤ Processing Industry
(Petrochemical, Ch

➤ Petroleum refinery

➤ Chemicals

➤ Clean Energy

➤ Price Inquiry

➤ Product Inquiry

➤ Technical Service /
Support Inquiry

➤ Other Inquiry