

CoreFocus
Report
No.362

GC

FID


AOC

SH Series

SH-PolarWax

Analysis of Residual Ethylene Oxide in Medical Devices


Keywords: EOG sterilization




1. Ethylene oxide (EO)
2. Ethylene chlorohydrin (ECH)
3. Ethylene glycol (EG)

Model	: Nexis™ GC-2030 / AOC-20i Plus
Detector	: FID-2030 Flame Ionization Detector
Column	: SH-PolarWax (30 m x 0.53 mm I.D. df = 1.00 μm), P/N: 221-75979-30
Column Temp.	: 60 °C (3 min) - 20 °C/min - 200 °C (10 min) Total 20 min
Injection Temp.	: 250 °C
Injection Mode	: Split
Split Ratio	: 1:3
Carrier Gas	: N ₂
Carrier Gas Controller	: Constant Linear Velocity
Linear Velocity	: 40 cm/s
Detector Temp.	: 250 °C
Detector Gas	: H ₂ 32 mL/min, Air 200 mL/min
Makeup Gas	: N ₂ 24 mL/min
Injection Volume	: 0.5 μL
Syringe	: Elastic Syringe, AOC, P/N: 221-49548*


* Using an elastic syringe for AOC (P/N: 221-49548) equipped with a plunger made of titanium enables stable sample introduction.



Chromatogram of
EO



Chromatogram of ECHx



Chromatogram of
EG

Source : Application News 01-00139 ([JP](#), [ENG](#))

CoreFocus and Nexis are trademarks of Shimadzu Corporation or its affiliated companies in Japan and/or other countries.

Shimadzu Corporation

www.shimadzu.com/an/

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu.

Company names, products/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol "TM" or "®".

Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.