


Application Note




Analysis of Acrylamide in Processed Foods in Japan


Dr. Hans-Ulrich Baier,
Shimadzu Deutschland GmbH.

Since announcement of presence of acrylamide in processed foods by University of Stockholm and National Food Administration of Sweden at a press conference in 24 April 2002, we have recognized this issue as a substantial matter of food safety and started analysis of acrylamide in processed foods in Japan by LC-MS/MS and GC-MS. Here we report the methods and results of the analyses on 63 samples covering 31 product types.

| | |
|-----------------------|--|
| GC-MS conditions | QP2010 (Shimadzu, Kyoto, Japan) |
| apparatus | CP-Sil 24 CB (0.25 mm i.d. x 30 m, 0.25 μ m film thickness, Varian, CA, USA) |
| column | |
| carrier gas | He |
| flow rate | 1.43 mL/min |
| injection volume | 1 μ L |
| injection temperature | 120 $^{\circ}$ C |
| temp. program | 85 $^{\circ}$ C (1 min) — (25 $^{\circ}$ C/min)—175 $^{\circ}$ C (6 min) — (40 $^{\circ}$ C/min)—250 $^{\circ}$ C (7.52 min) |
| retention time | 8.1 min |
| acquisition duration | 16 min |
| ionization | EI+ (70 eV) |
| detection | SIM (acrylamide m/z 150 & 152, internal standard m/z 153 & 155) |
| interface temp. | 280 $^{\circ}$ C |
| ion source temp. | 200 $^{\circ}$ C |
| limit of detection | 12 ng/mL (52 fmol) as 2,3-dibromopropanamide |
| limit of quantitation | 40 ng/mL (170 fmol) as 2,3-dibromopropanamide |



These analyses were carried out on single randomly selected samples from supermarkets. It is highly likely that there are variations in acrylamide concentration among production lots and among foods within a product type due to difference in the processing condition. Survey over wider range of foods with larger number of samples including home cooked foods is necessary for risk assessment.



References

- Ono H., Chuda Y., Ohnishi-Kameyama M., Yada H., Ishizaka M., Kobayashi H., and Yoshida M., Analysis of acrylamide by LC-MS/MS and GC-MS in processed Japanese foods, *Food Addit. Contam.*, 20, 215-220 (2003).
- Tareke, E., Rydberg, P., Karlsson, P., Eriksson, S., and Törnqvist M., Analysis of acrylamide, a carcinogen formed in heated foodstuffs. *J. Agric. Food Chem.*, 50, 4998-5006 (2002).
- Rosén, J. and Hellenäs, K., Analysis of acrylamide in cooked foods by liquid chromatography tandem mass spectrometry. *The Analyst*, 127, 880-882 (2002).
- Castle, L., Determination of acrylamide monomer in mushroom grown on polyacrylamide gel. *J. Agric. Food Chem.*, 41, 1261-1263 (1993).
- Castle, L., Campos, M., and Gilbert, J., Determination of acrylamide monomer in hydroponically grown tomato fruit by capillary gas chromatography-mass spectrometry. *J. Sci. Food. Agric.*, 54, 549-555 (1991).
- Food Standards Agency, Paper and board packaging: not likely to be a source of acrylamide in food. <http://www.foodstandards.gov.uk/multimedia/pdfs/27acryl.pdf>
- U. S. Food and Drug Administration, Center for Food Safety and Applied Nutrition, Detection and quantitation of acrylamide in foods (draft). <http://www.cfsan.fda.gov/~dms/acrylami.html>

The given specifications serve purely as technical information for the user. No guarantee is given on technical specifications of the described product and/or procedures.