

Application Note 159

Quantitative Analysis of Semivolatile Organic Compounds on Meridian MDN-5S

Due to the sensitivity of a mass spectrometric detector, it is imperative to use a column that will not bleed excessively when analyzing semivolatiles. The MDN™-5S offers both low bleed at high temperature and inertness.

Key Words:

● semivolatile ● GC/MS ● bleed

Analytes termed "semivolatiles" include acidic, basic and neutral compounds such as polynuclear aromatics and phthalates. Specific lists of these compounds have been identified as environmental pollutants and are routinely monitored for in soil, water, and air. A common method for the analysis of semivolatiles requires the use of a mass spectrometric detector (MSD) for positive identification and quantification. The analysis of semivolatiles requires that the column be brought to a final temperature in the range of 300°C. Due to the sensitivity of an MSD, it is imperative to use a column that will not bleed excessively at these temperatures and contaminate the ion source of the mass spectrometer.

The MDN-5S offers both low bleed at high temperature and inertness. These attributes are required to analyze semivolatile compounds. We used a 30m x 0.32mm ID, 0.50µm MDN-5S to generate a calibration curve typical for this type of analysis. Using a split injection allowed the calibration range to be run from 10 to 180ppm, and the constant flow option in the GCs EPC kept the analysis time to under 25 minutes. The TIC of an 80ppm standard is presented in Figure A. Very little column bleed was observed, even at 300 °C.

Figure B illustrates the resolution which was achieved between the difficult to separate pairs benzo(b)fluoranthene/benzo(k)fluoranthene and indeno(1,2,3-cd)pyrene/dibenzo(a,h)anthracene. The benzo (b) & (k) isomers have the same quantitation ion, which can lead to miss-identifications by many GC/MS data systems if they are not adequately resolved.

The results of the multipoint calibration are presented in Table 1. The retention times of each compound are listed along with the ions used to determine their relative response factors. Most compounds had RSD values of <15%, indicating that the column showed good linearity in this range. The inertness of the column provides adequate response for very active compounds such as 2,4-dinitrophenol, pentachlorophenol and benzidine despite the very low on-column amounts (2ng for the lowest standard and 36ng for the highest).

Figure A. Semivolatile Standard on a 30m x .32mm ID, 0.50µm Meridian™ MDN-5S Column

Sample: Semivolatile standard, 80ppm in methylene chloride
 Column: **MDN-5S, 30m x 0.32mm ID, 0.50µm**
 Cat. No.: **24393**
 Injection: 2µl, split 10:1
 Flow: 2.5 ml/min, He (constant)
 Oven: 50°C (1 min) to 120°C at 15°C/min, to 280°C at 20°C/min, to 300°C at 2 °C/min
 Inj. Temp.: 280°C
 Interface Temp.: 290°C
 Scan Range: HP 5973MSD, 45-450 amu,




Figure B. Key Semivolatile Resolutions on the 30m x 0.32mm ID, 0.50µm Meridian MDN-5S Column




Table 1. Results of Multipoint Calibration on the MDN-5S

Compound	RT	Quant. Ion	Avg. RRF	%RSD	Compound	RT	Quant Ion	Avg. RRF	%RSD
N-nitrosodimethylamine	2.48	74	0.848	4.0	3-Nitroaniline	9.68	138	0.357	4.8
Pyridine	2.54	79	1.54	4.2	Acenaphthene-d10	9.74	164		
2-Fluorophenol	3.90	112	1.247	1.8	Acenaphthene	9.78	153	1.202	6.6
Phenol	5.10	94	1.809	3.6	2,4-Dinitrophenol	9.81	184	0.152	19.8
Phenol-d6	5.09	99	1.76	3.1	4-Nitrophenol	9.87	109	0.248	14.8
Aniline	5.16	66	0.869	9.2	2,4-Dinitrotoluene	9.98	165	0.427	8.1
Bis(2-chloroethylether)	5.26	95	0.436	3.4	Dibenzofuran	10.00	168	1.788	7.4
2-Chlorophenol-d4	5.31	132	1.358	2.2	Diethyl phthalate	10.29	149	1.396	5.6
2-Chlorophenol	5.33	128	1.426	2.9	Fluorene	10.44	166	1.29	8.2
1,3-Dichlorobenzene	5.56	146	1.429	5.2	4-Chlorophenyl phenyl ether	10.44	204	0.740	8.1
1,4-Dichlorobenzene-d4	5.64	152			4-Nitroaniline	10.45	138	0.381	4.7
1,4-Dichlorobenzene	5.67	146	1.466	4.9	2-Methyl-4,6-dinitrophenol	10.49	198	0.101	37.3
Benzyl Alcohol	5.81	108	0.943	3.8	N-nitrosodiphenylamine	10.58	169	0.506	6.8
1,2-Dichlorobenzene-d4	5.86	152	0.919	6.1	Azobenzene	10.64	77	0.752	8.3
1,2-Dichlorobenzene	5.88	146	1.388	4.6	2,4,6-Tribromophenol	10.74	330	0.218	9.7
2-Methylphenol	5.96	108	1.285	4.1	4-Bromophenyl phenyl ether	11.06	248	0.215	3.6
Bis(2-chloroisopropyl)ether	6.02	45	1.214	4.1	Hexachlorobenzene	11.12	284	0.242	4.2
4-Methylphenol	6.18	108	1.401	3.9	Pentachlorophenol	11.37	266	0.129	26.9
N-nitroso-di-n-propylamine	6.21	70	1.082	4.2	Phenanthrene-d10	11.62	188		
Nitrobenzene-d5	6.42	82	0.428	4.2	Phenanthrene	11.65	178	1.028	7.8
Nitrobenzene	6.45	77	0.436	6.0	Anthracene	11.71	178	1.045	6.4
Isophorone	6.78	82	0.721	5.2	Carbazole	11.91	167	1.022	6.4
2-Nitrophenol	6.89	139	0.187	9.4	Di-n-butyl phthalate	12.35	149	1.095	7.3
2,4-Dimethylphenol	6.94	107	0.381	4.2	Fluoranthene	13.13	202	1.232	5.2
Benzoic Acid	7.03	105	0.262	22.6	Pyrene	13.41	202	1.191	6.9
bis(2-Chloroethoxy)methane	7.09	93	0.429	5.4	Benzidine	13.29	184	0.419	12.4
2,4-Dichlorophenol	7.21	162	0.297	4.9	Aramite #1	13.59	185	0.056	21.0
1,2,4-Trichlorobenzene	7.34	180	0.341	6.6	Terphenyl-d14	13.61	244	0.902	5.2
Naphthalene-d8	7.43	136			Aramite #2	13.68	185	0.075	19.7
Naphthalene	7.45	128	1.057	7.5	3,3'-Dimethylbenzidine	14.24	212	0.376	14.3
4-Chloroaniline	7.52	127	0.451	5.9	Butylbenzyl phthalate	14.26	149	0.477	6.0
Hexachlorobutadiene	7.63	225	0.208	4.6	3,3'-Dichlorobenzidine	15.10	252	0.388	11.8
4-Chloro-3-methylphenol	8.17	107	0.348	5.2	Benzo(a)anthracene	15.14	228	1.107	4.2
2-Methylnaphthalene	8.39	142	0.688	5.9	Chrysene-d12	15.15	240		
Hexachlorocyclopentadiene	8.60	237	0.377	10.9	Chrysene	15.21	228	1.028	6.1
2,4,6-Trichlorophenol	8.75	196	0.399	8.2	Bis(2-ethylhexyl)phthalate	15.21	149	0.603	9.4
2,4,5-Trichlorophenol	8.79	196	0.442	2.5	Di-n-octyl phthalate	16.62	149	1.089	22.6
2-Fluorobiphenyl	8.87	172	1.398	7.0	Benzo(b)fluoranthene	17.47	252	1.217	8.1
2-Chloronaphthalene	9.02	162	1.170	7.0	Benzo(k)fluoranthene	17.54	252	1.260	3.5
2-Nitroaniline	9.15	65	0.423	5.0	Benzo(a)pyrene	18.35	252	1.124	7.8
Dimethyl phthalate	9.39	163	1.347	5.6	Perylene-d12	18.50	265		
2,6-Dinitrotoluene	9.47	165	0.305	7.5	Indeno(1,2,3)cd-pyrene	22.01	276	1.010	12.0
Acenaphthylene	9.56	152	1.898	6.5	Dibenzo(a,h)anthracene	22.14	278	1.059	6.4
					Benzo(g,h,i)perylene	22.95	276	1.174	3.3

internal standards

Trademarks

Meridian – Sigma-Aldrich Co.
MDN – Sigma-Aldrich Co.

Fused silica columns manufactured under HP US Pat. No. 4,293,415.

Ordering Information:

Description	Cat. No.
Meridian MDN-5S Capillary Column	
30m x 0.32mm ID, 0.50µm film	24393

Contact our Technical Service Department
(phone 800-359-3041 or 814-359-3041, FAX 800-359-3044 or 814-359-5468) for expert answers to your questions.

Note 159

For more information, or current prices, contact your nearest Supelco subsidiary listed below. To obtain further contact information, visit our website (www.sigma-aldrich.com), see the Supelco catalog, or contact Supelco, Bellefonte, PA 16823-0048 USA.

ARGENTINA • Sigma-Aldrich de Argentina, S.A. • Buenos Aires 1119 AUSTRALIA • Sigma-Aldrich Pty. Ltd. • Castle Hill NSW 2154 AUSTRIA • Sigma-Aldrich Handels GmbH • A-1110 Wien
BELGIUM • Sigma-Aldrich N.V./S.A. • B-2880 Bornem BRAZIL • Sigma-Aldrich Quimica Brasil Ltda. • 01239-010 São Paulo, SP CANADA • Sigma-Aldrich Canada, Ltd. • 2149 Winston Park Dr., Oakville, ON L6H 6J8
CZECH REPUBLIC • Sigma-Aldrich s.r.o. • 186 00 Praha 8 DENMARK • Sigma-Aldrich Denmark A/S • DK-2665 Vallensbaek Strand FINLAND • Sigma-Aldrich Finland/YA-Kemia Oy • FIN-00700 Helsinki
FRANCE • Sigma-Aldrich Chimie • 38297 Saint-Quentin-Fallavier Cedex GERMANY • Sigma-Aldrich Chemie GmbH • D-82041 Deisenhofen GREECE • Sigma-Aldrich (o.m.) Ltd. • Ilioupoli 16346, Athens
HUNGARY • Sigma-Aldrich Kft. • H-1067 Budapest INDIA • Sigma-Aldrich Co. • Bangalore 560 048 IRELAND • Sigma-Aldrich Ireland Ltd. • Dublin 24 ISRAEL • Sigma Israel Chemicals Ltd. • Rehovot 76100
ITALY • Sigma-Aldrich s.r.l. • 20151 Milano JAPAN • Sigma-Aldrich Japan K.K. • Chuo-ku, Tokyo 103 KOREA • Sigma-Aldrich Korea • Seoul MALAYSIA • Sigma-Aldrich (M) Sdn. Bhd. • 58200 Kuala Lumpur
MEXICO • Sigma-Aldrich Química S.A. de C.V. • 50200 Toluca NETHERLANDS • Sigma-Aldrich Chemie BV • 3330 AA Zwijndrecht NORWAY • Sigma-Aldrich Norway • Torshov • N-0401 Oslo
POLAND • Sigma-Aldrich Sp. z o.o. • 61-663 Poznań PORTUGAL • Sigma-Aldrich Quimica, S.A. • Sintra 2710 RUSSIA • Sigma-Aldrich Russia • Moscow 103062 SINGAPORE • Sigma-Aldrich Pte. Ltd.
SOUTH AFRICA • Sigma-Aldrich (pty) Ltd. • Jet Park 1459 SPAIN • Sigma-Aldrich Quimica, S.A. • 28100 Alcobendas, Madrid SWEDEN • Sigma-Aldrich Sweden AB • 135 70 Stockholm
SWITZERLAND • Supelco • CH-9471 Buchs UNITED KINGDOM • Sigma-Aldrich Company Ltd. • Poole, Dorset BH12 4QH
UNITED STATES • Supelco • Supelco Park • Bellefonte, PA 16823-0048 • Phone 800-247-6628 or 814-359-3441 • Fax 800-447-3044 or 814-359-3044 • email: supelco@sial.com

Supelco is a member of the Sigma-Aldrich family. Supelco products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for a particular use. Additional terms and conditions may apply. Please see the reverse side of the invoice or packing slip. H DCE