

# TheReporter

Reprinted from Volume 15, No. 5, 1996

T296055

© 1999 Sigma-Aldrich Co.

For more information, or current prices, contact your nearest Supelco subsidiary listed below. To obtain further contact information, visit our website ([www.sigma-aldrich.com](http://www.sigma-aldrich.com)), see the Supelco catalog, or contact Supelco, Bellefonte, PA 16823-0048 USA.

**ARGENTINA** • Sigma-Aldrich de Argentina, S.A. • Av. Pueyrredon 2446/50 • Piso 5-B • Buenos Aires 1119  
**AUSTRALIA** • Sigma-Aldrich Pty. Ltd. • Unit #2, 14 Anella Avenue • Castle Hill NSW 2154  
**AUSTRIA** • Sigma-Aldrich Handels GmbH • Hebbelplatz 7 • A-1110 Wien  
**BELGIUM** • Sigma-Aldrich N.V./S.A. • K. Cardijnplein 8 • B-2880 Bornem  
**BRAZIL** • Sigma-Aldrich Quimica Brasil Ltda. • Rua Sabara, 566-Conj. 53 • 01239-010 São Paulo, SP  
**CANADA** • Sigma-Aldrich Canada, Ltd. • 2149 Winston Park Dr., Oakville, ON L6H 6J8  
**CZECH REPUBLIC** • Sigma-Aldrich s.r.o. • Pobrezni 46 • 186 21 Praha 8  
**DENMARK** • Sigma-Aldrich Denmark A/S • Vejlegaardsvej 65B • DK-2665 Vallensbaek Strand  
**FINLAND** • Sigma-Aldrich Finland/YA-Kemia Oy • Teerisuonkuja 4 • FIN-00700 Helsinki  
**FRANCE** • Sigma-Aldrich Chimie • Chromatographie Supelco • L'Isle d'Abeau Chesnes - B.P. 701 • 38297 Saint-Quentin Fallavier Cedex  
**GERMANY** • Sigma-Aldrich Chemie GmbH • Geschäftsbereich Supelco • Grünwalder Weg 30 • D-82041 Deisenhofen  
**GREECE** • Sigma-Aldrich (o.m.) Ltd. • 72 Argonafton Str. • 16346 Ilioupoli, Athens  
**HUNGARY** • Sigma-Aldrich Kft. • Nagy Diófa u. 7., IV fl. • H-1067 Budapest  
**INDIA** • Sigma-Aldrich Co. • Survey No. 31/1, Sitharamapalaya • Mahadevapura P.O. • Bangalore 560 048  
**IRELAND** • Sigma-Aldrich Ireland Ltd. • Airton Road • Tallaght • Dublin 24  
**ISRAEL** • Sigma Israel Chemicals Ltd. • Park Rabin • Rohovot 76100  
**ITALY** • Sigma-Aldrich s.r.l. • Via Gallarate, 154 • 20151 Milano  
**JAPAN** • Sigma-Aldrich Japan K.K. • Division Supelco • JL Nihonbashi Building • 1-10-15 Nihonbashi Horidome-cho, Chuo-ku • Tokyo 103  
**KOREA** • Sigma-Aldrich Korea • Samhan Camus Annex, 10<sup>th</sup> Floor • 17-26 Yoido-dong, Yungdeungpo-ku • Seoul  
**MALAYSIA** • Sigma-Aldrich (M) Sdn. Bhd. • 9-2, Jalan 2/128, Taman Gembira • Off Jalan Kuchai Lama • 58200 Kuala Lumpur • Selangor  
**MEXICO** • Sigma-Aldrich Quimica S.A. de C.V. • Calle 6 North No. 107 • Parque Industrial Toluca 2000 • 50200 Toluca  
**NETHERLANDS** • Sigma-Aldrich Chemie BV • Postbus 27 • 3330 AA Zwijndrecht  
**NORWAY** • Sigma-Aldrich Norway • Sandakerveien 102 • N-0483 Oslo  
**POLAND** • Sigma-Aldrich Sp. z o.o. • Szelagowska 30 • 61-626 Poznań  
**PORTUGAL** • Sigma-Aldrich Quimica, S.A. • P.O. Box 131 • Sintra 2710  
**RUSSIA** • Sigma-Aldrich Russia • TOO Techmedbiochem • Makarenko Str. 2/21 • Building 1, Flat 22 • Moscow 103062  
**SINGAPORE** • Sigma-Aldrich Pte. Ltd. • 102E Pasir Panjang Road • #08-01 Citilink Warehouse • Singapore 118529  
**SOUTH AFRICA** • Sigma-Aldrich (pty) Ltd. • CNR Kelly & Ackerman Streets • Southern Life Industrial Park Unit • Unit 16/17 • Jet Park 1459  
**SPAIN** • Sigma-Aldrich Quimica, S.A. • Apt. Correos 161 • 28100 Alcobendas, Madrid  
**SWEDEN** • Sigma-Aldrich Sweden AB • Solkraftsvägen 14C • 135 70 Stockholm  
**SWITZERLAND** • Supelco Switzerland • Industriestrasse 25 • P.O. Box 260 • CH-9471 Buchs  
**UNITED KINGDOM** • Sigma-Aldrich Company Ltd. • Supelco UK • Fancy Road, Poole • Dorset BH12 4QH  
**UNITED STATES** • Supelco • Supelco Park • Bellefonte, PA 16823-0048 • Phone 800-247-6628 or 814-359-3441 • Fax 800-447-3044 or 814-359-3044 • email:supelco@sial.com

H

This article is archived from a past issue of The Supelco Reporter. Information in the article was appropriate at the time of publication, but product specifications, catalog numbers, and availability may have changed over time.

If you have questions about applying methodology described in this article to a current application, please contact our technical service chemists.

 **SUPELCO**

ISO 9001  
REGISTERED

Supelco is a member of the Sigma-Aldrich family. Supelco products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for a particular use. Additional terms and conditions may apply. Please see the reverse side of the invoice or packing slip.

# Fast Screening of Wines Using SPME/GC

G. Vas, Research Institute for Viticulture & Enology of Agricultural Ministry, Eger, Hungary

*The ability to accurately determine the origin of a wine or detect the presence of additives is important to the wine industry. A new method of analysis was tested, consisting of sampling by headspace SPME and separation by GC using a PAG capillary GC column. This method characterized aroma patterns of two muscat-type wines from different vineyards, and revealed the presence of low concentrations of cold-pressed coriander seed oil in wine and distilled water. The resulting chromatograms showed noticeable differences between the two wines and were used to detect a 1ppm spike of additive.*


Aroma compounds of wines originate from grapes and are formed during fermentation. Their production is influenced by various factors: environment (soil, climate), grape variety, ripeness, fermentation conditions (pH, temperature, yeast flora), the wine production process (enological methods, treatment substances), and aging (bottle maturation). Wines contain aroma compounds in a wide concentration range — some components at high ppm

level, but most at low ppb or ppt levels. The low concentration of most volatile components in wine makes extraction and concentration necessary before analysis by gas chromatography (GC) or gas chromatography/mass spectrometry (GC/MS). The purpose of this study was to characterize aromas of muscat-type wines and to detect extraneous flavor addition, using the solid phase microextraction (SPME)/GC sample preparation and analysis method. Two wines from the same muscat grape variety — originating from the Eger wine region, Hungary, and the Trento wine region, Italy — were analyzed.

SPME<sup>®</sup> requires no solvents or complicated apparatus. It can be used to concentrate volatile and nonvolatile compounds in both liquid and gaseous samples. An SPME unit consists of a length of fused silica fiber coated with a phase film. The fiber is attached to a stainless steel plunger in a protective holder.

A 100µm PDMS (polydimethylsiloxane) SPME fiber was exposed to the headspace of the wine for 10 minutes at ambient temperature. After sampling, the fiber was retracted into the SPME needle, then inserted into a GC injection port and exposed for 5 minutes to desorb the analytes of interest.

**Figure A. Comparison of Muscat Wines by SPME/GC**




\*Peaks A, B, and C are characteristic of cold-pressed coriander seed oil (see Figures B and C).

Reprinted from Volume 15, No. 5, 1996


### Figure B. Muscat Wine from Eger Region Spiked with 1ppm Cold-Pressed Coriander Seed Oil

Conditions and callouts: see Figure A.



### Figure C. Distilled Water Spiked with 1ppm Cold-Pressed Coriander Seed Oil

Conditions and callouts: see Figure A.



A polyalkylene glycol (PAG) capillary GC column separated the analytes. The PAG phase has similar characteristics to the polyethylene glycol (PEG) phase, but retention indices are somewhat different. This can be advantageous for separating components that are not well resolved on a PEG column. Our results show that the PAG column is suitable for headspace SPME wine aroma characterization. Chromatograms of the two muscat wine samples show significant, characteristic differences between the respective wine regions (Figure A). This distinction is of major importance in wine analysis, providing a basis for analytical determination of the origin of a wine sample.

Using some additives (such as cold-pressed coriander seed oil) in wines is considered adulteration in many countries. Cold-pressed coriander seed oil is occasionally added to increase the muscat flavor of wine. The oil is typically used in the 10 to 20ppm concentration range—enologists easily recognize it by taste in this range. At lower concentrations it may still add a slight muscat

flavor, but is difficult to identify by taste. Sampling by headspace SPME with GC analysis, using a PAG column, can easily identify such adulteration. Figure B shows wine from the Eger region spiked with 1ppm cold-pressed coriander seed oil. The peaks marked with an asterisk (\*) are characteristic of this oil, but normally are found only in small amounts in wine. The chromatogram of pure, cold-pressed coriander seed oil in distilled water, studied under identical conditions, is shown in Figure C.

The headspace SPME sampling technique, in combination with GC analysis, easily distinguishes muscat wines of different origins and detects small concentrations of cold-pressed coriander seed oil additive. Use of this efficient and accurate sample preparation method will benefit analysts trying to ascertain or prove the origin of wines and the use of additives in wine.

#### Ordering Information:

Description	Cat. No.
<b>SPME Fiber, 100µm polydimethylsiloxane, pk. of 3</b>	
For manual sampling	57300-U
For Varian 8100/8200 AutoSampler or HPLC	57301
<b>SPME Holder<sup>■</sup></b>	
For manual sampling	57330-U
For Varian 8100/8200 AutoSampler or HPLC	57331
<b>PAG Fused Silica Capillary Column</b>	
30m x 0.25mm ID, 0.25µm film	24223

■ US patent #5,691,206. European patent # 0523092. Technology licensed exclusively to Supelco.

■ Initially you must order both holder and fiber assembly. Holder is reusable indefinitely. Use with AutoSampler requires Varian SPME upgrade kit (available from Varian).

Fused silica columns manufactured under HP US Pat. No. 4,293,415.

