



## Alcohols $C_1 - C_3$

# Separation of alcohols on a wide-bore fused silica column

## Application Note

Materials Testing & Research

### Authors

Agilent Technologies, Inc.

### Introduction

Gas chromatography using an Agilent CP-Sil 5 CB column separates 17  $C_1$  to  $C_3$  alcohols in five minutes.




**Agilent Technologies**

## Conditions

Technique : GC-capillary  
Column : Agilent CP-Sil 5 CB, 0.53 mm x 10 m fused silica  
WCOT CP-Sil 5 CB (5.0  $\mu\text{m}$ ) (Part no. CP7645)  
Temperature : 80  $^{\circ}\text{C}$   $\rightarrow$  250  $^{\circ}\text{C}$ , 10  $^{\circ}\text{C}/\text{min}$   
Carrier Gas :  $\text{N}_2$ , 10 kPa (0.1 bar), 52 cm/s  
Injector : direct  
T = 250  $^{\circ}\text{C}$   
Detector : FID, 100 x 10<sup>-12</sup> Afs  
T = 275  $^{\circ}\text{C}$   
Sample Size : 0.2  $\mu\text{L}$

## Peak identification

1. methanol
2. ethanol
3. isopropanol
4. tert. butanol
5. n-propanol
6. sec. butanol
7. isobutanol
8. 2-methyl-2-butanol
9. n-butanol
10. 3-methyl-2-butanol
11. 2-methyl-1 butanol + 3 methyl-1-butanol
12. n-pentanol + 4-methyl-2-pentanol
13. 3-hexanol
14. 2-methyl-1-pentanol
15. 2-ethyl-1-butanol
16. 2-heptanol
17. 1-heptanol



[www.agilent.com/chem](http://www.agilent.com/chem)

This information is subject to change without notice.

© Agilent Technologies, Inc. 2011

Printed in the USA

31 October, 2011

First published prior to 11 May., 2010

A00009



**Agilent Technologies**