



Buprenorphine in urine

Application Note

Forensic Toxicology

Authors

Agilent Technologies, Inc.

Introduction

An Agilent VF-DA GC column with GC/MS analyzes a sample of buprenorphine at +/- 5 ng/mL.



Agilent Technologies

Conditions



Technique : GC-capillary
Column : Agilent VF-DA, 0.20 mm x 12 m (df =optimized)
(Part no. CP8964)
Temperature : 70 °C, 1.2 min → 200 °C, 20 °C/min, → 27 °C,
7 °C/min, → 320 °C, 20 °C/min
Carrier Gas : Helium, ca. 1.0 mL/min
Pressure Program : 58.7 kPa, 2.2 min, → 97 kPa, 58 kPa/min →
132 kPa, 3 kPa/min → 180 kPa, 12 kPa/min
Injector : Splitless
Detector : MS in SIM mode
Sample Size : 1 µL
Solvent : methanol
Derivatization : as penta fluoro propionates
Courtesy : Joerg Szigan - Labor
Dr. lembke, Dr. Lempfrid - Cologne

Peak identification

1. buprenorphine -PFP

Real Sample

Method detection limit: ± 5 ng/ml



**Calibration curve 0.5-10 ng
buprenorphine**

www.agilent.com/chem

For Forensic Use only.

This information is subject to change without notice.

© Agilent Technologies, Inc. 2011

Printed in the USA

31 October, 2011

First published prior to 11 May, 2010

A01981



Agilent Technologies