


GC-MS/MS Handbook for the Analysis of Dioxins in Foods

Dioxins in
Foods



Procedures for the Analysis of Dioxins in

Extraction of Samples

SpeedExtractor E-914/E-916 (BÜCHI Labortechnik AG)



The SpeedExtractor is your best solution for fast Pressurized Solvent Extraction (PSE). Increase productivity by processing up to 6 samples in parallel. Streamline the workflow of your sample preparation with easy sample loading and ready-to-use extract collection.

Analysis of Samples

GCMS-TQ8050 (SHIMADZU)

(GC-MS/MS Method Package for Dioxins in Foods)



The method package consists of method files registered with the optimal conditions for the analysis of dioxins, as well as a report creation tool that can output the items required by EU regulations. This makes it possible to start an analysis without having to investigate analytical conditions.

Foods Using GC-MS/MS

Purification of Samples

GO-xHT (MIURA CO.,LTD.)



GO-xHT is an easy-to-operate automated system offering high throughput as well as the additional advantage of using less solvents and consumables. It provides labs with high-quality extraction as well as high return on investment thanks to an innovative flow path system.



Extraction of Samples

SpeedExtractor E-914/E-916



The SpeedExtractor is your best solution for fast Pressurized Solvent Extraction (PSE). Increase productivity by processing up to 6 samples in parallel. Streamline the workflow of your sample preparation with easy sample loading and ready-to-use extract collection.

■ Unrivaled throughput

- 6 times faster than other PSE instruments thanks to parallel extraction concept
- Extract up to 96 samples in an 8 hour shift
- Reduce number of replicates due to identical extraction conditions for up to 6 samples



Point 1

■ Ease of operation

- Unique design of extraction cells allows for ease of sample loading and cell assembly
- High level of automation saves time-consuming preparation and increases safety
- The automated sealing eliminates sources of error and guarantees safety



Point 2



Point 3

■ Low running costs

- Reduced solvent consumption and minimal need for consumables give lower running costs
- Save energy and costs with ECO mode for automated heater control



Point 1

Parallel workflow for high throughput

Significantly increase throughput with the concerted parallel extraction and parallel evaporation. Full compatibility of the SpeedExtractor collection glassware with Syncore™ Analyst and Multivapor™.

Point 2

Extraction made easy

The operation of the SpeedExtractor is easy and intuitive. Apply an established method at the push of a button and extract simultaneously up to 6 samples.



Point 3

Automated sealing


The fail-safe and reliable extraction cells seal automatically using a unique sealing principle and guarantee reproducible results.



Six Steps from the Preparation of Samples, Multi-Sample Extraction, to Multi-Sample Distillation



Step 1 Place the samples in the extraction cell.



Step 2 Attach the extraction cell to the heating block.




Step 4 After the extraction is completed, remove the collection unit.



Step 3 Close the protective shield and press the button. Six samples (or four samples) are extracted automatically.



Step 5 Transfer the samples from the collection vials to the parallel evaporator.



Step 6 Multiple samples are simultaneously distilled by the parallel evaporator.

Purification of Samples

GO-xHT



GO-xHT is an easy-to-operate automated system offering high throughput as well as the additional advantage of using less solvents and consumables. It provides labs with high-quality extraction as well as high return on investment thanks to an innovative flow path system.

■ Heating purification

• The quickness and selectivity of your current DXNs and PCBs analysis can be improved by using this technique.


■ Less than 100 mL of solvent per sample

• It uses only 90 mL of hexane and 2 mL of toluene, without dichloromethane. Thus, this system is both user and eco friendly.

■ No cross contamination in the valve-less system

• There are no solenoid valves on the DXNs and PCBs flow lines. All columns and tubings are disposable. Thus, there is no possibility of cross contamination.

*Please consult the relevant laws and regulations in your country before going for disposal.



- 1 Clean up with 90 mL of Hexane
- 2 Elute PCBs fraction with 0.9±0.1 mL of Toluene
- 3 Elute DXNs fraction with 1.3±0.1 mL of Toluene


Point 1 Heating purification

The benefit of the heated purification column is that DXNs and PCBs can be eluted quickly and a chemical reaction of sample matrices is facilitated with the chemically-modified silica gel.

A unique method for solvent flow switching with pinch valves.

Point 2 Less than 100 mL of solvent

The benefit of a heated concentration column is that the amount of toluene necessary for elution of DXNs and PCBs elution is decreased.



DXNs fraction









17 isomers of 2,3,7 & 8 substituted Chlorinated dibenzo-p-dioxins, dibenzofurans and 4 isomers of non-ortho PCBs in 1.3 ± 0.1 mL of Toluene

PCBs fraction


6 isomers of indicator PCBs, 8 isomers of mono-ortho PCBs in 0.9 ± 0.1 mL of Toluene

Four Purification Steps

Step 1 Select the column to use.

| ø18 column | Filter and Column Size | ø20 column | Filter and Column Size |
|---|--|---|--|
|  | AgNO ₃ -Silica gel ø15.6 × 81.3 1 |  | AgNO ₃ -Silica gel ø17.6 × 116.9 1 |
|  | H ₂ SO ₄ -Silica gel ø15.6 × 107 2 |  | H ₂ SO ₄ -Silica gel ø17.6 × 100 2 |
|  | Carbon ø6 × 38 3 |  | Carbon ø6 × 38 3 |
|  | Alumina ø6 × 38 4 |  | Alumina ø6 × 38 4 |

Step 2 Assemble the column and apply the sample to it.



Step 3 Set up the column to the system.



Step 4 Start purification using the method files for DXNs and PCBs. Note that the method files are prepared in advance and do not need to be created.



Analysis of Samples

GCMS-TQ8050 (Method Package for Dioxins in Foods)




The method package consists of method files registered with the optimal conditions for the analysis of dioxins, as well as a report creation tool that can output the items required by EU regulations. This makes it possible to start an analysis without the having to investigate analytical conditions.

Method files registered with the optimal conditions for the analysis of dioxins

Optimized analysis conditions (including transition and CE) are pre-registered in the method files. Additionally, the files are registered with retention times and retention indices, and the retention times can be adjusted automatically using the retention time adjustment function (AART: Automatic Adjustment of Retention Time), allowing analysis to start immediately.

Note that a comparison test (an analysis of 44 types and 201 samples of foods and feeds) with magnetic sector GC-MS was implemented using these method files and a Shimadzu TQ, and the performance was confirmed.




Report creation tool, capable of outputting items required by EU regulations

Reports of the analysis of dioxins in foods must contain the results of complicated calculations. A report creation tool is included in this product. It can automatically create reports showing items required by EU regulations. Additionally, it is capable of combining reports on DXNs and PCBs, changing the items displayed, and showing the Limit of Quantification (LOQ) calculation method, providing support to varied customer analyses.



Point 1 Automatic adjustment of retention times using the AART function


The retention times and retention indices for the target compounds are registered in the method files. The AART function adjusts the retention times for the target compounds using the retention indices for the target compounds and the results of the analysis of an n-alkane standard sample. The retention times and time programs are automatically adjusted even if the retention times for the target compounds change, such as for maintenance of the column tip.



Alkanes with different carbon numbers are investigated at a wide range of retention times, so the retention times of target compounds with a variety of boiling points can be adjusted with high accuracy.

Point 2 Report creation tool, capable of outputting items required by EU regulations

In the analysis of dioxins in foods, DXN analysis samples and PCB analysis samples are prepared from a single sample. However, depending on the pretreatment method, some of the PCBs can become mixed into the DXN sample, so the analysis results for PCBs are sometimes divided into two parts. With this product's report creation tool, even if the analysis results for PCBs are divided into two parts, they can be combined, enabling support for a variety of samples and pretreatment methods.



Reports Can Be Created Combining the Two Analysis Results

Four Steps from the Preparation for Analysis to Analysis of Samples

Step 1 Preparation for analysis

- Attach a suitable column* for the measurement compounds.
- Attach the insert liner**.
- Start up the instrument's vacuum system.

*DXNs analysis column

SH-Rxi™-5Sil MS (60 m, 0.25 mm I.D., 0.25 µm) P/N: 227-36036-02 (SHIMADZU)

PCBs analysis column

HT8-PCB (60 m, 0.25 mm I.D.) P/N: 054237 (Trajan Scientific)

BFRs analysis column


SH-Rtx™-1614 (15 m, 0.25 mm I.D., 0.1 µm) P/N: 227-36265-01 (SHIMADZU)

**Topaz® 3.5 mm I.D. Single Taper Inlet Liner w/ Wool, P/N: 23336 (Restek Corp.)



Step 2 Creation of method files

- Analyze the n-alkane standard sample.
- Adjust the retention times of the measurement compounds (AART function).
- Set the parameters for quantitation (including the calibration curve concentrations).



Step 3 Analysis of samples

- Measure and analyze the samples.



Export the results to Excel®.



Step 4 Creation of reports

- Create a report.

Analysis of PCDD/PCDF

| ID | Compound Name | LOQ | TEQ | | |
|----|---|--------|----------------|----------------|----------------|
| | | | Lowerbound | Mediumbound | Upperbound |
| 1 | 2,3,7,8-Tetrachlorodibenzo-p-dioxin | 0.0028 | 0.00000 | 0.00142 | 0.00284 |
| 2 | 1,2,3,7,8-Pentachlorodibenzo-p-dioxin | 0.0001 | 0.00500 | 0.00500 | 0.00500 |
| 3 | 1,2,3,4,7,8-Hexachlorodibenzo-p-dioxin | 0.0005 | 0.00070 | 0.00070 | 0.00070 |
| 4 | 1,2,3,6,7,8-Hexachlorodibenzo-p-dioxin | 0.0014 | 0.00080 | 0.00080 | 0.00080 |
| 5 | 1,2,3,7,8,9-Hexachlorodibenzo-p-dioxin | 0.0052 | 0.00170 | 0.00170 | 0.00170 |
| 6 | 1,2,3,4,6,7,8-Heptachlorodibenzo-p-dioxin | 0.0005 | 0.00265 | 0.00265 | 0.00265 |
| 7 | Octachlorodibenzo-p-dioxin | 0.0008 | 0.00039 | 0.00039 | 0.00039 |
| | Sum of PCDDs | | 0.01124 | 0.01266 | 0.01408 |
| 8 | 2,3,7,8-Tetrachlorodibenzofuran | 0.0034 | 0.00560 | 0.00560 | 0.00560 |
| 9 | 1,2,3,7,8-Pentachlorodibenzofuran | 0.0016 | 0.00054 | 0.00054 | 0.00054 |
| 10 | 2,3,4,7,8-Pentachlorodibenzofuran | 0.0009 | 0.00840 | 0.00840 | 0.00840 |
| 11 | 1,2,3,4,7,8-Hexachlorodibenzofuran | 0.0013 | 0.00170 | 0.00170 | 0.00170 |
| 12 | 1,2,3,6,7,8-Hexachlorodibenzofuran | 0.0012 | 0.00150 | 0.00150 | 0.00150 |
| 13 | 2,3,4,6,7,8-Hexachlorodibenzofuran | 0.0007 | 0.00090 | 0.00090 | 0.00090 |
| 14 | 1,2,3,7,8,9-Hexachlorodibenzofuran | 0.0023 | 0.00210 | 0.00210 | 0.00210 |
| 15 | 1,2,3,4,6,7,8-Heptachlorodibenzofuran | 0.0005 | 0.00041 | 0.00041 | 0.00041 |
| 16 | 1,2,3,4,7,8,9-Heptachlorodibenzofuran | 0.0004 | 0.00009 | 0.00009 | 0.00009 |
| 17 | Octachlorodibenzofuran | 0.0001 | 0.00001 | 0.00001 | 0.00001 |
| | Sum of PCDFs | | 0.02125 | 0.02125 | 0.02125 |

EU Regulation Compliant GC-MS/MS Method Package for Dioxins in Foods Specifications:

■ Product Contents

Method Files (DXNs analysis method file, PCBs analysis method file, and BFRs analysis method file)
Method File for Adjusting the Retention Times
Report Creation Tool

■ Supported Instruments

GC/MS: GCMS-TQ™8050

■ Operating Environment

Excel: Microsoft® Excel® 2016
Workstation: GCMSsolution™ Ver. 4.45 SP1 or later
LabSolutions™ Insight Ver. 3.2 SP1 or later

■ Recommended Consumables

Insert Liner: Topaz® 3.5 mm ID Single Taper Inlet Liner w/ Wool, P/N 23336 (Restek Corporation)
n-alkanes: C8-C40 Alkane calibration standard (SIGMA-ALDRICH, Cat#:40147-U)
Quantitative Retention Index Standard (C7-C33) (Restek Corporation, 31080)
Column: DXNs Analysis: SH-Rxi™-5Sil MS (60 m, 0.25 mm I.D., 0.25 m), P/N: 227-36036-02 (SHIMADZU)
PCBs Analysis: HT8-PCB (60 m, 0.25 mm I.D.), P/N: 054237 (Trajan Scientific)
BFRs Analysis: SH-Rtx™-1614 (15 m, 0.25 mm I.D., 0.1 m), P/N: 227-36265-01 (SHIMADZU)

Contact for Shimadzu and MIURA systems:



Shimadzu Europa GmbH

Albert-Hahn-Str. 6-10
D-47269 Duisburg
Tel.: +49 (0) 203 - 7687 0
Fax: +49 (0) 203 - 7666 25
info@shimadzu.eu
www.shimadzu.eu

Contact for BUCHI systems:



BÜCHI Labortechnik AG

Meierseggstrasse 40, Postfach
9230 Flawil 1, Switzerland
Tel.: +41 71 394 63 63
Fax: +41 71 394 64 64
buchi@buchi.com
www.buchi.com

Note: GO-xHT is not available in some countries/regions. Contact your Shimadzu sales representative for more details.

The n-alkanes mixed sample is used in analysis in order to adjust the measurement parameters and retention times in the method files using the AART function. BFR compounds with carbon Nos. C33 and higher are included.

P/N 31080: for customers who only analyze DXNs and PCBs.

P/N 40147-U: for customers who analyze BFRs.

GCMS-TQ, GCMSsolution, and LabSolutions Insight are trademarks of Shimadzu Corporation.

Topaz, Rxi, and Rtx are registered trademarks of Restek Corporation.

Syncore and Multivapor are trademarks of BUCHI Corporation.

Microsoft and Excel are either registered trademarks or trademarks of Microsoft Corp. in the United States and/or other countries.



Shimadzu Corporation

www.shimadzu.com/an/

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

Company names, products/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol "TM" or "®".

Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The contents of this publication are provided to you "as is" without warranty of any kind, and are subject to change without notice. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication.