



**New Application Specific Columns for Semi-volatile and ECD Methods  
Lead the Way to Higher Productivity in Environmental Laboratories**

Ken Lynam  
CSD Application Chemist

# Presentation Outline

- DB-UI 8270D columns for semi-volatile analysis
  - Testing procedures
  - Example chromatograms
  - Addition productivity enhancing tools
- DB-CLP1 and DB-CLP2 columns for dual ECD methods
  - Testing procedures
  - Example chromatograms
  - Additional productivity enhancing tools


# Agilent J&W DB-UI 8270D Design

Designed specifically for environmental laboratories performing regulated GC/MS semi-volatiles analysis needing improved resolution and peak shape for active target compounds, DB-UI 8270D columns are engineered and special tested to provide added assurance of optimal peak shape performance for trace level active semivolatile and priority pollutants to ensure confidence in analytical results.

Agilent DB-UI 8270D columns are manufactured with world class quality, column to column consistency and worldwide availability. Stringent testing included with each column provides proof of performance.




# DB-UI 8270D Example Test Chromatogram



1	Propionic acid	5	n-Octane	9	p-Xylene
2	Pyridine	6	1,2-Butanediol	10	2-Heptanone
3	1-Pentanol	7	1-Chloro-2-flourobenzene	11	n-Nonane
4	1-Octene	8	m-Xylene	12	Isopropylbenzene


## Testing conditions and probe significance.

- Ultra low volume solvent test mix eliminates solvent masking of activity
- Low concentration test probe mixture (~5-10 ng on column)
- 45°C low test temp more stringent test for inertness by allowing us to test small molecule (stronger) acid probe, smaller molecule (stronger) basic probe and smaller molecule (more demanding) primary alcohol plus a demanding –diol to measure silanol activity.
- m & p- xylene are surrogate representatives of aromatic isomers and optimized phase selectivity (RI measurement) for getting highest resolution (resolution measurement) for semi-volatile aromatic isomers.
- Halogenated aromatic to measure for reproducible selectivity of the column for similar semivolatiles found in typical analyses.





# 78 Semi-Volatile Components on a DB-UI 8270D 30m x 0.25 x 0.25




## 25 minute Semi-Volatile Analysis

Application note 5991-0250EN



# Resolution of benzo-b & k fluoranthene isomers




## Positional isomers

Column: Agilent DB-UI 8270D, 30 m x 0.25 mm, 0.25  $\mu$ m (p/n 122-9732)  
Liner: Dual taper direct connect liner (p/n G1544-80700)  
Inlet: MMI in non-pulsed splitless mode 1  $\mu$ L at 275  $^{\circ}$ C  
Carrier: He, 1.2 mL/min, constant flow  
Septum purge flow: 3 mL/min, purge time on 0.7 min 50 mL/min  
Oven program: 30  $^{\circ}$ C (1.0 min), 15  $^{\circ}$ C/min to 100  $^{\circ}$ C, 20  $^{\circ}$ C/min to 240  $^{\circ}$ C (0.5 min), 15  $^{\circ}$ C/min to 325  $^{\circ}$ C (6.7 min) Gas saver Off  
GC/MSD: Agilent 7890/5975C, 325  $^{\circ}$ C transfer line, 280  $^{\circ}$ C source, 150  $^{\circ}$ C quad, 35-500 AMU range  
Sampler: Agilent 7693, 10.0  $\mu$ L syringe (p/n G4513-80216)



# 10 ng/ul Semivolatile Checkout Standard on 20m x 0.18mm x 0.36um DB-8270D Capillary GC Column using an Ultra Inert Liner with Wool



1	N-Nitrosodimethylamine	11	4-nitrophenol	21	Phenanthrene-d10
2	Aniline	12	2,4-dinitrotoluene	22	Aldrin
3	1,4-Dichlorobenzene-d4	13	Flourene	23	Heptachlor epoxide
4	Isophorone	14	4,6-dinitro-2-methyl phenol	24	Endrin
5	1,3-dimethyl-2-nitrobenzene	15	Trifluralin	25	4,4'-DDT
6	Naphthalene	16	Simazine	26	3,3'-dichlorobenzidine
7	hexachlorocyclopentadiene	17	Atrazine	27	Chrysene d-12
8	Mevinphos	18	pentachlorophenol	28	benzo[b]flouranthene
9	Acenaphthene-d10	19	Terbufos	29	Perylene-d12
10	2,4-dinitrophenol	20	Chlorothanlonil		


## 16 minute Semi-Volatile Analysis

Checkout Standard mix Agilent part # 5190-0473





## Excellent Peak Shape and Response for 2,4 Dinitrophenol



Peak #	Compound
9	Acenaphthalene D10
10	2,4 Dinitrophenol
11	4 Nitro phenol
12	2,4 Dinitro toluene

Column 1: Agilent J&W 8270-D 20 m x 0.18mm x 0.36 um

Column 2: 1 m 0.15 mm ID deactivated 2 psi with 5ml/min bleeder

Carrier: 1 ml/min constant flow at 40°C, AUX EPC 2 PSI constant pressure

Backflush: postrun 75 psi AUX EPC 2 PSI inlet

Oven: 40°C (2.5min) to 320°C (25°C/min) 4.8 min hold

Injection: 280°C Pulsed splitless 44 psi until 0.3 min, purge flow 60ml/min on at 0.35 min

Liner: UI single taper

MSD: Transfer Line 290°C, Source 300°C, Quad 150°C 50 -550 AMU




9



# Additional productivity enhancing tools

- Free method translation software
- Ultra Inert Liners with “touch-less packaging”
- Gas Clean filters
- CFT fittings and devices
- Turn key analyzer solutions



# DB-UI 8270D

- World wide recognized EPA method for Semi-volatiles analysis
- Tested for optimal inertness, selectivity and low bleed
- Available in convenient 6-packs in the US for high volume labs
  - 20 m x 0.18mm x 0.36um format
  - 30m x 0.25mm x 0.25um format
  - Bulk discounting 6 for the price of 5 (US only)

DB-UI 8270D flyer



# Method Translation Software

GC Method Translation

Criterion:  Translate Only  Best Efficiency  Fast Analysis  None **Speed gain: 1.32526**

		Original Method	Translated Method				
<b>Column</b>							
Length,	m	30	<input type="checkbox"/> 20				
Internal Diameter,	µm	316	<input type="checkbox"/> 177				
<b>Film</b>							
Thickness,	µm	0.25	<input checked="" type="radio"/> Unlock				
Phase Ratio		316.0	<input type="radio"/> 0.18 <input type="radio"/> 245.8				
<b>Carrier Gas</b>							
Enter one Setpoint			<input checked="" type="checkbox"/> Helium				
Head Pressure,	psi	13.126	33.647				
Flow Rate,	mLn/min	2.0176	1.1301				
Outlet Velocity,	cm/sec	56.72	101.35				
Average Velocity,	cm/sec	38	43.16				
Hold-up Time,	min	1.31579	0.772393				
Outlet Pressure (absolute),	psi	14.696	<input type="checkbox"/> 14.696				
Ambient Pressure (absolute),	psi	14.696	<input type="checkbox"/> 14.696				
<b>Oven Temperature</b> 3-ramp Program							
		Ramp Rate	Final Temp.	Final Time	Ramp Rate	Final Temp.	Final Time
		°C/min	°C	min	°C/min	°C	min
Initial			120.00	1.17		120.00	0.883
Ramp 1		25	160	0	33.132	160	0.000
Ramp 2		10	260	0	13.253	260	0.000
Ramp 3		15	300	4	19.879	300	3.018
<b>Sample Information</b> None							

- Four translation modes
    - Translation only
    - Best efficiency
    - Fast analysis
    - None (unlock all carrier gas parameters)
- } Lock all carrier gas parameters, making the flow rate an independent parameter.

- If translating to a different ID column, **phase ratio should be maintained** for the **most reliable results**
- If there are significant differences in phase ratio, Method Translation Software can still be used but **elution order should be confirmed**.
- **Stationary phase** of a new column must be the **same** as the original – the Method Translation Software cannot account for differences in selectivity.



# Method Translation Made Simple

- Agilent Method Translation Software
  - A **FREE, stand-alone** software running on a PC
  - Method Translation allows easy “tweaking” of method parameters to speed up run time.
  - It preserves relative elution order by properly scaling gas velocity and temperature program to maintain equivalent elution pattern.
- Can Change
  - column dimensions (i.e. ID, column length, film thickness, or phase ratio)
  - carrier gas type (i.e. hydrogen, helium, nitrogen)
  - pneumatic set points (i.e. flow rates, head pressure, or holdup time).
  - The software then generate a translated method (new temperature program), which will attempt to maintain the resolution and selectivity of the original method.
- Benefits
  - reduces method development time
  - helps assess if GC method is compatible with hardware



## Useful References

Application Note: 5991-0250EN Semi-Volatile Analysis with Specifically Designed DB-UI 8270D Columns

Why Inertness Matter in Gas Phase Analysis, Poster at Pittcon 2012

Method Translation in Gas Chromatography LC/GC article July, 2011

**[GC Method Translation Software - Agilent - Agilent Technologies](#)**

*[www.chem.agilent.com/en-US/Support/.../GcMethodTranslation.aspx](http://www.chem.agilent.com/en-US/Support/.../GcMethodTranslation.aspx)*

*You tube video on Method Translation software:*


*[www.youtube.com/watch?v=bRupj5925\\_Y](http://www.youtube.com/watch?v=bRupj5925_Y)*





# Ultra Inert Deactivated Liners

- Agilent owns, manages, and does the deactivation in our Folsom facility
  - New manufacturing processes based on chemistry of **J&W Ultra Inert GC columns**
  - Elimination of fatty acid contamination – a common side effect of glass wool manufacturing
- Unequalled Manufacturing Quality Control
  - Lot testing ensures reproducible coverage of deactivation
  - QC test with probes selected to reveal activity
  - GC method tailored to test liner -- not column or system -- inertness
  - Lab notebook friendly **Certificate on a sticker** shipped with each liner
- New patent pending “Touchless” packaging...
  - Plasma treated Non-Stick **O-ring is preinstalled** on the liner
  - Packaging is Pharmaceutical grade PTEG tubing approved by GCMS extraction testing
  - Install new liner and O-ring without touching – or risk contaminating – the new, clean liner





# Reliability / Quality Assurance : Ultra Inert Liner Certificate of Performance

Lot to Lot Liner Reproducibility assured:

Each deactivation lot is Certified to ensure consistent and efficient coverage using both acidic and basic probes at trace (2 ng) levels on column

Certificate with every liner is printed on a label ready to peel and stick into analysts' laboratory notebooks for easier compliance.

Traceability:

Deactivation Lot number is on Certificate

Liner lot number (and part number) is permanently etched on glass







# Ultra Inert liners and packaging – to meet the needs of your lab

- 4 very popular liner configurations
- Splitless single taper without wool
- Splitless single taper With wool
- Split straight with wool
- Split low pressure drop single taper with wool



Convenient 5 pack in Touchless package




Individual liner in Touchless packaging



25 pack of liners in TouchLess Dispenser

# DB-CLP1 Example Test Chromatogram



1. BENZYLAMINE

2. TRIDECANE

3. 4-CHLOROPHENOL

4. TETRADECANE

5. DICHLORVOS

6. 1-UNDECANOL

7. EPTC


8. PENTADECANE

9. 1,2,3,4-TETRACHLORO BENZENE

10. BIPHENYL

11. HEXADECANE

# DB-CLP2 example Test chromatogram



- |                         |                        |                                |
|-------------------------|------------------------|--------------------------------|
| 1. 2-ETHYLHEXANOIC ACID | 5. TRIDECANE           | 9. TETRADECANE                 |
| 2. 1,6-HEXANEDIOL       | 6. 2-METHYLNAPHTHALENE | 10. 1,2,3,4-TETRACHLOROBENZENE |
| 3. 4-CHLOROPHENOL       | 7. EPTC                | 11. DICYCLOHEXYLAMINE          |
| 4. DICHLORVOS           | 8. 1-UNDECANOL         |                                |

## EPA Methods Covered by DB-CLP1 & DB-CLP2 columns

- **CLP Pesticides** 19 Targeted Organohalide pesticides
- **EPA 504.1** 1,2-Dibromoethane (EDB), 1,2-Dibromo-3-chloropropane (DBCP), and 1,2,3-Trichloropropane (123TCP)
- **EPA 505** Organohalide pesticides
- **EPA 508.1** Chlorinated pesticides, herbicides, and organohalides
- **EPA 551** Chlorinated solvents, Trihalomethanes (THMs), and Disinfection Byproducts (DBPs)
- **EPA 552.3** Haloacetic acids and Dalapon
- **EPA 8081B (extended analyte list)** Organochlorine pesticides
- **EPA 8082A** Polychlorinated biphenyls (PCBs) and Aroclors
- **EPA 8151A** Chlorophenoxyacid herbicides

# Agilent J&W DB-CLP1 and DB-CLP2 Columns

- **GC/ECD Analyses of Chlorinated Pesticides, Herbicides, Organohalides, and PCBs**
  - Contract Lab Program (CLP) targeted pesticides
  - Multiple EPA Methods
    - Testing contaminants in drinking water and groundwater
  - Dual Column analysis: Primary and Confirmation Column Sets
- **Issues**
  - Many methods use different column sets to achieve desired separation
  - Downtime of switching between different column sets for different methods
- **Solution**
  - **Agilent's J&W DB-CLP1 and DB-CLP2 column set**
    - Provides the selectivity to separate many different halogenated compounds
    - One column set is able to analyze via multiple methods



# Fast CLP Pesticides


## Chlorinated Pesticides

DB-CLP1 30m 0.32mm ID 0.25µm (cat.# 123-8232)  
 DB-CLP2 30m 0.32mm ID 0.5µm (cat.# 123-8336)  
 5m x 0.32mm ID deact. guard column  
 Inert Tee CFT device (cat.# G3184-60065) 1:1 Split


Instrument: Agilent 7890 GC with dual µECD  
 Sampler: Agilent 7693  
 Sample: 50 ng/mL CLP Pesticides  
 Inj. Vol.: 1 µL splitless  
 Liner: Ultra Inert liner, single taper splitless  
 (cat.#5190-2292)  
 Inj. Temp.: 250°C  
 Oven Temp: 150°C (hold 0.2 min),  
 45°C/min to 250°C, 18°C/min to 300°C,  
 30°C/min to 330°C, hold 2.5 min  
 Carrier Gas: Helium, constant flow 3.5 mL/min  
 Detector: µECD @340 °C

- |                          |                  |                         |
|--------------------------|------------------|-------------------------|
| 1. Tetrachloro-m-xylene* | 9. γ-Chlordane   | 16. Endosulfan II       |
| 2. α-BHC                 | 10. α-Chlordane  | 17. 4,4'-DDT            |
| 3. γ-BHC                 | 11. Endosulfan I | 18. Endrin aldehyde     |
| 4. β-BHC                 | 12. 4,4'-DDE     | 19. Endosulfan sulfate  |
| 5. Heptachlor            | 13. Dieldrin     | 20. Methoxychlor        |
| 6. δ-BHC                 | 14. Endrin       | 21. Endrin ketone       |
| 7. Aldrin                | 15. 4,4'-DDD     | 22. Decachlorobiphenyl* |
| 8. Heptachlor epoxide    |                  | *surrogate standard     |

Agilent DB-CLP1



Agilent DB-CLP2





# EPA Method 504.1


## EDB, DBCP, and 123TCP

DB-CLP1 30m 0.32mm ID 0.25µm (cat.# 123-8232)  
 DB-CLP2 30m 0.32mm ID 0.5µm (cat.# 123-8336)  
 5m x 0.32mm ID deact. guard column  
 Inert Tee CFT device (cat.# G3184-60065) 1:1 Split


Instrument: Agilent 7890 GC with dual µECD  
 Sampler: Agilent 7693  
 Sample: 100 ng/ mL 504.1 analytes  
 100 ng/ mL Chlorinated Solvents + Trihalomethanes  
 Inj. Vol.: 2 µL splitless  
 Liner: Ultra Inert liner, single taper splitless  
 (cat.#5190-2292)  
 Inj. Temp.: 200°C  
 Oven Temp: 50°C (hold 1.5 min) to 95°C at 20°C/ min,  
 40°C/ min to 175°C, hold 1.25 min  
 Carrier Gas: Helium, constant flow 3.75 mL/ min  
 Detector: µECD @300 °C

- |                          |   |
|--------------------------|---|
| 1. Chloroform            | 7. 1,1,2-Trichloroethane                      |
| 2. 1,1,1-Trichloroethane | 8. Dibromochloromethane                       |
| 3. Carbon tetrachloride  | <b>9. 1,2-Dibromoethane (EDB)</b>             |
| 4. Trichloroethene       | 10. Bromoform                                 |
| 5. Bromodichloromethane  | <b>11. 1,2,3-Trichloropropane (123TCP)</b>    |
| 6. Tetrachloroethene     | <b>12. 1,2-Dibromo-3-chloropropane (DBCP)</b> |

Agilent DB-CLP1



Agilent DB-CLP2





# EPA Method 505


## Organohalide Pesticides

DB-CLP1 30m 0.32mm ID 0.25µm (cat.# 123-8232)  
DB-CLP2 30m 0.32mm ID 0.5µm (cat.# 123-8336)  
5m x 0.32mm ID deact. guard column  
Inert Tee CFT device (cat.# G3184-60065) 1:1 Split


Instrument: Agilent 7890 GC with dual µEOD  
Sampler: Agilent 7873B  
Sample: 100 ng/ mL 505 analytes  
Inj. Vol.: 2 µL splitless  
Liner: Ultra Inert Liner double taper, splitless (cat.# 5190-3983)  
Inj. Temp.: 250°C  
Oven Temp: 90°C (hold 0.5 min), 35°C/ min to 175°C, 12°C/ min to 300°C, hold 1.75 min  
Carrier Gas: Helium, constant flow 2.5 mL/ min  
Detector: µEOD @325 °C

- |                               |                              |
|-------------------------------|------------------------------|
| 1. Hexachlorocyclopentadienol | Heptachlor epoxide           |
| 2. Hexachlorobenzene          | 10. γ-Chlordane              |
| 3. Atrazine (25 µg/mL)        | 11. trans-Nonachlor          |
| 4. Simazine (25 µg/mL)        | 12. α-Chlordane              |
| 5. γ-BHC                      | 13. Dieldrin                 |
| 6. Heptachlor                 | 14. Endrin                   |
| 7. Alachlor (1 µg/mL)         | 15. cis-Nonachlor            |
| 8. Aldrin                     | 16. Methoxychlor (0.5 µg/mL) |

Agilent DB-CLP1



Agilent DB-CLP2







# EPA Method 508.1


## Chlorinated Pesticides and Herbicides

DB-CLP1 30m 0.32mm ID 0.25µm (cat.# 123-8232)  
 DB-CLP2 30m 0.32mm ID 0.5µm (cat.# 123-8336)  
 5m x 0.32mm ID deact. guard column  
 Inert Tee CFT device (cat.# G3184-60065) 1:1 Split


Instrument: Agilent 7890 GC with dual µEOD  
 Sampler: Agilent 7873B  
 Sample: 100 ng/ mL 508.1 analytes  
 100 ng/ mL Pesticide Surrogate Mix  
 Inj. Vol.: 2 µL splitless  
 Liner: Ultra Inert liner, single taper splitless  
 (cat.# 5190-2292)  
 Inj. Temp.: 250°C  
 Oven Temp: 80°C (hold 0.5 min) to 175°C at 26°C/ min,  
 6.5°C/ min to 235°C, 15°C/ min to 300°C,  
 hold 6 min  
 Carrier Gas: Helium at 35 cm/ sec (constant flow)  
 Detector: µEOD @340 °C

1. Hexachlorocyclopentadiene	11. γ-BHC	21. Heptachlor epoxide	31. 4,4'-DDD
2. Etriazole	12. β-BHC	22. Cyanazine	32. Endosulfan II
3. Chloroneb	13. Heptachlor	23. Butachlor	33. 4,4'-DDT
4. Trifluralin	14. Alachlor	24. γ-Chlordane	34. Endrin aldehyde
5. TCMX*	15. δ-BHC	25. α-Chlordane	35. Endosulfan sulfate
6. Propachlor	16. Chlorothalonil	26. Endosulfan I	36. Methoxychlor
7. Hexachlorobenzene	17. Aldrin	27. 4,4'-DDE	37. cis-Permethrin
8. α-BHC	18. Metribuzin	28. Dieldrin	38. trans-Permethrin
9. Atrazine	19. Metolachlor	29. Chlorobenzilate	39. Decachlorobiphenyl*
10. Simazine	20. DCPA	30. Endrin	<i>*surrogate standard</i>

Agilent DB-CLP1



Agilent DB-CLP2



# EPA Method 551


## Chlorinated Solvents, THMs, and DBPs

DB-CLP1 30m 0.32mm ID 0.25µm (cat.# 123-8232)  
 DB-CLP2 30m 0.32mm ID 0.5µm (cat.# 123-8336)  
 5m x 0.32mm ID deact. guard column  
 Inert Tee CFT device (cat.# G3184-60065) 1:1 Split


Instrument: Agilent 7890 GC with dual µECD  
 Sampler: Agilent 7693  
 Sample: 100 ng/ mL Chlorinated Solvents, THMs, and DPBs  
 Inj. Vol.: 2 µL splitless  
 Liner: Ultra Inert liner, single taper splitless  
 (cat.# 5190-2292)  
 Inj. Temp.: 200°C  
 Oven Temp.: 35°C (hold 5.75 min) to 95°C at 20°C/ min,  
 40°C/ min to 200°C, hold 1.25 min  
 Carrier Gas: Helium at 45 cm/ sec (constant flow)  
 Detector: µECD @300 °C

- |                             |                                 |
|-----------------------------|---------------------------------|
| 1. Chloroform               | 11. Tetrachloroethene           |
| 2. 1,1,1-Trichloroethane    | 12. 1,1,2-Trichloroethane       |
| 3. Carbon tetrachloride     | 13. Dibromochloromethane        |
| 4. Trichloroacetonitrile    | 14. 1,2-Dibromoethane           |
| 5. Trichloroethene          | 15. 1,1,1-Trichloro-2-propanone |
| 6. Chloral hydrate          | 16. Bromochloroacetonitrile     |
| 7. Bromodichloromethane     | 17. Bromoform                   |
| 8. 1,1-Dichloro-2-propanone | 18. 1,2,3-Trichloropropane      |
| 9. Dichloroacetonitrile     | 19. Dibromoacetonitrile         |
| 10. Chloropicrin            | 20. 1,2-Dibromo-3-chloropropane |

Agilent DB-CLP1



Agilent DB-CLP2





# EPA Method 552.3


## Haloacetic Acids and Dalapon

DB-CLP1 30m 0.32mm ID 0.25µm (cat.# 123-8232)  
DB-CLP2 30m 0.32mm ID 0.5µm (cat.# 123-8336)  
5m x 0.32mm ID deact. guard column  
Inert Tee CFT device (cat.# G3184-60065) 1:1 Split


Instrument: Agilent 7890 GC with dual µEOD  
Sampler: Agilent 7693  
Sample: 10-100 ng/ mL Haloacetic acids and Dalapon (methyl esters)  
Inj. Vol.: 1 µL splitless  
Liner: Ultra Inert liner, single taper splitless (cat.# 5190-2292)  
Inj. Temp.: 180°C  
Oven Temp.: 40°C (hold 0.5 min) to 95°C at 10°C/ min, 30°C/ min to 200°C, hold 1 min  
Carrier Gas: Helium at 54.79 cm<sup>3</sup>/ sec (constant flow)  
Detector: µEOD @340 °C

- |  |   |
|--|---|
| 1. Methyl chloroacetate (30ng/mL)        | 7. Methyl bromochloroacetate (20ng/mL)    |
| 2. Methyl bromoacetate (20ng/mL)         | 8. Methyl bromodichloroacetate (20ng/mL)  |
| 3. Methyl dichloroacetate (30ng/mL)      | 9. Methyl dibromoacetate (10ng/mL)        |
| 4. Dalapon methyl ester (20ng/mL)        | 10. Methyl dibromochloroacetate (50ng/mL) |
| 5. Methyl trichloroacetate (10ng/mL)     | 11. Methyl tribromoacetate (100ng/mL)     |
| 6. 1,2,3-Trichloropropane (IS) (50ng/mL) |   |

Agilent DB-CLP1



Agilent DB-CLP2





# EPA Method 8081B (extended)


## Organochlorine Pesticides

DB-CLP1 30m 0.32mm ID 0.25µm (cat.# 123-8232)  
 DB-CLP2 30m 0.32mm ID 0.5µm (cat.# 123-8336)  
 5m x 0.32mm ID deact. guard column  
 Inert Tee CFT device (cat.# G3184-60065) 1:1 Split


Instrument: Agilent 7890 GC with dual µECD  
 Sampler: Agilent 7693  
 Sample: 50 ng/ mL 8081B analytes  
 Inj. Vol.: 2 µL splitless  
 Liner: Ultra Inert liner, single taper splitless (cat.#5190-2292)  
 Inj. Temp.: 250°C  
 Oven Temp: 80°C (hold 0.5 min) to 150°C at 20°C/ min, 5°C/ min to 235°C, 15°C/ min to 300°C, hold 5 min  
 Carrier Gas: Helium at 43.5 cm/ sec (constant flow)  
 Detector: µECD @325 °C

1. 1,2-Dibromo-3-chloropropane	11. α-BHC	21. DCPA	31. Chlorobenzilate (250ng/mL)	41. Captafol
2. Hexachlorocyclopentadiene	12. Pentachloronitrobenzene	22. Isodrin	32. Perthane (250ng/mL)	42. Methoxychlor
3. 1-Bromo-2-nitrobenzene	13. γ-BHC	23. Kelthane	33. Chloropropylate (250ng/mL)	43. Endrin ketone
4. Etriazole	14. β-BHC	24. Heptachlor epoxide	34. Endrin	44. Mirex
5. Chloroneb	15. Heptachlor	25. γ-Chlordane	35. Nitrofen	45. cis-Permethrin
6. Trifluralin	16. Dichlone	26. trans-Nonachlor	36. 4,4'-DDD	46. trans-Permethrin
7. TCMX*	17. Alachlor	27. α-Chlordane	37. Endosulfan II	47. Decachlorobiphenyl*
8. Propachlor	18. δ-BHC	28. Endosulfan I	38. 4,4'-DDT	
9. Diallate isomers (250ng/mL)	19. Chlorothalonil	29. 4,4'-DDE	39. Endrin aldehyde	
10. Hexachlorobenzene	20. Aldrin	30. Dieldrin	40. Endosulfan sulfate	

Agilent DB-CLP1



Agilent DB-CLP2





# EPA Method 8082A


## PCB Congeners

DB-CLP1 30m 0.32mm ID 0.25µm (cat.# 123-8232)  
 DB-CLP2 30m 0.32mm ID 0.5µm (cat.# 123-8336)  
 5m x 0.32mm ID deact. guard column  
 Inert Tee CFT device (cat.# G3184-60065) 1:1 Split


Instrument: Agilent 7890 GC with dual µEOD  
 Sampler: Agilent 7693  
 Sample: 100 ng/ mL PCB congeners  
 100 ng/ mL Pesticide Surrogate Mix  
 Inj. Vol.: 2 µL splitless  
 Liner: Ultra Inert liner, single taper splitless  
 (cat.#5190-2292)  
 Inj. Temp.: 250°C  
 Oven Temp: 125°C (hold 0.25 min),  
 20°C/ min to 210°C, hold 0.5 min  
 7°C/ min to 235°C, hold 0.75 min,  
 25°C/ min to 325°C, hold 2 min,  
 Carrier Gas: Helium, constant flow 3 mL/ min  
 Detector: µEOD @340 °C

- |           |             |                         |
|-----------|-------------|-------------------------|
| 1. BZ #1  | 8. BZ #66   | 15. BZ #138             |
| 2. TCMX*  | 9. BZ #101  | 16. BZ #187             |
| 3. BZ #5  | 10. BZ #87  | 17. BZ #183             |
| 4. BZ #18 | 11. BZ #110 | 18. BZ #180             |
| 5. BZ #31 | 12. BZ #151 | 19. BZ #170             |
| 6. BZ #52 | 13. BZ #153 | 20. BZ #206             |
| 7. BZ #44 | 14. BZ #141 | 21. Decachlorobiphenyl* |
- \*surrogate standard

Agilent DB-CLP1



Agilent DB-CLP2





# EPA Method 8082A


## Aroclors

DB-CLP1 30m 0.32mm ID 0.25µm (cat.# 123-8232)  
 DB-CLP2 30m 0.32mm ID 0.5µm (cat.# 123-8336)  
 5m x 0.32mm ID deact. guard column  
 Inert Tee CFT device (cat.# G3184-60065) 1:1 Split

Instrument: Agilent 7890 GC with dual µEOD  
 Sampler: Agilent 7693  
 Sample: 100 ng/ mL Individual Aroclor  
 100 ng/ mL Decachlorobiphenyl  
 Inj. Vol.: 2 µL splitless  
 Liner: Ultra Inert liner, single taper splitless  
 (cat.#5190-2292)  
 Inj. Temp.: 250°C  
 Oven Temp.: 160°C (hold 0.25 min),  
 35°C/ min to 330°C, hold 3 min  
 Carrier Gas: Helium, constant flow 3.75 mL/ min  
 Detector: µEOD @340 °C

### Agilent DB-CLP1

### Agilent DB-CLP2





# EPA Method 8151A


## Chlorophenoxyacid Herbicides

DB-CLP1 30m 0.32mm ID 0.25µm (cat.# 123-8232)  
 DB-CLP2 30m 0.32mm ID 0.5µm (cat.# 123-8336)  
 5m x 0.32mm ID deact. guard column  
 Inert Tee CFT device (cat.# G3184-60065) 1:1 Split


Instrument: Agilent 7890 GC with dual µEOD  
 Sampler: Agilent 7693  
 Sample: 50 ng/ mL Chlorophenoxyacid Herbicides (methyl esters)  
 50 ng/ mL DCAA methyl ester  
 50 ng/ mL 4,4'-Dibromooctafluorobiphenyl  
 Inj. Vol.: 2 µL splitless  
 Liner: Ultra Inert liner, single taper splitless (cat.#5190-2292)  
 Inj. Temp.: 250°C  
 Oven Temp: 85°C (hold 0.25 min) to 135°C at 25°C/ min, 11°C/ min to 200°C, 20°C/ min to 300°C, hold 1.5 min  
 Carrier Gas: Helium, constant flow 3 mL/ min  
 Detector: µEOD @325 °C

- |  |                              |
|--|------------------------------|
| 1. Dalapon methyl ester                  | 11. Pentachloroanisole       |
| 2. 3,5-dichlorobenzoic acid methyl ester | 12. 2,4,5-TP methyl ester    |
| 3. 4-Nitroanisole                        | 13. 2,4,5-T methyl ester     |
| 4. DCAA methyl ester (SS)                | 14. Chloramben methyl ester  |
| 5. Dicamba methyl ester                  | 15. Dinoseb methyl ester     |
| 6. MCPP methyl ester                     | 16. 2,4-DB methyl ester      |
| 7. MCPA methyl ester                     | 17. Bentazon methyl ester    |
| 8. 4,4'-DBOB (IS)                        | 18. DCPA methyl ester        |
| 9. Dichlorprop methyl ester              | 19. Picloram methyl ester    |
| 10. 2,4-D methyl ester                   | 20. Acifluorfen methyl ester |

Agilent DB-CLP1



Agilent DB-CLP2



# Summary


- DB-UI8270D columns for semi-volatile analysis
  - Specifically tested for optimal semi-volatile selectivity and inertness performance
- DB-CLP1 and DB-CLP2 columns for dual ECD methods
  - Specifically tested for optimal CLP method performance, cover 9 EPA methods

Additional productivity enhancing tools work for both semi-volatiles and CLP methods



# Additional productivity enhancing tools

- Free method translation software
- Ultra Inert Liners with “touch-less packaging”
- Gas Clean filters
- CFT fittings and devices
- Turn key analyzer solutions



## More Useful References

Application Note: 5990-0207EN Organohalide Pesticides in Water by GC/ $\mu$ ECD with Agilent J&W DB-CLP1 and DB0CLP2

Application Note: 5991-0541EN Rapid Analysis of CLP Pesticides Using High-Temperature Agilent J&W DB-CLP-1 and DB-CLP2 Columns

Application Note: 5991-0615EN Evaluating CLP and EPA Methods for Pesticides in Water Using Agilent J&W DB-CLP1/DB-CLP2 GC Columns

Agilent J&W Column Selection Guide 5990-9867EN

Stop back to our website often for the latest updates



# New Tools for Success

- **Agilent's J&W DB-8270D**
  - Designed and tested for optimal semi-volatile analysis performance
- **Agilent's J&W DB-CLP1 and DB-CLP2 column set**
  - Analyze nine EPA methods along with CLP pesticides without changing column sets
  - Increase sample throughput with faster analysis times
- **Flow-path productivity enhancing for all**

# Questions?

