

Waters

Multi-Residue
Pesticide Analysis of
Food Matrices using
GC/MS/MS and
LC/MS/MS

For Complete  Confidence

- Introduction to pesticide residue analysis
- GC/MS/MS based methods
 - Quattro micro GC for complex matrices
 - 100 pesticide residues in a single analysis
- LC/MS/MS based methods
 - 100 pesticide residues using the Quattro Premier API
 - Enhancing the method
 - ACQUITY UPLC™
- TargetLynx Data Management Software

Waters

GC/MS/MS Analysis of Pesticides in Complex Food Matrices

For Complete  Confidence

- Hundreds of Pesticide Compounds in Use
 - Requirement for Surveillance Monitoring of Foodstuffs
 - Permitted Concentrations Regulated by Legislation
- Reporting Levels
 - Limits of Detection below the regulation are Required
 - target LOD of 0.01mg / kg (10 ppb)
- Need for the Targeting of Multiple Analytes per Analysis
 - Wide Range of Analyte Chemistries
 - polarity, thermal stability, volatility
 - Wide Range of Sample Matrices
- Common Sample Matrix Extraction and Minimal Sample Prep
 - Selectivity obtained from Chromatographic and Determinative Step
 - Analyte Specific Multiple Reaction Monitoring (MRM)
 - Triple Quadrupole Mass Spectrometry
 - Combined GC, chromatographic selectivity, and MS/MS is recommended

- Develop a rapid multi-residue pesticide screening method
 - 100 Residue Analytes in a Single Analysis
 - Reduce overall method cycle time, injection to injection
 - To increase efficiency and sample throughput
- Able to quantify and confirm the analyte at MRL tolerance levels
 - Include an additional MRM transition for the most commonly observed residues
 - Confirmatory method
 - Potentially reducing overall analysis time for a batch
 - TargetLynx

- Increasing selectivity and/or sensitivity
 - Sample Preparation (unlikely for multi-residue methods)
 - Chemical Ionization (narrows target range)
 - High Resolution MS
 - Magnetic Sector (expensive)
 - Dioxin/Furans, Brominated Flame Retardants, PCB
 - EPA approved methods
 - Time of Flight (limited dynamic range)
 - Tandem MS/MS
 - Ion Trap (limited dynamic range, complex matrices)
 - Tandem Quadrupole MS/MS

Waters

System Solution



Quattro micro GC/MS/MS

- 100 residue compounds in a single analysis
- Multiple compound classes
 - OPs, OCs, Pyrethroids
- Run time approximately <45 minutes
- Chromatographic and mass selectivity to compensate for non-selective sample preparation
- Less than 10 ppb sensitivity to minimize injection volumes of “dirty” matrix

- Column
 - DB5 MS, 30m x 0.25mm i.d., 0.25 μ m
- Injector
 - Cryo cooled PTV in solvent vent mode, 2 μ L injected
- Constant flow rate of 1.0 mL/min Helium
- Temperature program

0 min	50°C (2min)
6 min	150°C @ 25 °C/min
32 min	280°C @ 5°C/min
37 min	280°C


Dichlorvos
 cis-Mevinphos
 trans-Mevinphos
 Ethoprophos
 Chlorpropham
 Cadusafos
 Phorate
 α -HCH
 Hexachlorobenzene
 Dimethoate
 Carbofuran
 β -HCH
 Dimethipin
 Quintozene
 Lindane
 Diazinone
 Chlorothalonil
 Pentachloroaniline
 Phosphamidon
 Chlorpyriphos-methyl
 Vinclozolin
 Parathion-methyl
 Alachlor
 Carbaryl
 Heptachlor


Pirimiphos-methyl
 MPCPS
 Malathion
 Chlorpyriphos
 Aldrin
 Fenthion
 Parathion
 Fenpropimorf
 Triadimefon
 Oxychlordane
 Fipronil
 Heptachlor epoxide
 Chlorfenvinphos
 Phenthoate
 Quinalphos
 Triadimenol
 Methidathion
 α -Chlordane
 Tetrachlorvinphos
 Endosulfan α
 γ -Chlordane
 Prothiofos
 Profenofos
 p,p'-DDE


Dieldrin
 o,p'-TDE
 Cinerin I
 Endrin
 Nitrofen
 Endosulfan β
 Fensulfothion
 p,p'-TDE
 o,p'-DDT
 Ethion
 Jasmolin I
 Triazophos
 Carbophenothion
 Endosulfan sulphate
 p,p'-DDT
 Propargite
 Piperonylbutoxide
 Pyrethrin I
 Phosmet
 Bromopropylate
 Fenpropathrin
 Tetradifon
 Phosalone
 Phenothrin

Azinphos-methyl
 Cinerin-II
 λ -Cyhalothrin
 Azinphos-ethyl
 Acrinathrin
 Jasmolin II
 cis-Permethrin
 trans-Permethrin
 Prochloraz
 Cyfluthrin I
 Cyfluthrin II
 Cyfluthrin III
 Cyfluthrin IV
 Cypermethrin I
 Cypermethrin II
 Cypermethrin III
 Flucythrinate I
 Cypermethrin IV
 Etofenprox
 Flucythrinate II
 Fenvalerate I
 Fluvalinate
 Fenvalerate II
 Deltamethrin

- MRM functions arranged into 14 time windows
 - More flexible use of dwell times
 - Improved signal to noise (S/N) ratios for some analytes
 - Overlapping allows for small changes in retention time









2 μ l injected = 20 pg on column


Compound name: Dichlorvos
 Correlation coefficient: $r = 0.998799$, $r^2 = 0.997600$
 Calibration curve: $12.7295 * x + -2.86235$
 Response type: External Std, Area
 Curve type: Linear, Origin: Exclude, Weighting: 1/x, Axis trans: None




Compound name: Tetradifon
 Correlation coefficient: $r = 0.999367$, $r^2 = 0.998735$
 Calibration curve: $6.93326 * x + 21.1316$
 Response type: External Std, Area
 Curve type: Linear, Origin: Exclude, Weighting: 1/x, Axis trans: None




Compound name: Pirimiphos-methyl
 Correlation coefficient: $r = 0.999511$, $r^2 = 0.999022$
 Calibration curve: $20.2691 * x + -27.0132$
 Response type: External Std, Area
 Curve type: Linear, Origin: Exclude, Weighting: 1/x, Axis trans: None



Compound name: Deltamethrin
 Correlation coefficient: $r = 0.998918$, $r^2 = 0.997838$
 Calibration curve: $14.0112 * x + 125.819$
 Response type: External Std, Area
 Curve type: Linear, Origin: Exclude, Weighting: 1/x, Axis trans: None







Limit of Detection = S/N ratio \geq 3:1

- Sample matrix
 - Strawberry, Baby Food, Dried Tea, Tobacco, Spices, Ginko, Cannabis
- Homogenize using high speed blending
- Extraction with ethyl acetate (Turrax)
- Cleanup with dispersive SPE
 - Primary Secondary Amine followed by Graphitised Carbon Black
- Final extract
 - Fresh produce, baby food = 0.5 g/mL
 - Dried herbs, spices, tobacco = 0.1 g/mL







Pesticides Spiked
in Baby Food
Quattro micro GC
0.01 mg / kg
10 ppb




Herbal Tea = 0.2 mg matrix






Pesticides Spiked
in Herbal Tea
Quattro micro GC
0.3 ppm dry
0.03 ppm fresh

Name: TNC633
Description: Herbal Tea 15ppM spike
Date: 22-Apr-2004



Curry Spice Mix = 0.2 mg matrix on column


TIC Scan EI+





Spiked Pesticides in
Curry Spice Mix
Quattro micro GC
0.15 ppm

Name: TNC655
Description: Curr: 15 ppb/1 spike
Date: 22-Apr-2004



- High % of pesticides correctly identified & quantified
- Highly complex matrices
- Generic, limited sample preparation

Food Matrix	Pesticide Fortification	# Pesticides Spiked	# Pesticides Identified	% Detection
Baby Food	0.01 mg/kg	93	81	87
Baby Food	0.05 mg/kg	93	92	99
Herbal Tea	0.03 mg/kg	93	88	95
Herbal Mix	0.03 mg/kg	93	87	94
Curry	0.15 mg/kg	93	90	97
Massala	0.15 mg/kg	93	87	94
Tobacco	0.15 mg/kg	93	87	94
Ginko	0.02-4.0 mg/kg	52	44	85
Cannabis	0.02-4.0 mg/kg	52	41	79

Waters

A Rapid Multi-Residue
LC-MS/MS Method for the
Analysis of Pesticides

For Complete  Confidence

- Validated a multi-residue screening method in 2003, suitable for surveillance monitoring of fruits and vegetables using a Waters Alliance 2795 with Quattro Micro API MS
 - 81 Residue Compounds in a Single Analysis
 - Benzimidazoles, Carbamates, Organophosphorous etc.
 - All LODs below those required for monitoring purposes
 - Overall cycle time of the method was 40 min
 - Analyses performed from a variety of food matrices (avocado, tomato, raisin, lemon, flour)
- Method transferred to Waters Alliance 2795 with Waters micromass **Quattro Premier API MS**
 - Confirmatory method for the 20 most commonly observed residues
 - Overall cycle time of the method was 25 min

- Develop a rapid multi-residue pesticide screening method
 - 100 Residue Compounds in a Single Analysis
 - Reduce overall cycle time of the method
 - To increase efficiency and sample throughput
- Able to quantify and confirm at MRL tolerance levels
 - Include an additional transition for a selection of the most commonly observed residues
 - Confirmatory method
 - Potentially reducing overall analysis time for a batch

Waters

System Solution ACQUITY UPLC/MS/MS

Quattro Premier





- Column: BEH C₁₈, 2.1x100 mm, 1.7 μm
- Gradient: 20% to 90% MeOH containing 5 mM NH₄OAc
Linear for 8.5 mins, hold for 2 min, re-equilibrate
- Flow: 450 μL/min
- Col Temp: 40° C
- Inj Volume: 20 μL
- Waters Micromass Quattro Premier
 - Ionization mode: ES+ / ES-
 - Capillary voltage: 0.8 kV (ES+ and ES-)
 - Gas Flow: 800 L/hr
 - Source Temp: 120° C
 - Desolv'n Temp: 80° C
 - Cone Voltage: Analyte specific
 - MS/MS: operated in MRM mode
 - Collision Voltage: Analyte specific

5-Hydroxy-clethodim-sulfon	Clethodim-sulfon	Imazalil	Propoxur
6-chloro-4-hydroxy-3-phenyl-pyridazin	Clethodim-sulfoxid	Imidacloprid	Prosulfuron
Acephate	Cyprodinil	Indoxacarb	Pymetrozin
Aldicarb	Daminozid	loxynil	Pyridate
Aldicarb-sulfoxid	Demeton-S-methyl-sulfon	lprovalicarb	Pyrimethanil
Aldoxycarb	Desmedipham	Isoproturon	Quinmerac
Amidosulfuron	Difenzoquat methylsulfate	Isoxaflutole	Quizalofop-ethyl
Atrazin	Diflubenzuron	Linuron	Rimsulfuron
Azoxystrobin	Dimethoat	Metalaxyl	Spiroxamine
Bendiocarb	Diuron	Metamitron	Tebuconazol
Bensulfuron-methyl	Ethiofencarb	Methamidophos	Tebufenozide
Bromoxynil	Ethiofencarbsulfon	Methiocarb	Teflubenzuron
Butocarboxim	Ethiofencarbsulfoxid	Methomyl	Thiabendazol
Butoxycarboxim	Fenhexamid	Metolachlor	Thiacloprid
Butoxycarboxim-sulfoxid	Fenoxycarb	Metsulfuron-methyl	Thifensulfuron-methyl
Carbaryl	Fenpropimorph	Monocrotophos	Thiodicarb
Carbendazim	Flazasulfuron	Nicosulfuron	Thiofanox
carbofuran	Florasulam	Omethoat	Thiofanox-sulfon
Carbofuran-3-hydroxy	Fluazifop-P-butyl	Oxamyl	Thiofanox-sulfoxid
Chlorsulfuron	Fludioxonil	Oxydemeton-methyl	Triasulfuron
Cinosulfuron	Flufenoxuron	Phenmedipham	Triflumuron
Clethodim	Formetanate	Pirimicarb	Triflusulfuron-methyl
Clethodim-imin-sulfon	Furathiocarb	Primisulfuron methyl	Vamidothion
	Haloxypop-ethoxyethyl	Promecarb	


10 ppb matrix matched standard
acquityraisin201204_021

142MRRM010 Channels ES+
TIC
5.38e6








Quattro micro API, Alliance 2695
81 residues
40 min cycle time





Quattro Premier, Alliance 2695
100 residues
25 min cycle time



Quattro Premier, ACQUITY UPLC
100 residues
14.5 min cycle time








Butoxycarboxim Sulfoxide
and Aldicarb Sulfoxide
have the same MRM transition

m/z 207.1 > 89

And

Similar retention time: $\Delta = 12$ seconds



Requires Enhanced
Chromatographic Resolution
for Correct Analyte
Identification and Quantification


For accurate quantitation a minimum of 10 data points are needed across each peak.

Due to the

- Narrow peak widths
- Number and type of compounds to be analysed
 - Multi-analyte method
 - Inclusion of confirmatory transitions for some analytes
 - Some analytes ionise preferentially in negative ion mode


Significant demands are placed on the mass spectrometer


- The T-Wave is used to rapidly clear the cell between transitions, even when using a very short inter-channel delay (5ms) – minimising crosstalk.
- The T-Wave also ensures that the cell is completely refilled before acquisition commences allowing very short dwell times without loss of intensity.
- This enables quantification of narrower chromatographic peaks or quantification of more MRM transitions in the same time window.



“Travelling Wave Ion Propulsion in Collision Cells” K. Giles, S. Pringle, K. Worthington and R. Bateman—Presented at the 51st ASMS Conference, Montreal, Canada 2003. The travelling wave device described here is similar to that described by Kirchner in US Patent 5,206,506 (1993).

- MRM functions arranged into 26 retention time windows
 - More flexible use of dwell times
 - Positive and negative switching
 - Improved signal to noise (S/N) ratios






- Dimethoate and vamidothion determined under positive ion
- Bromoxynil determined under negative ion


Pesticide Residue	Retention Time (min)	Precursor Ion (m/z)	Product Ion (m/z)	Cone Voltage (V)	Collision Voltage (V)	Dwell Time (ms)	LOD (pg/uL)
Daminozid	0.50	161.1	143.1	18	12	200	0.01
Methamidophos	0.79	141.8	93.8	22	14	80	0.02
			124.9	22	13	80	
Acephate	0.89	184.1	143.0	16	8	40	0.04
Butoxycarboxim-sulfoxide	1.00	207.1	132.1	17	6	30	0.05
Omethoate	1.01	214.0	183.0	20	12	30	0.01
			154.9	20	15	30	
Aldicarb-sulfoxide	1.11	207.1	132.0	16	10	30	0.1
			89.0	16	14	30	
Butoxycarboxim	1.20	240.1	106.1	10	14	30	0.04
Aldoxycarb	1.26	240.1	86.0	15	20	30	0.005
Oxamyl	1.32	237.1	71.9	12	10	30	0.003
Propamocarb	1.36	189.1	102.0	25	17	30	0.01
			144.0	25	12	30	
Oxydemeton-methyl	1.49	247.0	169.0	20	13	10	0.002
Pymetrozin	1.57	218.0	105.0	25	17	10	0.02
6-chloro-4-hydroxy-3-phenyl-pyridazin	1.60	207.1	77.0	35	30	10	0.04
			104.0	35	21	10	
Methomyl	1.60	162.9	87.8	15	8	10	0.01
			105.9	15	10	10	
Demeton-S-methyl-sulfon	1.61	263.1	169.1	28	16	10	0.02
			121.2	28	16	10	
Quinmerac	1.69	222.0	141.0	22	33	10	0.008
Monocrotophos	1.78	224.0	126.9	20	15	10	0.005
Bendiocarb	1.78	224.1	109.0	18	18	10	0.01
			167.1	18	9	10	
Nicosulfuron	1.80	411.0	182.1	22	18	10	0.05
Amidosulfuron	1.84	370.0	261.2	18	14	10	0.02
Metsulfuron-methyl	2.00	382.0	167.0	22	15	10	0.02
Thifensulfuron-methyl	2.00	388.0	167.1	22	15	10	0.02
Ethiofencarbsulfon	2.04	275.1	107.1	10	20	10	0.006
Rimsulfuron	2.05	431.9	182.1	30	22	10	0.02
Ethiofencarbsulfoxide	2.13	242.1	107.0	18	18	10	0.003
Thiofanox-sulfoxide	2.14	252.1	104.0	10	12	10	0.3
Imidacloprid	2.14	256.1	209.2	22	16	10	0.02
			175.1	22	20	10	
Florasulam	2.28	360.1	129.0	30	20	10	0.09
5-Hydroxy-clethodim-sulfon	2.29	408.2	204.2	22	16	10	0.1
Thiofanox-sulfon	2.32	268.1	76.0	10	10	10	0.02
Clethodim-imin-sulfon	2.35	302.2	98.1	35	30	10	0.04
Metamitron	2.37	203.0	175.1	28	16	10	0.02
Cinosulfuron	2.42	414.1	183.1	25	18	10	0.05
Chlorsulfuron	2.43	358.1	141.1	25	16	10	0.08
			167.1	25	16	10	
Bromoxynil*	2.45	273.9	78.9	40	25	30	0.2

LOD at 3x S/N
in the sub-ppb range
0.003 to 2 ppb



- Sample: Californian Sun-dried Seedless Raisins (Thompson variety)
- Chopped
- 5g Aliquot homogenised
- 9 mL water added
- After 10 min, 20 mL of methanol added
- 6mL of extract was mixed with 2mL of sodium chloride solution
- 5mL aliquot of the mixture was transferred to a ChemElut column containing 5 mL of diatomaceous earth
- After 5min the column was eluted with 16 mL DCM
- Elute evaporated
- Reconstituted in methanol (250 μ L)
- Addition of 1 mL water and filtering into glass vial

- Blank matrix: organically grown Sun-dried Seedless Raisins (Thompson variety)
- Same sample preparation, extraction and clean-up procedure as for sample
- All analytes spiked into the matrix at a range of concentrations equivalent to 1, 2, 5, 10, 20 ppb





Compound name: Carbendazim
Correlation coefficient: $r = 0.998134$, $r^2 = 0.996272$
Calibration curve: $5128.64 * x + 66.1767$
Response type: External Std, Area
Curve type: Linear, Origin: Exclude, Weighting: 1/x, Axis trans: None




1 ppb Chromatogram

acquityraisin201204_013
1 ppb matrix matched standard



F18:MRM of 11 channels, ES+
226.2 > 93.1
1.952e+004

Compound name: Cyprodinil
Correlation coefficient: $r = 0.997715$, $r^2 = 0.995435$
Calibration curve: $807.68 * x + 136.057$
Response type: External Std, Area
Curve type: Linear, Origin: Exclude, Weighting: 1/x, Axis trans: None



acquityraisin201204_013
1 ppb matrix matched standard


F18:MRM of 11 channels, ES+
226.2 > 108.1
1.532e+004



1 ppb Chromatogram

acquityraisin201204_002
1 ppb matrix matched standard



F21:MRM of 2 channels, ES-
379 > 339.1
7.478e+002



Compound name: Teflubenzuron
Correlation coefficient: $r = 0.992141$, $r^2 = 0.984344$
Calibration curve: $19.9838 * x + 12.801$
Response type: External Std, Area
Curve type: Linear, Origin: Exclude, Weighting: 1/x, Axis trans: None

acquityraisin201204_002
1 ppb matrix matched standard





F21:MRM of 2 channels, ES-
379 > 196
2.412e+002



1 ppb chromatogram

Imidacloprid 2.2 ppb


Tebufenozide 1.7 ppb






- Comprehensive Tool for Data Processing
 - Quantitative Data Generation and Evaluation
 - QA/QC System and Method suitability
- When injecting a series of samples, the analyte concentrations are reported as follows:
 - Analytes which are below the reporting limit are not reported
 - user defined values, i.e. MRL, LOQ, S/N, etc
 - Analytes between the reporting limit and the maximum allowed MRL are reported
 - Analytes **above** the MRL are reported in **red Bold**
 - The method QA/QC parameters are flagged in red if exceeds user defined limits, and the explanation appears as tool tip
 - RT, S/N, r^2 , plate count, calib checks, spiked duplicates, etc

- Carbaryl detected as positive in a raisin matrix blank using the pesticide multi-residue screening method
- 6 additional analytes also detected in matrix blank
 - Confirmatory method created, monitoring 3 MRM transitions per compound
 - TargetLynx method created by exporting the compounds from the screening method
 - Matrix blank, matrix spiked and 'real' raisin samples acquired and processed





- Multi-residue method for surveillance monitoring of 100 pesticide residues in food and natural products, complex matrices
 - GC/MS/MS and UPLC/MS/MS
- Chromatographic and MS/MS Selectivity to allow generic sample preparation
 - Enhanced chromatographic resolution and shorter analysis times
 - Ability to group MRM functions into time windows, enabling the incorporation of confirmatory MRM traces
 - Ability to switch rapidly between MRM channels and between positive and negative ionization modes
- Good sensitivity (<5 ppb generally) for pesticide residues
 - Sufficient to limit the amount of matrix introduced into the chromatographic inlet

- Waters Corporation, Manchester, UK
 - David Douce
 - Gordon Kearney
 - Peter Hancock
 - Catherine Ryan
- TNO, Zeist, The Netherlands
 - Hans Mol
- FIRA, Berlin, Germany
 - Lutz Alder
 - Jeannette Klein