

GC/MS Mass Spectra Library

Polymer Additives Library

The analysis of additives contained in polymer materials is essential for managing and improving the quality of these materials, and for complying with regulations dealing with chemical substances.

The Polymer Additives Library is a GCMS mass spectral library containing information on a wide range of additives utilized in polymer materials. In addition to mass spectra and retention indices, it contains information on the classification of additives to enable detailed analysis of additives.



GCMS-QP2020 NX

Contains Approximately 4,900 Mass Spectra

In addition to the existing library*, which contains 4,804 mass spectra for 494 types of polymer material additives and GC/MS pyrolysates from the pyrolysis GC/MS analysis of additives, a library has been added containing the mass spectra for 65 compounds often targeted for analysis, selected based on their usage in the commercial market, and chemical substance regulatory information. In total, this product contains 4,869 mass spectra, which provides strong support for the analysis of polymer additives.

*ADD-MS16B F-Search Additives Library from Frontier Laboratories

Ready-to-Use Methods

This library provides GCMS and PY methods; both contain all necessary conditions and measurement parameters. Anyone can easily start analysis with optimal analysis conditions.

Library Specifications

Registered compounds*: 4,804+65

Registered information: Mass spectrum, retention index for each analytical condition, compound name, molecular weight, compositional formula, structural formula, classification of additives

Applicable models: GCMS-QP series + GCMSsolution™ ver. 2.6 or later

GCMS-TQ™ series + GCMSsolution ver. 4.0 or later (This library does not include MS/MS spectra.)

Pyrolyzer (Pyrolysis analysis system): PY-2020D, PY-2020iD, EGA/PY-3030D

Includes Retention Indices and Classification Information of Additives

Retention indices are registered for all compounds. Compounds can be identified with a high degree of accuracy by reducing the library search results using retention indices. In addition, information on the classification of additives is included, so even without detailed knowledge of additives, users can see which additive is associated with a compound found in the library.

Can Be Used with a Variety of GC/MS Systems

The library can be used for a variety of GC/MS applications, including pyrolysis GC/MS, which is widely used for the analysis of additives in polymer materials, and liquid sample injection GC/MS. It supports a wide range of additive analysis for customers.

*: This consists of two types of libraries. The library with 4,804 compounds consists of the ADD-MS16B F-Search Additives Library (Ver. 16B.02) developed by Frontier Laboratories, which has been converted to the Shimadzu mass spectral library format, and prepared for use with GCMSsolution. The library with 65 compounds registered is an additive library developed exclusively by Shimadzu. It contains additives often targeted for analysis, selected based on their usage in the commercial market, and chemical substance regulatory information. Using these two libraries in combination enables the analysis of an even wider range of additives.

Polymer Additives Library

Filtering with Retention Index

Multiple compounds with similar mass spectra are listed as candidates when performing a library search using only the mass spectrum. Filtering with the retention index sorts the candidates by retention index, thereby providing highly accurate identification results.

Results of similarity search using mass spectrum

Hit#	Similarity	Register	Ret. Index	Compound Name	Mol Wt	Formula	Library
1	85	<input type="checkbox"/>	1732	Ethyl ethyl benzoate [Original Additive] B	232	C18H22O2	LAB_ADD1
2	84	<input type="checkbox"/>	2645	Undecyl benzoate [Original Additive] C	276	C18H28O2	LAB_ADD1
3	79	<input type="checkbox"/>	1795	Octyl benzoate [Original Additive] T	254	C18H26O2	LAB_ADD1
4	79	<input type="checkbox"/>	2103	Dodecyl benzoate [Original Additive] D	330	C22H34O2	LAB_ADD1
5	66	<input type="checkbox"/>	792	1-Octyl-1-Propyl-Additive-Additive	172	C18H34	LAB_ADD1
6	55	<input type="checkbox"/>	791	Isomer of C18H34 [Original Additive] A	172	C18H34	LAB_ADD1




Results sorted using retention index filtering

Hit#	Similarity	Register	Ret. Index	Compound Name	Mol Wt	Formula	Library
1	85	<input type="checkbox"/>	1732	Ethyl ethyl benzoate [Original Additive] B	232	C18H22O2	LAB_ADD1

Confirming the Additive Classification Information

The additive classification information (such as plasticizers and flame retardants) registered in the library plays a role in confirming the type of additive associated with the compounds included in search results.



GCMSsolution, GCMS-TQ and Py-Screener are trademarks of Shimadzu Corporation or its affiliated companies in Japan and/or other countries.

Screening System for Phthalate Esters

Can be Used in Combination with Py-Screener™

In addition to screening for phthalate esters, using this product in combination with Py-Screener enables the analysis of a wide range of other additives.



Example of System Configuration*

Pyrolysis Analysis System



Liquid Injection Analysis System



*: The accuracy of the retention indexes differs depending on the sample introduction method. Some of the pyrolysates may only be seen using a pyrolysis GC/MS analysis.



Shimadzu Corporation

www.shimadzu.com/an/

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

Company names, products/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol "TM" or "®".

Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The contents of this publication are provided to you "as is" without warranty of any kind, and are subject to change without notice.

Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication.