



# Agilent Oil Analyzer: customizing analysis methods

## Application Note

### Author

Alexander Rzhevskii &  
Mustafa Kansiz  
Agilent Technologies, Inc.

### Introduction

Traditionally, the analysis of used oils has been conducted by physical and wet chemical methods. FTIR spectroscopy has become a routinely used technique to analyze used oils, providing the following major advantages<sup>1</sup>:

- Ability to simultaneously determine several parameters from a single experiment
- Increase in speed of analysis
- More cost effective than traditional techniques
- Mobility and portability allowing remote on-site analysis

The Agilent FTIR Oil Analyzer is designed to meet the requirements of the US Department of Defense Joint Oil Analysis Program (JOAP)<sup>2</sup> for use in their condition monitoring program as well as commercial applications. It is optimized for monitoring relative changes in various indicators of oil conditions (oil failure symptoms) using a standardized protocol developed by the Joint Oil Analysis Program Technical Support Center (JOAP-TSC). This protocol sets the data extraction algorithm for several types of petroleum and synthetic-based lubricants and hydraulic fluids, and eliminates the need for reference samples as spectral subtraction is no longer required.

The Agilent Oil Analyzer software allows users to readily customize existing methods as well as create new methods to measure other parameters and properties of lubricants defined by the user. The methods can be easily adjusted for performing analysis of samples where spectral subtraction is required.



**Agilent Technologies**

This application note describes the tools available with the Agilent FTIR Oil Analyzer and procedures that a user should follow to customize analysis methods, while reinforcing the importance of reliable calibration in quantitative spectral analysis.

## Analysis methods

The sampling and analyzing procedures available in the Agilent FTIR Oil Analyzer conform to the ASTM E 2412-04 "Standard practice for condition monitoring of used lubricants by trending analysis using Fourier Transform Infrared (FTIR) Spectrometry"<sup>3</sup>. These methods provide a generalized protocol for condition monitoring of contaminants and breakdown products in used lubricants including water, ethylene glycol, fuels, incorrect oil, soot, oxidation, nitration and sulfonation. The methods are based on calculating trends and distributions from mid-IR absorption measurements, and encompass both direct and differential (spectral subtraction) trend analysis approaches.

The Agilent Oil Analyzer software is configured to run twelve predefined analysis methods that correspond to different classes of lubricating oils or hydraulic fluids, and their applications with differing limits. The methods are:

- Aircraft hydraulic (Mil-H-83282)
- Aircraft hydraulic (Mil-H-83282\_350 ppm limit for water)
- Dextron transmission fluid
- Engine crankcase (Diesel\_gasoline\_natural gas)
- Fire retardant hydraulic (Mil-H-46170)
- Gas turbine or Helo Gbx (Mil-L-23699)
- Ground equipment hydraulic (Mil-L-2104\_10W)
- Ground equipment synthetic hydraulic (Mil-H-5606)

- Marine diesel crankcase (Mil-L-9000)
- Conostan IR OTS fluid
- Steam turbine (Mil-L-17331)
- Generic or undetermined (Unknown lubricant type)

Each of the methods measures numerical indicators (parameters) that are related to the oil's condition. The software then generates a report that contains thirteen measurement parameters, as listed below:

- Water in EP fluids
- Antioxidant reading
- Ester breakdown
- Water in petroleum
- Soot value
- Oxidation by-products
- Nitration by-products
- Antiwear reading
- Gasoline dilution
- Diesel/JP8 dilution
- Sulfate by-products
- Ethylene glycol
- Other fluid contamination

Additionally, a separate procedure for predicting Total Base Number (TBN) is available and can be integrated into existing methods.

The parameters are reported in the units of spectral absorbance (peak areas or heights) rather than in physical concentrations, such as ppm, wt.% or mg of KOH. Figure 1 shows an example of a typical standard Oil Analysis report.

<b>Oil Analysis</b>	
Date: 7/27/2005	
Time: 05:09 PM	
Software Version: 4.2.8	
Sample ID: Preview	
TEC: XXXX	
Component Model Number: XXXXXX	
Component Serial Number: XXXXXX	
End Item: XXXXX	
End Item Serial Number: XXXXXX	
Time Since Fluid Change: 0	
Total Component Hours: 0	
Matched Spectra Name:	
Matched Spectra Comment:	
Lube Analysis Type: TEST	
Water in EP Additive Fluids... (N/A).....	1.
Antioxidant Reading.....	1.
Ester Breakdown I... (N/A).....	0.
Water Petroleum Lube... (Normal 10 to 40)...65 = 2000 ppm...	264.
Soot Value... (Normal 0).....	0.
Oxidation By-Products... (Normal 10 to 12).....	514.
Nitration By-Products.....	965.
Antiwear Reading... (Normal 8 to 12).....	1.
Gasoline Dilution... (N/A).....	1.
Diesel/JP8 Dilution... (N/A).....	1.
Sulfate By-Products... (Normal 10 to 14).....	736.
Ethylene Glycol (Antifreeze)... (N/A).....	487.
Other Fluid Contamination... (Normal 100).....	679.
<b>Notes and Warnings</b>	

Figure 1. Typical standard Oil Analysis report

## Calibration

All analysis methods in the Agilent FTIR Oil Analyzer consist of a set of calibration models (procedures) in the form of corresponding files with an indication of the calibration model's type (univariate, or multivariate, or a combination). The analysis method may be composed of one or several calibration files.

The construction of calibration models in quantitative spectral analysis is a two-step procedure: calibration and validation. In the calibration step, indirect instrumental measurements (spectra) are obtained from standard samples in which the value of the parameter of interest has been determined by a standard reference method (an accurate direct measurement method). The set of spectra and results from the reference method, referred to as the calibration set or training set, is used to construct a model that relates parameter values to the spectra. Before the calibration model is accepted and used for prediction, it should be validated by a set of independent (not used in the calibration set) samples of known parameter concentrations (validation set). If parameters from the validation set fall within acceptable accuracy limits using the model derived

from the calibration set, an acceptable model has been constructed that can be used to predict for new "unknown" samples.

To build a univariate calibration model, it is necessary to specify a single measurement from a spectrum, such as peak area or height that demonstrates the most distinctive spectral response for the parameter of interest. The univariate calibration and prediction procedures are available as a standard part of Resolutions/Resolutions Pro software and are defined as a simple quantitative analysis. The analysis is described in detail in the Resolutions online help and the corresponding system reference manuals for previous software versions (Win-IR Pro and Merlin). The user must generate a quantitative calibration document and save it as \*.BSQ file using Resolutions/Resolutions Pro (Win-IR Pro or Merlin) software.

Where spectral responses attributed to different parameters overlap and the selective spectral measurements for the parameter of interest is very difficult, univariate models may not be reliable. Multivariate methods such as Principal Component Regression (PCR) and Partial Least Squares (PLS) allow multiple responses at the selected wavenumbers to be used. These methods are better suited to extracting spectral information where bands overlap and it is difficult to discern the relevant spectral regions attributable to a particular parameter. The main advantage of multivariate methods is the ability to calibrate for a parameter of interest when it correlates in a complicated (non-specific) way with multiple spectral regions, while minimizing background matrix interferences in the lubricants.

The Agilent Oil Analysis software allows multivariate calibration models created with the use of third party software to be incorporated in analysis methods. The PLSplus IQ package available as an additional application in the Galactic GRAMS/AI (GRAMS/32) software suite must be used. The "PLSplus IQ User's guide" gives step-by-step instructions on how to construct and validate a multivariate calibration model

as well as theory of advanced statistical analysis in spectroscopic quantitative analysis. The user must build an accurate calibration model and save it into a \*.CAL file using PLSplus IQ.

The validity of empirically-built calibration models depends heavily on how well the standard samples (calibration set) represents the unknown samples to be analyzed (prediction set). In all cases, the selection of standard samples to be used for calibration must adequately cover the expected range of measurement parameters in the prediction set. This means that the expected extreme values for each parameter of interest in unknown samples must be included in the calibration set, as extrapolation outside the calibrated value range can be unreliable. It is important to ensure that any phenomena that influence the spectral measurements (e.g., not only the total amount of soot but its particle size distribution) also vary in the calibration set over ranges that span the levels of the phenomena occurring in the prediction set. It is also very important to minimize the errors in the standard sample parameters that are used to construct the empirical calibration model, as any calibration model can only be as accurate as the reference measurements from which it was constructed.

Many conditions can affect the results obtained from FTIR lubricant monitoring such as lubricant type, engine type, operational conditions, environmental conditions, etc. When the conditions are changed significantly, new calibration models and methods may be required to ensure accurate prediction of oil properties. For instance, new calibrations may be required when a new oil type with a different base stock and additive chemistries comes for the analysis.

Care must be taken when measuring overall oil quality parameters such as Total Acid Number (TAN) and Total Base Number (TBN) using FTIR spectroscopy. The secondary formation of acidic products in lubricants is characterized by TAN or indirectly by TBN, which assesses the consumption

of basic reserve additives in the oil. While the various acids or bases present in a lubricant could, in principle, be individually quantified based on their characteristic absorption bands, no unique absorption bands can be directly related to TAN or TBN. Thus, only indirect FTIR spectroscopic methods for TAN and TBN have been standardized to date. In addition, there is a large discrepancy in new lubricant TAN values, from less than 0.1 mg KOH/g for R&O type oils to 9 or higher for some synthetic oils in industrial applications. On the other hand, the incremental decrease in TBN used to indicate that a product is failing, varies in broad ranges: some oils may have a new TBN value of 12, but rapidly decrease to a value of 3, whereas other synthetic oils may have the beginning TBN of 40.

A calibration model for TBN is currently available in Agilent Oil Analyzer. The calibration is intended for prediction of the values in gasoline and diesel engine oils having typical baseline numbers not higher than 12 mg KOH/g.

Note that in many individual cases, in order to estimate TAN and TBN satisfactorily the user needs to construct a multivariate calibration model that would cover the higher range of values as well as take into account any other factors that could influence the accuracy and the reproducibility of spectral measurements.

## Method editor

Once the univariate or multivariate calibration models are built, the corresponding \*.BSQ or \*.CAL files must be moved or copied into the directory

C [Local Disk]:\ Program

Files\Varian\Resolutions\Oil Analyzer\Methods.

This is the storage location for the available calibration and method files. Then, log in as Administrator to the Agilent Oil Analysis software and enter the Method editor. Follow the Chapter 11 "Method Editor" in "Agilent Oil Analyzer operational manual" to incorporate the calibrations to an existing method or to develop a new method.

Note that spectral subtraction is available in the Agilent Oil Analyzer but was not utilized in JOAP protocol. It is not considered to be practical in view of the deployability aspect of many JOAP laboratories and that the required sample volume would increase because of the necessity of new oil samples to act as references. In order to apply the spectral subtraction procedure, the user needs to select "Use spectral subtraction" option in the Sampling method group in the General option dialog and edit the relevant analysis method, by clearing the "Zero less than Zero" check box in all the associated calibration models. Refer to Chapter 4 "General Options–Setup" and Chapter 11 "Method Editor" of "Agilent Oil Analyzer operational manual" for more information.

## Conclusion

FTIR spectroscopy has been gaining increased acceptance as a method of choice for used oil analysis. Designed and optimized as a complete system for predictive maintenance programs, according to JOAP standards, the Agilent FTIR Oil Analyzer combines specific capabilities with the flexibility to be successfully used in any oil analysis laboratory.

The Oil Analyzer software allows new and improved analysis methods to be built and ensures that new types of lubricating oils and fluids used in a variety of different machinery are timely and reliably monitored and tested.

The software allows the user to include PCR/PLS methods to measure oil parameters and convert the units of spectral absorbance into physical results (ppm, wt.%, cSt, mg KOH/g oil, etc.) applying spectral subtraction if needed.

## References

- <sup>1</sup> Larry A. Toms, "Machinery Oil Analysis. Methods, Automation & Benefits", 2<sup>nd</sup> ed., Coastal Skills Training, Virginia Beach, VA, 1998.
- <sup>2</sup> Allison M. Toms, "FTIR for the Joint Oil Analysis Program", in Proc. 1994 Joint Oil Analysis Program International Condition Monitoring Conference, Squalls, M., ed., JOAP-TSC, Pensacola, FL (1994), pp.387-419.
- <sup>3</sup> Available from [www.astm.org](http://www.astm.org)

[www.agilent.com/chem](http://www.agilent.com/chem)

© Agilent Technologies, Inc., 2005, 2011  
Published March, 2011  
Publication Number ftir127

