

# Application Data Sheet

## No.72

### GC-MS

Gas Chromatograph Mass Spectrometer

## Scan/MRM Analysis of Residual Pesticides in Foods Using GC-MS/MS (3)

The GCMS-TQ8030 is a triple quadrupole GC-MS/MS system equipped with scan/MRM mode to allow simultaneous scan and MRM data measurements. This data sheet introduces the results of an investigation using the scan/MRM mode, where target pesticides were quantitatively determined using the MRM data, and concentrations of the untargeted pesticides were estimated by applying the scan data to the Compound Composer Database Software Ver.2.

### Experimental

For the evaluation, analytical standards (0.001 mg/L to 0.1 mg/L) were used, as well as samples (0.01 mg/L) created by pretreating paprika with the QuEChERS method, and then adding pesticides to the obtained solution. The pesticides specified as targets and their transitions were those recommended in the results of the validity evaluations by the European Reference Laboratory<sup>1)</sup>. The analysis conditions are shown in Table 1.

Table 1 Analytical Conditions

GC-MS	:GCMS-TQ8030	[MS]	
Column	:Rxi-5Sil MS (30 m length, 0.25 mm I.D., df=0.25 μm)	Interface Temp.	:300 °C
Glass Liner	:Sky Liner, Splitless Single Taper Gooseneck w/Wool (Restek Corporation, catalog # 567366)	Ion source Temp.	:200 °C
[GC]		Data Acquisition Mode	:Scan/MRM
Injection Temp.	:250 °C	Event Time	:0.1 sec (Scan), 0.3 sec (MRM)
Column Oven Temp.	:40 °C (2 min) → (8 °C /min) → 310 °C (5 min)	Scan Mass Range	:m/z 50 – 500
Injection Mode	:Splitless	Scan Speed	:5,000 u/sec
Flow Control Mode	:Linear velocity (40.0 cm/sec)		
Injection Volume	:1 μL		

### MRM Monitoring m/z

Compound Name	Quantitative Transition		Qualitative Transition		Compound Name	Quantitative Transition		Qualitative Transition					
	Precursor>Product	CE (V)	Precursor>Product	CE (V)		Precursor>Product	CE (V)	Precursor>Product	CE (V)				
Diphenylamine	169.10>77.00	26	169.10>115.10	30	Bupropion	172.10>57.10	18	105.10>104.10	4				
Ethoprophos	200.00>157.90	6	200.00>114.00	14	200.00>97.00	26	273.10>193.20	8	273.10>108.00	18			
Chlorpropham	213.10>171.10	6	213.10>127.10	18	beta-Endosulfan	240.90>205.90	14	238.90>203.90	14				
Trifluralin	306.10>264.00	8	264.10>206.10	8	264.10>160.10	18	Oxadixyl	163.10>132.10	10	163.10>117.10	24		
Dicloran	206.00>176.00	12	206.00>124.00	26	176.00>148.00	12	Ethion	231.00>174.90	14	231.00>128.90	26		
Propyzamide	172.90>144.90	16	172.90>109.00	26	Triazophos	161.10>134.10	8	161.10>106.10	14				
Chlorothalonil	265.90>230.90	14	265.90>167.90	24	263.90>167.90	24	Endosulfan sulfate	386.90>252.90	10	386.90>216.90	26		
Diazinon	304.10>179.10	12	179.20>137.20	18	Propiconazole-1	259.10>190.90	8	259.10>172.90	18	259.10>69.10	12		
Pyrimethanil	199.10>184.10	14	199.10>158.10	14	Propiconazole-2	259.10>190.90	8	259.10>172.90	18	259.10>69.10	12		
Tefluthrin	197.10>141.10	26	177.10>127.10	32	Tebuconazole	252.10>127.00	24	250.10>125.10	24				
Pirimicarb	238.20>166.10	10	166.10>96.00	14	Iprodione	314.10>244.90	12	314.10>56.10	24				
Chlorpyrifos-methyl	285.90>270.90	12	285.90>93.00	22	Bromopropylate	340.90>184.90	18	182.90>154.90	16				
Vinclozolin	212.10>172.00	14	212.10>144.90	26	212.10>109.00	30	Bifenthrin	181.10>166.10	16	181.10>165.10	22	181.10>153.10	10
Parathion-methyl	263.10>109.00	18	263.10>81.00	26	Fenpropathrin	265.10>210.10	12	181.10>152.10	24	181.10>127.10	26		
Tolclofos-methyl	265.00>249.90	12	265.00>93.00	24	Fenazaquin	160.20>145.10	8	145.20>115.10	24	145.20>91.10	24		
Metalaxyl	206.20>162.10	8	206.20>132.10	18	Tebuconazole	333.20>276.10	8	333.20>171.00	22				
Fenitrothion	277.10>125.00	18	277.10>109.00	18	Tetradifon	355.90>158.90	12	353.90>159.00	12	228.90>200.90	14		
Pirimiphos-methyl	305.10>290.10	12	290.10>125.00	24	Phosalone	182.00>138.00	8	182.00>111.00	18	182.00>102.10	18		
Dichlofuanid	332.00>167.10	6	224.00>123.00	12	Pyriproxyfen	136.10>96.00	12	136.10>78.00	24				
Malathion	173.10>117.00	12	173.10>99.00	18	Cyhalothrin	181.10>152.10	24	163.10>127.00	14	163.10>91.00	22		
Chlorpyrifos	196.90>168.90	14	196.90>107.00	26	Fenarimol	251.00>139.00	18	139.10>111.00	16				
Fenthion	278.10>125.00	22	278.10>109.00	18	Acrinathrin	289.10>93.10	12	181.10>152.10	24	208.10>181.10	8		
Parathion	291.10>109.00	14	291.10>81.00	26	Permethrin-1	183.10>168.10	12	183.10>153.10	18	183.10>115.10	24		
Tetraconazole	336.10>218.00	18	336.10>204.00	26	Pyridaben	147.20>132.10	14	147.20>117.10	22				
Pendimethalin	252.20>162.10	12	252.20>161.10	12	Permethrin-2	183.10>168.10	12	183.10>153.10	18	183.10>115.10	24		
Cyprodinil	225.20>224.10	6	224.20>208.10	18	Cyfluthrin-1	206.10>151.20	24	163.10>127.10	6	163.10>91.00	14		
(E)-Chlorfenvinphos	323.10>266.90	14	267.00>159.00	18	Cyfluthrin-2	206.10>151.20	24	163.10>127.10	6	163.10>91.00	14		
Tolyfluanid	137.10>91.00	18	137.10>65.00	26	Cyfluthrin-3	206.10>151.20	24	163.10>127.10	6	163.10>91.00	14		
Fipronil	367.00>227.90	26	367.00>212.90	26	Cyfluthrin-4	206.10>151.20	24	163.10>127.10	6	163.10>91.00	14		
Captan	79.00>77.00	8	79.00>51.00	22	Cypermethrin-1	181.10>152.10	24	163.10>127.10	6	163.10>91.00	14		
(Z)-Chlorfenvinphos	323.10>266.90	14	267.00>159.00	18	Cypermethrin-2	181.10>152.10	24	163.10>127.10	6	163.10>91.00	14		
Phenthoate	274.10>125.00	18	274.10>121.10	12	Cypermethrin-3	181.10>152.10	24	163.10>127.10	6	163.10>91.00	14		
Folpet	147.10>103.10	10	147.10>76.00	26	Cypermethrin-4	181.10>152.10	24	163.10>127.10	6	163.10>91.00	14		
Procymidone	283.10>96.10	12	283.10>67.10	24	Ethofenprox	163.20>135.00	10	163.20>107.10	18				
Methidathion	145.10>85.00	8	145.10>58.00	18	Fenvalerate-1	125.10>99.00	22	125.10>89.00	22				
alpha-Endosulfan	240.90>205.90	14	238.90>203.90	16	tau-Fluvalinate-1	250.10>200.10	16	250.10>55.00	18				
Mepanipyrim	222.20>220.10	8	222.20>193.10	26	Fenvalerate-2	125.10>99.00	22	125.10>89.00	22				
Profenofos	337.10>266.80	16	207.90>63.00	26	tau-Fluvalinate-2	250.10>200.10	16	250.10>55.00	18				
Myclobutanil	179.10>152.00	8	179.10>125.00	16	Deltamethrin-1	252.90>93.10	18	181.10>152.10	24				
Flusilazole	233.10>165.10	18	233.10>152.10	18	Deltamethrin-2	252.90>93.10	18	181.10>152.10	24				



## Results

An example of the results of the analysis of the analytical standards (0.001 mg/L to 0.1 mg/L) and the pesticide-spiked samples (0.01 mg/L) in scan/MRM mode are shown in Fig. 1. As with Procymidone, shown in Fig. 1, strict quantification of the targeted pesticides could be performed by creating calibration curves from the MRM data. Furthermore, since the scan data was sampled simultaneously, the pesticides could be confirmed from the mass spectrum.

For the untargeted pesticides, data analysis was performed utilizing the Compound Composer Database Software (P/N: 225-13106-92). The simultaneous analysis database software contains information (mass spectra, retention times, and calibration curves) on more than 450 pesticides. As a result, it is possible to identify pesticides from their estimated retention times and mass spectra without using analytical standards, and then calculate semi-quantitative values from the calibration curves. Also in this investigation, it was possible to confirm the detection and semi-quantification of untargeted pesticides such as Quinoxifen.

### Analysis Results for Target Compounds


Procymidone: Quantitative value 0.0098 mg/L



Mass spectral confirmation is possible with the scan data.  
(Upper: Sampled mass spectrum; Lower: Mass spectrum in the library)


Quantitative determination by creating calibration curves from the MRM data

Measurement Results for Pesticide-Spiked Samples




### Screening Results for Other Pesticides

Quinoxifen: Quantitative value 0.0134 mg/L



Compound identification and semi-quantitative calculations can be performed utilizing the information contained in the simultaneous analysis database software.



Reference

1) EURL-FV Multiresidue Method using QuEChERS followed by GC-QqQ/MS/MS and LC-QqQ/MS/MS for Fruits and Vegetables (European Reference Laboratory, 2010-M1)

Fig. 1 Scan/MRM Analysis Results

First Edition: Jan, 2013