

Separation of Fatty Acid Methyl Esters (FAME) on an Agilent J&W Select CP-Sil 88 for FAME GC Column

Application Note

Authors

Frau Thomae and
Frau Dr. Schwendig-Radke
Office for Consumer Protection,
Mettmann

Introduction

The routine, detailed separation of fatty acid methyl esters (FAME) requires high polarity liquid phases, which will differentiate between the multiple FAME isomers. The CP-Sil 88 is among the GC columns frequently used for the profiling of complex FAME mixtures. It is based on a stabilized, highly substituted cyanopropyl siloxane phase and due to its highly polar properties is able to effectively separate on small structural differences of many positional FAME isomers. This application note shows the routine separation of fatty acid methyl esters (FAME) using a Select CP-Sil 88 column and GC-FID resulting in excellent selectivity; all 37 components were resolved in one run.



Agilent Technologies

Conditions

Technique: GC-FID
Column: Select CP-Sil 88 for FAME,
100 m x 0.25 mm df=0.2 µm
(part number CP7489)
Oven: 80 °C, 4 °C/min to 220 °C, 5
min, 4 °C/min to 240 °C, 10
min
Carrier Gas: Helium, constant flow 1.0
mL/min
Injection: Splitflow 20 mL/min,
Temperature = 250 °C
Injection Volume: 0.5 µL
Detection: FID, 270 °C
Sample: 37 component FAME mix
Sample Solvent: Dichloromethane
Sample Conc: 2-4% (w/w)

Peak Identification

1.	C4:0	20.	C18:2 cis 9,12
2.	C6:0	21.	C20:0
3.	C8:0	22.	C18:3 cis 6,9,12 gamma
4.	C10:0	23.	C20:1 cis 11
5.	C11:0	24.	C18:3 cis 9,12,15 alpha
6.	C12:0	25.	C21:0
7.	C13:0	26.	C20:2 cis 11,14
8.	C14:0	27.	C22:0
9.	C14:1	28.	C20:3 cis 8,11,14
10.	C15:0	29.	C22:1 cis 13
11.	C15:1	30.	C20:3 cis 11,14,17
12.	C16:0	31.	C20:4 cis 5,8,11,14
13.	C16:1	32.	C23:0
14.	C17:0	33.	C22:2 cis 13,16
15.	C17:1	34.	C24:0
16.	C18:0	35.	C20:5 cis 5,8,11,14,17 EPA
17.	C18:1 trans 9	36.	C24:1 cis 15
18.	C18:1 cis 9	37.	C22:6 cis 4,7,10,13,16,19 DHA
19.	C18:2, trans 9, 12		




Figure 1a. Separation of a 37 component FAME mix (peaks 1-14)

Peak Identification

1. C4:0	14. C17:0	27. C22:0
2. C6:0	15. C17:1	28. C20:3 cis 8,11,14
3. C8:0	16. C18:0	29. C22:1 cis 13
4. C10:0	17. C18:1 trans 9	30. C20:3 cis 11,14,17
5. C11:0	18. C18:1 cis 9	31. C20:4 cis 5,8,11,14
6. C12:0	19. C18:2, trans 9, 12	32. C23:0
7. C13:0	20. C18:2 cis 9,12	33. C22:2 cis 13,16
8. C14:0	21. C20:0	34. C24:0
9. C14:1	22. C18:3 cis 6,9,12 gamma	35. C20:5 cis 5,8,11,14,17 EPA
10. C15:0	23. C20:1 cis 11	36. C24:1 cis 15
11. C15:1	24. C18:3 cis 9,12,15 alpha	37. C22:6 cis 4,7,10,13,16,19 DHA
12. C16:0	25. C21:0	
13. C16:1	26. C20:2 cis 11,14	




Figure 1b. Separation of a 37 component FAME mix (peaks 15-37)

www.agilent.com/chem

This information is subject to change without notice.

© Agilent Technologies, Inc. 2010

Published in UK, October 08, 2010

SI-02178



Agilent Technologies