

Poster Reprint

**ASMS 2022**  
**Poster number TP090**

# Novel Algorithms and Visualizations for Discovering Polymers and Chemical Series in Industrial and Environmental Applications

Jeremy P Koelmel<sup>1</sup>, David Weil<sup>2</sup>, Paul Stelben<sup>1</sup>, David Godri<sup>3</sup>, Jiarong Qi<sup>1</sup>, Emma Rennie<sup>2</sup>, Elizabeth Lin<sup>1</sup>, John A Bowden<sup>4,5</sup>, Carrie McDonough<sup>6</sup>, Doug Postl<sup>2</sup>, and Krystal Pollitt<sup>1</sup>

<sup>1</sup>School of Public Health, Yale University, New Haven, CT, United States;

<sup>2</sup>Agilent Technologies Inc., Santa Clara, CA, United States;

<sup>3</sup>3rd Floor Solutions;

<sup>4</sup>Center for Environmental and Human Toxicology & Department of Physiological Sciences, University of Florida, Gainesville, FL, United States;


<sup>5</sup>Department of Chemistry, University of Florida, Gainesville, FL, United States;

<sup>6</sup>Department of Civil Engineering, Stony Brook University, Stony Brook, NY, United States

Polymers and chemicals that form series are essential in every day products including plastics, emulsifiers, lubricants, and detergents. The physicochemical properties and health impacts of a product are determined by the properties of molecules in each mixture, such as number of repeating units, chemical moieties and degree of unsaturation. Thus, it is essential to be able to characterize these mixtures at the molecular level. The repeating nature of chemical moieties means that mass spectra and chromatographic data of polymers and chemical series have obvious patterns, which when discerned can be used to comprehensively characterize a sample. After the success of the LipidMatch and FluoroMatch workflows for non-targeted analysis which take advantage of these unique spectral patterns, we present a new workflow for polysorbates and related species.

## Libraries


Class based libraries with fragmentation rules were generated for over 150,000+ species across 16 different types of polysorbate, polysorbide, and PEG classes, including with and without fatty acids.



Libraries are also generated for data that just contain MS/MS fragments corresponding to PEG fragments, forming over 64,450 PEG species. Below is an example of the naming conventions and identifications in PolySorbate 80.

**Simulated PEGs: 64,450 species**

**Structures and Naming Conventions:**



PEG(n=3\_NH4+\_CH2=6\_DBE=4\_(4-OH))

**Polysorbate 80:**

**580+ unique structures**

**217 with at least one NH<sub>3</sub>**

**74 with ketones**

**332 as NH<sub>4</sub><sup>+</sup> adducts...**




PEG(n=4\_NH4+\_CH2=0\_DBE=0\_)

Example	Definition	Built-in Range
n=3	Number of PEG units	3 - 44
NH4+	Adduct	NH4+, Na+, H+
CH2=6	Number of carbons (non-PEG)	0 - 19
DBE=4	Double bond equivalents	0 - 6
(4-OH)	Functional groups	OH (0-4), NH3 (0-2), =O (0-2)

**PolyMatch Flow and PolyMatch Modular (acquisition and data-processing workflow)**

The PolyMatch software data analysis workflow starts by importing data collected using MS, and MS/MS data dependent (DDA), iterative exclusion MS/MS, or targeted MS/MS modes from individual, pooled and blank samples. PolyMatch algorithms cover file -conversion, blank filtering, feature annotation, and visualization via PolyVI (visualization software). One can also generate your own libraries using PolyMatch Generator from your own series of compounds. The PolyMatch Software also directly imports data processed initially using Agilent's Mass Profiler software or other peak picking software.



**Annotation Evidence Compiled by PolyMatch**

PolyMatch integrates a wide range of evidence to classify polymers: mass defect, retention time, and exact mass can be used alongside homologous series to compile groups of chemicals that likely belong to the same class. MS/MS evidence can give structural information pertaining to class or species level assignments. Evidence used for annotation is shown bottom left and a homologous series retention time vs *m/z* is shown for Polysorbate-PEG bottom right.

