

# Biotage® Extrahera™ Software (SW) Release 2.0.0

Here at Biotage, we want to keep your Extrahera updated with the most efficient, user-friendly and easy to program software. So, every now and then, we issue software updates to stamp out any pesky software bugs or we introduce minor cosmetic changes.

This time we're introducing features based on your feedback and it's at no cost to you. No matter your current SW version, we can help you upgrade to 2.0.0.

## What's in the 2.0.0 SW release?

You'll get all the features shown on the right, plus some minor system improvements. If you're interested in learning more about these features, you can find details in our release notes (P/N 414368-H) or in our Extrahera user's manual (P/N 414158-L). Both can be download from our website ([www.biotage.com](http://www.biotage.com)).

## How do I know what SW version I have?

With your Extrahera powered up, select "About" on the splash screen. Note the "Application SW version" and, for good measure, write down the serial number (just in case you need it).

## How do I upgrade to 2.0.0?

Once you've determined your SW version, simply contact your local Biotage account manager or sales representative and request SW 2.0.0. If your Extrahera SW predates version 1.4.0, you will require our dual flow upgrade (P/N: 415039) before you can upgrade to the 2.0.0.

## I have questions?

If you don't find the information you are looking for, or have additional questions, please contact your Biotage account manager or sales representative.

**We're here to support you.**



Visual Alarms



Audible Alarms



Run Individual Samples



Volatile Solvent Designation



ISOLUTE® HYDRO DME+ (Dual Mode Extraction)



Process 48 x 3 mL Tableless SPE Cartridges



Pipette Pump Calibration Notifications



QR and Barcoding (with purchase of reader)

Literature Number: biomark140

© 2020 Biotage. All rights reserved. No material may be reproduced or published without the written permission of Biotage. Information in this document is subject to change without notice and does not represent any commitment from Biotage. E&OE. A list of all trademarks owned by Biotage AB is available at [www.biotage.com/legal](http://www.biotage.com/legal). Other product and company names mentioned herein may be trademarks or registered trademarks and/or service marks of their respective owners, and are used only for explanation and to the owners' benefit, without intent to infringe.