


Terpenes in Gin by GC-TOFMS

Life Science and Chemical Analysis Centre
LECO Corporation, Saint Joseph, Michigan USA



Spectra of terpenes in gin obtained by spectral deconvolution of TOFMS data.

Column:

10 m x 0.2 mm x 0.2 μ m Rtx-5

TOFMS:

35 to 500 m/z @ 5 spectra/s

For questions on this analysis e-mail us at:
life_science@leco.com

For a complete listing of snapshots and application notes visit us on the web at www.leco.com

