


GC-TOFMS Rapid Drug Screen Analysis of 12 Drugs of Abuse in 4.9 Minutes

Life Science and Chemical Analysis Centre
LECO Corporation, Saint Joseph, Michigan USA



GC-TOFMS analysis of 12 drugs of abuse in 4.9 minutes.

Column:
20 m x 0.18 mm x 0.18 μ m Rxi-5ms

TOFMS:
40 to 450 m/z @ 10 spectra/s

For questions on this analysis e-mail us at:
life_science@leco.com

For a complete listing of snapshots and application notes visit us on the web at www.leco.com

