


Shim-pack GIST Amide Application

Analysis of Tetrodotoxin



Analyte:

1. Tetrodotoxin 100 ng/mL

LC Conditions

Column : Shim-pack GIST Amide (3 μ m, 150 x 2.1 mm I.D.)
P/N: 227-30818-06

Mobile Phase : A) CH₃CN
B) 0.1 % HCOOH in H₂O

Time (min)	A (%)	B (%)
0.00	85	15
8.00	50	50
8.10	85	15
15.1	85	15

Flow Rate : 0.3 mL/min

Col. Temp. : 40 °C

Detection : LC/MS/MS

Injection Vol. : 2 μ L

Sample : Standard